



Science Together

 **KNAUER**

Product Selection Guide

2021/2022



Get in touch

Sales

If you want to learn more about our products and services or get a quote, the experts from our Sales team are happy to assist you with your request.

Phone: +49 30 809727-0 (workdays 9-17h CET)

Fax : +49 30 8015010

Email: sales@knauer.net

Support

Do you have questions about the installation or the operation of your device or software?

International Support:

Contact your local KNAUER partner for support:

www.knauer.net/en/Support/Distributors-worldwide

Support in Germany

Phone: +49 30 809727-111 (workdays 9-17h CET)

Fax : +49 30 8015010

Email: support@knauer.net



Disclaimer

Technical data or prices are subject to change without notice. Prices may vary by country and do not include taxes, customs duties or delivery. All trademarks are the property of their respective owners. Our general terms and conditions apply: www.knauer.net/terms.

Welcome to KNAUER



About KNAUER

Based in Berlin, KNAUER is a medium-sized company that has been serving the sciences since 1962. We develop and manufacture scientific instruments of superior quality for liquid chromatography and other laboratory tasks. The product range includes systems and components for analytical HPLC/UHPLC, preparative HPLC, fast protein liquid chromatography (FPLC), multi-column chromatography/simulated moving bed (SMB), as well as equipment for high pressure dosing and osmometry.

Sustainability & ecological commitment

We are committed to protect the environment for ourselves and our children. KNAUER contributes to the conservation of a healthy environment by basing our work on an environmental management system according to DIN EN ISO 14001. The KNAUER quality management system according to DIN EN ISO 9001 makes sure that we continuously manufacture products in the best quality possible. As a family business with around 160 employees, KNAUER focuses on sustainability and takes responsibility for our future.

Some of our ecological activities:

- The regular creation of an input and output balance for the determination and evaluation of energy and resource flows
- Environmentally friendly product development, energy-efficient production, and shipping with biodegradable packaging materials and re-usable packaging with local suppliers
- Fixed specifications for the development of new products according to ecological aspects such as low solvent consumption, repairability, and longevity of the products
- Complete modernization of the company building included thermal insulation, new windows, electric blinds, and a green rooftop, which resulted in a 50 % heating energy saving
- 100 % green electricity and generation of solar power with our photovoltaic system on the roof
- Guidelines for business travel from an environmental, economic, and social perspective
- Tips and instructions for clients to reduce solvent consumption during instrument use
- Environmentally compatible working and manufacturing of HPLC instruments and accessories, e.g. by using energy-efficient working equipment and reducing the use of solvents and harmful substances
- A life cycle assessment to optimize the manufacturing process and concentrate on electricity saving components

Sustainability: #KNAUERforFuture

Many KNAUER employees have good ideas for sustainability, and so we all get better together every year. We would like to inspire YOU to implement sustainability in many areas of your company, too. May these short videos keep you entertained and invite you to act! www.knauer.net/sustainable.

Table of contents

Instruments

| | |
|---------------------------------|----|
| Pumps | 1 |
| Assistants | 11 |
| Autosampler | 13 |
| Column thermostat | 15 |
| Detectors | 21 |
| Fraction collectors | 27 |
| Degasser | 30 |
| Valves unifier and valves | 31 |
| Osmometer | 33 |

Accessories

| | |
|---|----|
| Pump accessories | 34 |
| Autosampler accessories - spare parts | 41 |
| Detector accessories | 44 |
| Valve accessories - spare parts | 48 |
| Osmometry accessories | 52 |
| Purification accessories | 53 |
| Consumables | 54 |
| Lab equipment | 70 |

PC Hardware & software

| | |
|---|----|
| Mobile Control (Chrom) for Windows 10 | 75 |
| ClarityChrom® | 76 |
| PurityChrom® | 78 |
| Purity Chrom® MCC / MCC PLUS | 79 |
| OpenLAB® CDS EZChrom Edition | 80 |
| Chromeleon™ 7.2 Drivers | 81 |
| PC Hardware & periphery | 82 |

KNAUER Services 83

| | |
|--|----|
| Application Services | 83 |
| Chiral Column Screening Services | 84 |
| KNAUER Academy | 84 |
| Research | 86 |
| Compliance | 86 |
| Support | 88 |
| Development Services | 89 |
| Storage of instruments and systems | 89 |
| Configuration of your PC | 89 |

Power cable overview 90

System configurator 92

| | |
|-------------------------------|----|
| Analytic HPLC/UHPLC | 92 |
| Preparative HPLC | 93 |
| FPLC | 94 |

Bluishadow Devices 95

KNAUER GMP Services 96

Conversion tables 97

Terms & Conditions 98

KNAUER Brochures 102

AZURA® Pump P 6.1L

The AZURA® Pump P 6.1L uses technology to overcome the challenges of pumping LC solvents at high pressure with low pulsation. This pump is designed to fulfil the needs for high pressure and low pressure mixing tasks. The pump can deliver flow in the range of 0.001 - 50 ml/min at pressures up to 1000 bar (depending on model and flow rate). The AZURA® binary pump contains two identical high-pressure pumps, a 2 × 2-channel solvent selection valve and the new developed AZURA® mixer, a low-volume microfluidic mixing device. The AZURA® quaternary pump contains one high-pressure pump and an integrated LPG mixing block with a 4 channel valve and mixer. The integrated degasser and AZURA® inline filter are completing the analytical AZURA® HPLC pump and turn this pump into a working horse in the lab. This pump is also available with wetted materials made from ceramic, PEEK and titanium for bio-compatible applications.



Specifications

Solvent delivery

| | |
|------------------------------------|--|
| Pump type | Analytical HPLC pump |
| Delivery system | Double Serial Piston Pump |
| Pulsation compensation | Active Pulsation Compensation |
| Piston seal washing | Active Wash |
| Flow rate accuracy | < 1%, measured at 5 - 50 % of max. flow range using ethanol/water 10:90 |
| Flow rate precision | ≤ 0.1% RSD |
| System protection | Soft start, Pmin und Pmax are programmable |
| Gradient range | 0 - 100 % in 0.1 % increments |
| Solvent selection valve | HPG only |
| Gradient formation | LPG / HPG |
| Liquid temperature range | 4–60 °C (39.2–140 °F) |
| HPG: gradient accuracy | ± 0.3 % at 1 ml/min, 150 bar (ethanol/caffeine tracer) ± 1 % (5 - 95 %, measured at 0.1 - 10 ml/min, water/caffeine tracer) |
| HPG: gradient precision | < 0.1 % RSD at 1 ml/min, 0.3% RSD overall, based on retention time at constant room temperature |
| LPG: gradient accuracy | ± 0.3 % at 1 ml/min, 150 bar (ethanol/caffeine tracer) ± 2 % (1 - 99 %, measured at 5 - 50 % of the flow range, water/caffeine tracer) |
| LPG: gradient precision | < 0.1 % RSD at 1 ml/min, 0.5 % RSD overall, based on retention time at constant room temperature |
| Pump head inlet (standard) | UNF 1/4-28 Thread (for 1/8" tubing) |
| Pump head outlet (standard) | UNF 10-32 Thread (for 1/16" capillary) |

Degasser module

| | |
|---------------------------------|---|
| Degasser channels | 4 channels (LPG Versions), 2 / 4 channels (HPG Versions); optional |
| Max. flow rate/channel | 10 ml/min |
| Degassing method | Gas Permeation through Teflon® AF amorphous fluoropolymer membrane |
| Degassing efficiency | < 0.5 ppm dissolved O ₂ at 1 ml/min |
| Degassing chamber volume | 280 µl volume per channel |
| Solvent applicability | Universal, with exception of hydrochloric acid and halogenated hydrocarbons |
| Wetted materials | PEEK, Tefzel® (ETFC), Systec AF™ |

Communication

| | |
|-----------------------------|--|
| Display | Mobile Control (optional) |
| Inputs | LAN, Pin header connectors (Analog IN, Start IN, Error IN) |
| Analog inputs | Flow rate, 0 - 10 V via pin header connectors |
| Analog control input | Flow Rate |
| Level/event outputs | 8 event outputs (TTL, OC, Relais) and 24 V |
| Control | LAN, Analog and event control, Mobile Control |

Technical parameters

| | |
|---------------------------|---|
| Leak sensor | Yes |
| Special features | Pump Head is detected automatically using Radio frequency identification (RFID) |
| Ambient conditions | 4-40 °C (39.2-104 °F) Air humidity below 90 %, non-condensing |

General

| | |
|---------------------|---|
| Power supply | 100 - 240 V; 50 - 60 Hz; Maximum power consumption 100 Watt |
| Dimensions | 361 mm x 208.2 mm x 523 mm (W × H × D) |
| Weight | 14.1 kg |



KNAUER offers various software control options:

www.knauer.net/softwarecontrol



For pump accessories see p. 34

AZURA® Pump P 6.1L with 5 ml pump head

Pump specifications

| | |
|-------------------------------|--|
| Pump head | 5 ml |
| Continuous working conditions | 0.1 - 4 ml/min |
| Best working conditions | 0.02 - 5 ml/min |
| Flow rate increment | 0.001 ml/min |
| Mixing volume | 100 µl (HPG) |
| Wetted materials | GFP, Stainless Steel, FKM, PEEK, Sapphire, Aluminiumoxide, Ruby, Zirconiumoxide |
| Maximum delivery pressure | 14500 psi / 1000 bar / 100 MPa up to 2 ml/min, 10150 psi / 700 bar / 70 MPa up to 5 ml/min |
| Flow rate range | 0.001 - 5 ml/min |
| Pump head material | Stainless Steel |

Ordering details:

Device

| | |
|---------|--|
| APH34GA | AZURA® Pump P 6.1L (LPG), with 5 ml pump head (stainless steel), degasser and mixer (100 µl) |
| APH35GA | AZURA® Pump P 6.1L (HPG), with 5 ml pump head (stainless steel), degasser and mixer (100 µl) |

AZURA® Pump P 6.1L with 10 ml pump head

Pump specifications

| | |
|-------------------------------|--|
| Pump head | 10 ml |
| Continuous working conditions | 0.1 - 4.0 ml/min |
| Best working conditions | 0.1 - 8.0 ml/min |
| Flow rate increment | 0.001 ml/min |
| Mixing volume | 100 µl HPG, 200 µl LPG |
| Wetted materials | GFP, Stainless Steel, FKM, PEEK, Sapphire, Aluminiumoxide, Ruby, Zirconiumoxide |
| Maximum delivery pressure | 12500 psi / 862 bar / 86 MPa up to 2 ml/min; 5800 psi / 400 bar / 40 MPa up to 10 ml/min |
| Flow rate range | 0.001 - 10 ml/min |
| Pump head material | Stainless Steel |

Ordering details:

| | |
|---------|---|
| APH30EA | AZURA® Pump P 6.1L isocratic, without degasser, with 10 ml pump head (stainless steel) |
| APH31EA | AZURA® Pump P 6.1L isocratic, with degasser, with 10 ml pump head (stainless steel) and solvent selection valve |
| APH30ED | AZURA® Pump P 6.1L isocratic, without degasser, with 10 ml NP pump head (stainless steel) |
| APH34EA | AZURA® Pump P 6.1L (LPG), with 10 ml pump head (stainless steel), degasser and mixer (200 µl) |
| APH35EA | AZURA® Pump P 6.1L (HPG), with 10 ml pump head (stainless steel), degasser and mixer (100 µl) |
| APH35ED | AZURA® Pump P 6.1L (HPG), with 10 ml NP pump head (stainless steel), degasser and mixer (100 µl) |
| APH38EA | AZURA® Pump P 6.1L (HPG), without degasser, with 10 ml pump head and mixer (100 µl) |
| APH38ED | AZURA® Pump P 6.1L (HPG), without degasser, with 10 ml NP pump head and mixer (100 µl) |
| APH39EA | AZURA® Pump P 6.1L (LPG), without degasser, with 10 ml pump head (stainless steel) and mixer (200 µl) |

AZURA® Pump P 6.1L with 50 ml pump head

Pump specifications


| | |
|-------------------------------|--|
| Pump head | 50 ml |
| Continuous working conditions | 0.1 - 20 ml/min |
| Best working conditions | 0.1 - 40 ml/min |
| Flow rate increment | 0.001 ml/min |
| Mixing volume | 200 µl |
| Wetted materials | GFP, FKM, PEEK, Sapphire, Aluminiumoxide, Ruby, Zirconiumoxide |
| Maximum delivery pressure | 4350 psi / 300 bar / 30 MPa up to 10 ml/min; 2900 psi / 200 bar / 20 MPa up to 50 ml/min |
| Flow rate range | 0.01 - 50 ml/min |
| Pump head material | Stainless Steel |

Ordering details:

| | |
|---------|---|
| APH30FA | AZURA® Pump P 6.1L isocratic, without degasser, with 50 ml pump head (stainless steel) |
| APH30FD | AZURA® Pump P 6.1L isocratic, without degasser, with 50 ml normal phase pump head (stainless steel) |
| APH38FA | AZURA® Pump P 6.1L (HPG), without degasser, with 50 ml pump head (stainless steel) and mixer (200 µl) |

AZURA® Pump P 6.1L Biocompatible

Pump specifications

| | |
|--|--|
| Pump head | 10 ml / 50 ml |
| Continuous working conditions | for 10 ml pump heads: 0.1 - 4.0 ml/min; for 50 ml pump heads: 0.1 - 20 ml/min |
| Best working conditions | for 10 ml pump heads: 0.1 - 8.0 ml/min; for 50 ml pump heads: 0.1 - 40 ml/min |
| Flow rate increment | 0.001 ml/min |
| Mixing volume | 250 µl |
|  Wetted materials | UHMW PE, PEEK, Sapphire, Aluminiumoxide, Ruby |
| Maximum delivery pressure | 5800 psi / 400 bar / 40 MPa for 10 ml head, 2900 psi / 200 bar / 20 MPa for 50 ml head |
| Flow rate range | 0.001 - 10 ml/min / 0.01 - 50 ml/min |
| Pump head material | Ceramic |

Ordering details:

| | |
|---------|---|
| APH60EB | AZURA® Pump P 6.1L, isocratic, without degasser, with 10 ml pump head (ceramic) |
| APH60FB | AZURA® Pump P 6.1L, isocratic, without degasser, with 50 ml pump head (ceramic) |
| APH64EB | AZURA® Pump P 6.1L (LPG), with 10 ml pump head (ceramic), degasser and mixer (250 µl) |
| APH69EB | AZURA® Pump P 6.1L (LPG), without degasser, with 10 ml pump head (ceramic) and mixer (250 µl) |
| APH65EB | AZURA® Pump P 6.1L (HPG), with degasser, with 10 ml pump head (ceramic) and mixer (250 µl) |
| APH68EB | AZURA® Pump P 6.1L (HPG), without degasser, with 10 ml pump head (ceramic) and mixer (250 µl) |
| APH68FB | AZURA® Pump P 6.1L (HPG), without degasser, with 50 ml pump head (ceramic) and mixer (250 µl) |

AZURA® Pump P 2.1L

AZURA® preparative HPLC pump P 2.1L covers wide flow rate range and pressure capabilities. It has been designed for purification of mg and gram samples. The pump can deliver flow in the range of 0.01 - 1000 ml/min at pressures up to 400 bar (depending on model). The integrated automatic recognition of the pump head with RFID technology allows fast adaptations of the pump for various applications.



Specifications

Solvent delivery

| | |
|---------------------------------|---|
| Pump type | Preparative HPLC pump |
| Delivery system | Dual Piston Pump with pistons parallel |
| Pulsation compensation | Yes, with compressibility factor |
| Piston seal washing | Active Wash |
| Flow rate accuracy | ± 2 %, measured at 5 - 50 % of flow range using ethanol/water 10:90 |
| Flow rate precision | < 0.1% RSD |
| System protection | Soft start, Pmin und Pmax are programmable |
| Gradient range | 0 - 100 % |
| Gradient formation | LPG / HPG |
| Liquid temperature range | 4–60 °C (39.2–140 °F) |
| HPG: gradient accuracy | ± 2 % (5 - 95 %, measured at 5 - 50 % of flow range, water/caffeine tracer) |
| Leak management | Yes |
| HPG: gradient precision | < 1 % RSD based on retention time at constant room temperature |
| LPG: gradient accuracy | ± 3 % (5 - 95 %, measured at 5 - 50 % of flow range, water/caffeine tracer) |
| LPG: gradient precision | 2 % RSD, based on retention time at constant room temperature |

Communication

| | |
|-----------------------------|--|
| Display | Mobile Control (optional) |
| Inputs | LAN, Pin header connectors (Analog IN, Start IN, Error IN) |
| Analog inputs | Flow rate, 0 - 10 V via pin header connectors |
| Analog control input | Flow rate |
| Level/event outputs | 8 event outputs (TTL, OC, Relais) and 24 V |
| Control | LAN, Analog and event control, Mobile Control |

Technical parameters

| | |
|---------------------------|---|
| Leak sensor | Yes |
| Special features | Pump Head is detected automatically using Radio frequency identification (RFID) |
| Ambient conditions | 10 - 40 °C (50-104 °F), Air humidity below 90%, non-condensing |

General

| | |
|-----------------------------|---|
| Power supply | 100 - 240 V; 50 - 60 Hz; Maximum power consumption 320 Watt |
| Dimensions | 361 mm x 208.2 mm x 523 mm (W x H x D) |
| Weight | 19 kg |
| Optional accessories | Ternary low pressure gradient valve block, 10 - 220 ml/min, binary low pressure gradient valve block, 10 - 800 ml/min, pump head heating and cooling device |



KNAUER offers various software control options:

www.knauer.net/softwarecontrol



For pump accessories see p. 34

AZURA® Pump P 2.1L with 100 ml pump head

Pump specifications

| | |
|-------------------------------|--|
| Pump head | 100 ml |
| Continuous working conditions | 1 - 40 ml/min |
| Best working conditions | 1 - 80 ml/min |
| Flow rate increment | 0.01 ml/min |
| Wetted materials | Zirconium oxide (ZrO ₂), FFKM, graphite fiber reinforced PTFE, PEEK, sapphire, ruby, stainless steel, titanium |
| Maximum delivery pressure | 5800 psi / 400 bar / 40 MPa |
| Flow rate range | 0.01 - 100 ml/min |
| Pump head material | Stainless steel / titanium |

Ordering details:

| | |
|---------|--|
| APE20KA | AZURA® Pump P 2.1L with 100 ml pump head (stainless steel) |
| APE20KB | AZURA® Pump P 2.1L with 100 ml pump head (titanium) |

AZURA® Pump P 2.1L with 250 ml pump head

Pump specifications

| | |
|-------------------------------|--|
| Pump head | 250 ml |
| Continuous working conditions | 2.5 - 100 ml/min |
| Best working conditions | 2.5 - 200 ml/min |
| Flow rate increment | 0.01 ml/min |
| Wetted materials | Zirconium oxide (ZrO ₂), FFKM, graphite fiber reinforced PTFE, PEEK, sapphire, ruby, stainless steel, titanium |
| Maximum delivery pressure | 3260 psi / 225 bar / 22.5 MPa up to 100 ml/min, 2900 psi / 200 bar / 20 MPa up to 250 ml/min |
| Flow rate range | 0.01 - 250 ml/min |
| Pump head material | Stainless steel / titanium |

Ordering details:

| | |
|---------|--|
| APE20LA | AZURA® Pump P 2.1L with 250 ml pump head (stainless steel) |
| APE20LC | AZURA® Pump P 2.1L with 250 ml pump head (titanium) |

AZURA® Pump P 2.1L with 500 ml pump head

Specifications

Pump specifications

| | |
|-------------------------------|--|
| Pump head | 500 ml |
| Continuous working conditions | 5 - 200 ml/min |
| Best working conditions | 5 - 400 ml/min |
| Flow rate increment | 0.1 ml/min |
| Wetted materials | Zirconium oxide (ZrO ₂), FFKM, graphite fiber reinforced PTFE, PEEK, sapphire, ruby, stainless steel, titanium |
| Maximum delivery pressure | 1450 psi / 100 bar / 10 MPa |
| Flow rate range | 0.01 - 500 ml/min |
| Pump head material | Stainless steel / titanium |

Ordering details:

| | |
|---------|--|
| APE20MA | AZURA® Pump P 2.1L with 500 ml pump head (stainless steel) |
| APE20MC | AZURA® Pump P 2.1L with 500 ml pump head (titanium) |

AZURA® Pump P 2.1L with 1000 ml pump head

Pump specifications

| | |
|--------------------------------------|--|
| Pump head | 1000 ml |
| Continuous working conditions | 10 - 400 ml/min |
| Best working conditions | 10 - 800 ml/min |
| Flow rate increment | 0.1 ml/min |
| Wetted materials | Zirconium oxide (ZrO ₂), FFKM, graphite fiber reinforced PTFE, PEEK, sapphire, ruby, stainless steel, titanium |
| Maximum delivery pressure | 1080 psi / 75 bar / 7.5 MPa up to 350 ml/min, 720 psi / 50 bar / 5 MPa up to 1000 ml/min, |
| Flow rate range | 1 - 1000 ml/min |
| Pump head material | Stainless steel / titanium |

Ordering details:

| | |
|---------|---|
| APE20NA | AZURA® Pump P 2.1L with 1000 ml pump head (stainless steel) |
| APE20NB | AZURA® Pump P 2.1L with 1000 ml pump head (titanium) |

LPG Modules

Ordering details:

| | |
|---------|---|
| AZZ00AA | LPG module for Pump P 2.1L binary up to 800 ml/min (stainless steel) |
| AZZ00AB | LPG module for Pump P 2.1L ternary up to 220 ml/min (stainless steel) |
| AZZ10AB | LPG module for Pump P 2.1L ternary up to 220 ml/min (PEEK) |

AZURA® Pump P 4.1S

The AZURA® Pump P 4.1S was developed for high-pressure dosing applications of up to 400 bar and for flow rates of up to 50 ml/min. Whenever a compact and easy-to-integrate pump is required, this pump is a perfect choice.

The pump contains a manual purge valve with a built-in pressure sensor. The pump automatically stops the flow when minimum or maximum pressure limits are reached. The exchangeable pump heads are compatible with a wide range of chemicals and the versatile control options allow easy remote and standalone operation.



Specifications

Solvent delivery

| | |
|------------------------------------|---|
| Pump type | Ultra-compact high pressure pump |
| Delivery system | Dual piston pump with one working piston, one auxillary piston |
| Pulsation compensation | No |
| Piston seal washing | Passive Wash |
| Flow rate accuracy | ± 2 %, measured at 5 - 50 % of flow range using ethanol/water 10:90 |
| Flow rate precision | ≤ 0.5% RSD, measured at 1/5 ml/min using ethanol/water 10:90 |
| System protection | Pmin und Pmax are programmable |
| Liquid temperature range | 4–60 °C (39.2–140 °F) |
| Pump head inlet (standard) | 1/8" OD, 2.1 mm ID FEP tubing (UNF 1/4-28 thread, flat bottom) |
| Pump head outlet (standard) | UNF 10-32 Thread (for 1/16" capillary) |

Communication

| | |
|-----------------------------|--|
| Display | Yes |
| Inputs | LAN, Pin header connectors (Analog IN, Start IN, Error IN), RS-232 |
| Analog inputs | 0 - 10 V |
| Analog control input | Flow rate |
| Control | LAN, RS-232, analog, standalone |

Technical parameters

| | |
|---------------------------|--|
| Display | Yes |
| Ambient conditions | 10-40 °C (50-104 °F) Air humidity below 90 %, non-condensing |

General

| | |
|---------------------|---|
| Power supply | 100 - 240 V; 50 - 60 Hz; Maximum power consumption 100 Watt |
| Dimensions | 121 x 129 x 220 mm (W x H x D) |
| Weight | 2.4 kg |



KNAUER offers various software control options:
www.knauer.net/softwarecontrol



For pump accessories
 see p. 34

AZURA® Pump P 4.1S with 10 ml pump head

Pump specifications

| | |
|-------------------------------|--|
| Pump head | 10 ml |
| Flow rate range | 0.001 - 10 ml/min |
| Maximum delivery pressure | 5800 psi / 400 bar / 40 MPa up to 10 ml/min |
| Wetted materials | Graphite fiber reinforced PTFE, FKM (FFKM for APG20EC), PEEK (PCTFE for APG20EC), sapphire, ruby, zirconium oxide, titanium and pump head material |
| Maximum viscosity | 100 mPa s (at reduced max. flow) |
| Flow rate increment | 0.001 ml/min |
| Best working conditions | 0.1 - 8.0 ml/min |
| Continuous working conditions | 0.1 - 4.0 ml/min |
| Pump head material | Stainless steel / ceramic / Hastelloy® C |

Ordering details:

| | |
|---------|---|
| APG20EA | AZURA® Pump P 4.1S compact HPLC pump with 10 ml/min stainless steel pump head, stainless steel connections |
| APG20EB | AZURA® Pump P 4.1S compact HPLC pump with 10 ml/min ceramic pump head, PEEK connections |
| APG20EC | AZURA® Pump P 4.1S compact HPLC pump with 10 ml/min Hastelloy® C pump head, Titanium pressure sensor, Hastelloy® C connections |
| APG20EF | AZURA® Pump P 4.1S compact HPLC pump with 10 ml/min ceramic pump head, Ti connections. |
| APG20EG | AZURA® Pump P 4.1S compact HPLC pump with 10 ml/min stainless steel pump head, stainless steel connections, recommended for aqueous solutions |
| APG20EH | AZURA® Pump P 4.1S compact HPLC pump with 10 ml/min ceramic pump head, Titanium connections, recommended for aqueous solutions |

AZURA® Pump P 4.1S with 50 ml pump head

Pump specifications

| | |
|-------------------------------|--|
| Pump head | 50 ml |
| Flow rate range | 0.01 - 50 ml/min |
| Maximum delivery pressure | 2180 psi / 150 bar / 15 MPa up to 50 ml/min |
| Wetted materials | Graphite fiber reinforced PTFE, FKM (FFKM for APG20FC), PEEK (PCTFE for APG20FC), sapphire, ruby, zirconium oxide, titanium and pump head material |
| Maximum viscosity | 100 mPa s (at reduced max. flow) |
| Flow rate increment | 0.01 ml/min |
| Best working conditions | 0.5 - 40.0 ml/min |
| Continuous working conditions | 0.5 - 20 ml/min |
| Pump head material | Stainless steel / ceramic / Hastelloy® C |

Ordering details:

| | |
|---------|---|
| APG20FA | AZURA® Pump P 4.1S compact HPLC pump with 50 ml/min stainless steel pump head, stainless steel connections |
| APG20FB | AZURA® Pump P 4.1S compact HPLC pump with 50 ml/min ceramic pump head, PEEK connections |
| APG20FC | AZURA® Pump P 4.1S compact HPLC pump with 50 ml/min Hastelloy® C pump head, Hastelloy® C connections |
| APG20FG | AZURA® Pump P 4.1S compact HPLC pump with 50 ml/min stainless steel pump head, stainless steel connections, recommended for aqueous solutions |
| APG20FI | AZURA® Pump P 4.1S compact HPLC pump with 50 ml/min ceramic pump head, PEEK connections, recommended for aqueous solutions |

AZURA® Pump P 2.1S

The AZURA® Pump P 2.1S was developed for high-pressure dosing applications of up to 400 bar and for flow rates of up to 50 ml/min. Whenever a compact and easy-to-integrate pump is required, this pump is a perfect choice.

The exchangeable pump heads are compatible with a wide range of chemicals and the versatile control options allow easy remote and standalone operation. For aggressive liquids, a Hastelloy® C version is available.



KNAUER offers various software control options:
www.knauer.net/softwarecontrol



For pump accessories
see p. 34

Specifications

Solvent delivery

| | |
|------------------------------------|--|
| Pump type | Ultra-compact high pressure pump |
| Delivery system | Dual piston pump with one working piston, one auxiliary piston |
| Pulsation compensation | No |
| Piston seal washing | Passive Wash |
| Flow rate accuracy | ± 5%, measured at 5 - 50% of flow range using ethanol/water 10:90. ± 2 % at calibration point (one point calibration), measured at 5 - 50% of flow range |
| Flow rate precision | ≤ 0.5 % RSD, measured at 1/5 ml/min using ethanol/water 10:90 |
| System protection | Imin und Imax are programmable (I ~ pressure) |
| Liquid temperature range | 4–60 °C (39.2–140 °F) |
| Pump head inlet (standard) | 1/8" OD, 2.1 mm ID FEP tubing (UNF 1/4-28 thread, flat bottom) |
| Pump head outlet (standard) | UNF 10-32 Thread (for 1/16" capillary) |

Communication

| | |
|-----------------------------|--|
| Display | Yes |
| Inputs | LAN, Pin header connectors (Analog IN, Start IN, Error IN), RS-232 |
| Analog inputs | 0 - 10 V |
| Analog control input | Flow rate |
| Control | LAN, RS-232, analog, standalone |

Technical parameters

| | |
|---------------------------|--|
| Display | Yes |
| Ambient conditions | 10-40 °C (50-104 °F) Air humidity below 90 %, non-condensing |

General

| | |
|---------------------|---|
| Power supply | 100 - 240 V; 50 - 60 Hz; Maximum power consumption 100 Watt |
| Dimensions | 121 x 129 x 220 mm (W x H x D) |
| Weight | 2.3 kg |

AZURA® Pump P 2.1S with 10 ml pump head

Pump specifications

| | |
|--------------------------------------|--|
| Pump head | 10 ml |
| Flow rate range | 0.001 - 10 ml/min |
| Maximum delivery pressure | 5800 psi / 400 bar / 40 MPa up to 10 ml/min |
| Wetted materials | Graphite fiber reinforced PTFE, FKM (FFKM for APG90EC), PEEK (PCTFE for APG90EC), sapphire, ruby, zirconium oxide and pump head material |
| Maximum viscosity | 100 mPa s (at reduced max. flow) |
| Flow rate increment | 0.001 ml/min |
| Best working conditions | 0.1 - 8.0 ml/min |
| Continuous working conditions | 0.1 - 4.0 ml/min |
| Pump head material | Stainless steel / ceramic / Hastelloy® C |

Ordering details:

| | |
|---------|--|
| APG90EA | AZURA® Pump P 2.1S compact HPLC pump with 10 ml/min stainless steel pump head |
| APG90EB | AZURA® Pump P 2.1S compact HPLC pump with 10 ml/min ceramic pump head |
| APG90EC | AZURA® Pump P 2.1S compact HPLC pump with 10 ml/min Hastelloy® C pump head |
| APG90EG | AZURA® Pump P 2.1S compact HPLC pump with 10 ml/min stainless steel pump head, recommended for aqueous solutions |

AZURA® Pump P 2.1S with 50 ml pump head

Pump specifications

| | |
|--------------------------------------|--|
| Pump head | 50 ml |
| Flow rate range | 0.01 - 50 ml/min |
| Maximum delivery pressure | 2180 psi / 150 bar / 15 MPa up to 50 ml/min |
| Wetted materials | Graphite fiber reinforced PTFE, FKM (FFKM for APG90FC), PEEK (PCTFE for APG90FC), sapphire, ruby, zirconium oxide and pump head material |
| Maximum viscosity | 100 mPa s (at reduced max. flow) |
| Flow rate increment | 0.01 ml/min |
| Best working conditions | 0.5 - 40.0 ml/min |
| Continuous working conditions | 0.5 - 20 ml/min |
| Pump head material | Stainless steel / ceramic / Hastelloy® C |

Ordering details:

| | |
|---------|--|
| APG90FA | AZURA® Pump P 2.1S compact HPLC pump with 50 ml/min stainless steel pump head |
| APG90FB | AZURA® Pump P 2.1S compact HPLC pump with 50 ml/min ceramic pump head |
| APG90FC | AZURA® Pump P 2.1S compact HPLC pump with 50 ml/min Hastelloy® C pump head |
| APG90FG | AZURA® Pump P 2.1S compact HPLC pump with 50 ml/min stainless steel pump head, recommended for aqueous solutions |

AZURA® Assistant ASM 2.2L

Docking station for pumps, valves and detectors

The Assistant ASM 2.2L is a docking station for up to three compact devices. Valves, pumps and UV detectors can be combined in one housing.


The plug-in modules are removed by loosening four screws allowing the user to exchange modules in case of service within minutes. Likewise, the configuration of the LC system can be adapted to new requirements. Routine maintenance work e.g. replacing the lamp of a detector are easily performed by the user.

Depending on the integrated modules the assistant fulfills many different tasks like eluent delivery, detection, sample and solvent selection, sample injection, column switching or fraction collection. An assistant including a pump, injection valve, and detector features a complete, compact chromatographic system. As a part of a larger system, the ASM 2.2L is extremely versatile in analytical, preparative and continuous liquid chromatography.

Select your desired plug-in modules for the left, middle and right position in the assistant and you will get your perfect assistant for chromatography and beyond.

www.knauer.net/assistants

Freely combine pumps, valves and detectors in one housing



KNAUER offers various software control options:
www.knauer.net/softwarecontrol

Specifications

General

| | |
|--------------------|---|
| Power supply | 100 - 240 V, 50 - 60 Hz, maximum 130 W |
| Dimensions | 361 x 208 x 523 mm (W x H x D) |
| Weight | About 17 kg (depending on integrated modules) |
| Leak sensor | Yes |
| Ambient conditions | Temperature range: 4 - 40 °C, 39,2 - 104 °F Humidity: 10 - 90 % non-condensing |

Communication

| | |
|---------------|--|
| Interfaces | LAN |
| Control | Mobile Control, Software |
| Inputs | Error (IN), Start (IN), Autozero, 0-10 V Analog (IN) |
| Outputs | Event 1-2, Error (OUT) (OC), + 5 V, + 24 V |
| Analog inputs | Integrator output (detector signal) |

Software functions

Assistant configuration: The ASM 2.2L is supported as complete device. Modules are addressed via the assistant.

| | ClarityChrom® | OpenLAB® | Mobile Control (version 6) |
|-------------------------|---------------|----------|----------------------------|
| Two pumps (independent) | yes | no | yes |
| Fraction valve | one | one | yes, one valve |
| Injection module* | no | no | yes, but part of a method |

Single device configuration: The ASM 2.2L is not supported as device. Integrated modules are addressed as separate devices via IP port.

| | ClarityChrom® | OpenLAB® | PurityChrom® |
|-------------------------|---------------------------|---|---------------------------|
| Two pumps (independent) | no | yes | yes |
| Fraction valve | no | cascading (Multi valve fraction collector) | one |
| Injection module* | yes, but part of a method | yes, fully automatic module with trigger for data acquisition | yes, but part of a method |

* An injection module is a combination of one pump and one 6 port 2 position valve.

Plug-in modules

The Assistant ASM 2.2L can be equipped with following plug-in modules. Use the web-based configurator to customize your assistant: www.knauer.net/assistantconfigurator.

Select the chosen plug-in modules for the left, middle and right position of the assistant to receive the article number of your assistant variant.

Configuration note

An assistant with following configuration is not allowed:

- more than two pump modules - a high-pressure gradient is not supported
- more than one UV detector
- without a plug-in module

Basic device

| | |
|-----------------------|-----|
| ASM 2.2L basic device | AY* |
|-----------------------|-----|

Plug-in modules

| Basic plug-in modules | Article number for ordering individual modules (without the assistant housing)* |
|--|---|
| Empty module | AG2022 |
| AZURA® Valve Unifier VU 4.1** | AWA04 |
| AZURA® UV Detector UVD 2.1S | ADA03XA |
| AZURA® UV Detector UVD 2.1S, fiber optics | ADA07XA |
| Compact pump without pressure sensor | |
| AZURA® Pump P 2.1S, 10 ml, stainless steel | APG92EA |
| AZURA® Pump P 2.1S, 10 ml, Hastelloy C | APG92EC |
| AZURA® Pump P 2.1S, 10 ml, ceramic | APG92EB |
| AZURA® Pump P 2.1S, 50 ml, stainless steel | APG92FA |
| AZURA® Pump P 2.1S, 50 ml, Hastelloy C | APG92FC |
| AZURA® Pump P 2.1S, 50 ml, ceramic | AGP92FB |
| Compact pump with pressure sensor | |
| AZURA® Pump P 4.1S, 10 ml, stainless steel | APG22EA |
| AZURA® Pump P 4.1S, 10 ml, stainless steel, normal phase | APG22ED |
| AZURA® Pump P 4.1S, 10 ml, ceramic | APG22EB |
| AZURA® Pump P 4.1S, 50 ml, stainless steel | APG22FA |
| AZURA® Pump P 4.1S, 50 ml, stainless steel, normal phase | APG22FD |
| AZURA® Pump P 4.1S, 50 ml, ceramic | APG22FB |
| AZURA® Pump P 4.1S, 50 bar, 10 ml, stainless steel | APG12EA |
| AZURA® Pump P 4.1S, 50 bar, 10 ml, ceramic | APG12EB |
| AZURA® Pump P 4.1S, 50 bar, 50 ml, stainless steel | APG12FA |
| AZURA® Pump P 4.1S, 50 bar, 50 ml, ceramic | APG12FB |

* Use the assistant configurator to find your desired plug-in module combination configured in the assistant housing: www.knauer.net/assistantconfigurator.

** Note that valves V 4.1 must be ordered in addition to the valve drive VU 4.1. For valves, see p. 32.

Accessories



| | |
|--|---------|
| Column holder - replacing empty module | AG2022B |
|--|---------|



AZURA® Valve Unifier VU 4.1**



AZURA® Detector UVD 2.1S





Column holder


AZURA® Autosampler AS 6.1L

The Autosampler AS 6.1L can inject from up to 768 positions when equipped with microtiter plates (either high or low formats) or from up to 108 standard 1.5 ml sample vials. The sample carryover is significantly minimized thanks to a highly-effective interior and exterior needle wash procedure. This autosampler is also fast and flexible: one complete sample injection cycle takes less than one minute, including needle wash. Three different injection modes are supported; "full loop filling" (highest precision and reproducibility), "partial loop filling" (variable volumes, e.g. for dilution series) and "µl pickup" (loss-free injection of extremely small sample volumes), allowing the user to optimize sample usage. The headspace pressure function prevents bubbles from forming in the vial during sample uptake. Precolumn derivatization is supported.

For high-pressure injections, the autosampler is equipped with a so-called ILDTM valve (Intermediate Loop Decompression). This valve consists of a rotor-stator combination and includes a central port for depressurizing. For high-pressure applications, the sample loop is depressurized prior to receiving the sample. This way, the sample is not diluted with a solvent. Because the valve is switched extremely fast, pressure spikes are reduced. Analyses are more precise and wear of the column is reduced.



 KNAUER offers various software control options: www.knauer.net/softwarecontrol

 For autosampler accessories see p. 41

Specifications

Sample injection

| | |
|---------------------------|--|
| Autosampler Flow Path | Analytical |
| Maximum back pressure | See device versions |
| Vial/plate dimensions | Well plate dimensions according to ANSI SLAS 4-2004 (formerly ANSI/SBS 4-2004) max. plate/vial height: 47 mm (incl. septa or capmat) |
| Injection volume range | 0.1 µl - 10 ml depending on sample loop |
| Headspace pressure | Built-in compressor, only for sample vials with septum |
| Switching time inj. valve | < 100 ms |
| Piercing needle precision | ± 0.6 mm |
| Sample tray cooling | Optional (4 - 40 °C) |
| Vial detection | Missing vial/well plate detection by sensor |
| Wetted materials | ETFE (buffer & needle tubing), stainless steel (sample needle, valve stator), Vespel (rotor seal), Kel-F (syringe valve), glass (syringe), PTFE (tip of syringe plunger) |

Analytical performance

| | |
|----------------------|--|
| Injection modes | Full loop filling, partial loop filling and microliter pickup; PASA™ (pressure-assisted sample aspiration) |
| Injection precision | Full loop filling: < 0.3 % RSD partial loop injection at injection volumes >5 µl: <0.5 % RSD microliter pickup at injections >5 µl: <1.0 % RSD |
| Injection accuracy | 0.2 µl for 250 µl injection syringe |
| Sample carryover | < 0.0015 % for partial loop (chlorhexidine) < 0.0003 % with extended needle wash (s. Technical Note VTN0004) |
| Injections per vial | Max. 9 injections |
| Injection cycle time | Min. 7 s from the same vial, 14 s from different vials; < 60 s for >100 µl sample injection in all injection modes, incl. 300 µl needle wash |
| Analysis time | Max. 9 h, 59 min, 59 s |

Communication

| | |
|------------|---|
| Inputs | 2 programmable TTL inputs (next injection, freeze, stop) |
| Outputs | 1 programmable relay output (inject marker, auxiliary, alarm) |
| Control | Ethernet (LAN) |
| Interfaces | LAN, analog |

Technical parameters

| | |
|--------------------|---|
| Ambient conditions | Temperature range: 10-40 °C; 50-104 °F; air humidity: 20 - 80 % |
|--------------------|---|

NEW

General

| | |
|--------------|--------------------------------|
| Power supply | 95-240 V AC |
| Dimensions | 361 x 375 x 570 mm (W x H x D) |
| Weight | 30 kg |

Device versions

| | HPLC+ | UHPLC | Bio | Prep |
|---------------------------------|-------------------|-------------------|-------------------|----------|
| Maximum back pressure | 862 bar | 1240 bar | 345 bar | 200 bar |
| Sample needle | 15 µl | 15 µl | 15 µl | 60 µl |
| Dispenser syringe | 250 µl | 250 µl | 250 µl | 2500 µl |
| Buffer tubing | 500 µl | 500 µl | 1000 µl | 2000 µl |
| Sample loop | 100 µl, 0.4 mm ID | 10 µl, 0.18 mm ID | 100 µl, 0.4 mm ID | 10 ml |
| Order number | AAA50AA | AAA10AA | AAA20AA | AAA40AA |
| Order number (cool/heat option) | AAA51AA | AAA11AA | AAA21AA | AAA41AA* |

* also available as biocompatible version: AAA31AA

Ordering details:

| | |
|---------|---|
| AAA50AA | AZURA® Autosampler AS 6.1L analytical HPLC autosampler, 862 bar |
| AAA51AA | AZURA® Autosampler AS 6.1L analytical HPLC autosampler, 862 bar, with sample cooling/heating |
| AAA10AA | AZURA® Autosampler AS 6.1L analytical HPLC autosampler, 1240 bar |
| AAA11AA | AZURA® Autosampler AS 6.1L analytical HPLC autosampler, 1240 bar, with sample cooling/heating |
| AAA20AA | AZURA® Autosampler AS 6.1L analytical HPLC autosampler, 345 bar, with biocompatible flow path |
| AAA21AA | AZURA® Autosampler AS 6.1L analytical HPLC autosampler, 345 bar, with biocompatible flow path and sample cooling/heating |
| AAA31AA | AZURA® Autosampler AS 6.1L preparative HPLC autosampler, 200 bar, with biocompatible flow path and sample cooling/heating |
| AAA40AA | AZURA® Autosampler AS 6.1L preparative HPLC autosampler, 200 bar |
| AAA41AA | AZURA® Autosampler AS 6.1L preparative HPLC autosampler, 200 bar, with sample cooling/heating |

AZURA® Column thermostat CT 2.1

The AZURA® CT 2.1 is a price attractive basic column thermostat. It allows temperature control in the range of 5 °C and 85 °C and thus is appropriate for most HPLC applications. For advanced purification and analysis purposes, the oven can optionally be equipped with an eluent pre-heating cartridge. This ensures even more constant separation conditions leading to higher selectivity and an improved peak shape.

The instrument operates with a microprocessor controlled Peltier element for precise temperature settings. In combination with its high temperature stability, this allows programming of linear as well as non-linear temperature gradients.



Specifications

Thermostating

| | |
|----------------------------|--|
| Heating and cooling system | microprocessor controlled Peltier element for heating and cooling, fan supported 2-way air circulation |
| Temperature range | 5-85 °C |
| Heating/cooling rate | 2 °C/min |
| Temperature accuracy | ± 0.2 °C |
| Temperature stability | ± 0.1 °C |



KNAUER offers various software control options:
www.knauer.net/softwarecontrol

Column compartment

| Column dimensions | max. number | max. length* | max. outer diameter* | matching column |
|-------------------|-------------|--------------|----------------------|-----------------------------------|
| | 8 | 160 mm | 12 mm | 125 mm x 4.6 mm ID with precolumn |
| | 4 | 325 mm | 12 mm | 300 mm x 4.6 mm ID |
| | 1 | 325 mm | 35 mm | 300 mm x 16 mm ID |

*total outer dimensions of the column including screw caps.

| | |
|----------------------|--|
| Dimensions, internal | 90 x 390 x 47 mm (W x H x D) |
| Safety | self-check and auto-calibration at power-on, selectable turn-off temperature |
| Leak sensor | gas sensor, adjustable sensitivity, acoustic signal, turn-off switch |

Communication

| | |
|------------|--|
| Control | optional for stand-alone functionality: Mobile Control |
| Interfaces | LAN Interface |

General

| | |
|--------------|--------------------------------|
| Power supply | 90-230 V, 50-60 Hz, 100 W |
| Dimensions | 150 x 470 x 310 mm (W x H x D) |
| Weight | 8.4 kg |

Other

| | |
|----------------------|---|
| Optional accessories | Cartridge for eluent pre-heating for capillary with an ID of 0.1 or 0.18 mm |
|----------------------|---|

Ordering details:

Device

| | |
|-------|--|
| ATC00 | Column Thermostat AZURA® CT 2.1 for constant temperatures and reproducible results |
|-------|--|

Accessories

| | |
|----------|---|
| A05852-3 | Cartridge for eluent Pre-heating ID 0.1 mm, ~5.5 µl |
| A05852-2 | Cartridge for eluent Pre-heating ID 0.18 mm, ~18 µl |

AZURA® Detector DAD 6.1L

The AZURA® DAD 6.1L is a high-end diode array detector (DAD) which combines outstanding performance with easy handling.

A wide range of easily exchangeable flow cells make this device the right choice for fast, standard analytical, semi-preparative and preparative separations with bio-compatible or stainless steel wetted parts.

State-of-the-art total reflection flow cells (LightGuide technology) are available for this detector providing maximum light throughput (due to total internal reflection) with minimal peak dispersion (due to small cell volume) to guarantee an optimized signal to noise (S/N) ratio.

An optional fiber optics adapter offers the possibility to separate the flow cell spatially from the device and thus provides enhanced security for hazardous, explosive or toxic work processes, as well as protecting the device from leakages at high flow rates.

The newly developed optical unit with KNAUER Polka-Dot technology and intelligent temperature management ensure maximum sensitivity combined with minimal baseline drift over the whole spectrum.

Furthermore, easy frontal access and improved safety features enable effortless lamp replacement. This eases maintenance and guarantees short downtimes.

The DAD 6.1L comes installed with a high brightness deuterium and tungsten halogen lamps, which cover a wavelength range from 190 to 1000 nm.



Key features

- Wide application range
- Large choice of flow cells
- Fiber optics adapter available
- Attractively priced
- Made in Germany



KNAUER offers various software control options:

www.knauer.net/softwarecontrol



For detector accessories see p. 44

Specifications

| Detection | |
|-------------------------|---|
| Detector type | Diode array detector |
| Number of diodes | 1024 |
| Pixel pitch | 0.8 nm/diode |
| Detection channels | 8 (Digital)/4 (Analog) |
| Light source | High brightness deuterium (D ₂) lamp and halogen lamp with integrated GLP chip |
| Wavelength range | 190 - 1000 nm |
| Spectral bandwidth | < 3.5 nm at H _α line (FWHM) /Note: digital bandwidth 1 - 32 nm |
| Slit width | 70 μm |
| Wavelength accuracy | ± 1 nm |
| Noise | ± 3.5 μAU at 254 nm (ASTM E1657-98) |
| Drift | 300 μAU/h at 254 nm (ASTM E1657-98) |
| Linearity | > 2.5 AU at 274 nm (ASTM E1657-98) |
| Maximum data rate | 100 Hz (LAN)/12.5 Hz (analog) |
| Flow cell | Not included (see Accessories / Spare parts) |
| Time constants | 0.00 / 0.01 / 0.02 / 0.05 / 0.1 / 0.2 / 0.5 / 1.0 / 2.0 / 5.0 / 10.0 s |
| Integration time | Automatic |
| Wavelength verification | Internal holmium filter and deuterium lines |
| Leak sensor | Yes |
| Communication | |
| Inputs | Error (IN), Start (IN), Autozero |
| Outputs | Events 1 - 2 (Relay and TTL compatible, respectively), Error (OUT), + 5 V, Valve + 24 V, Valve (OUT) |
| Analog outputs | 4 x 0 - 5 V, 20 bit, offset adjustable |
| Control | Mobile Control, software, event control, analog, terminal protocol |
| Interfaces | LAN (RJ-45), USB (service only), multi-pin connector, analog (RCA cinch connector) |
| Technical parameters | |
| GLP | Detailed report including lamp recognition, operating hours, lamp operating hours, number of lamp ignitions |
| Display | Mobile Control (optional) |
| Ambient conditions | Temperature range: 4 - 40 °C, 39.2 - 104 °F, Humidity: below 90 % non condensing |

General

| | |
|--------------|--------------------------------|
| Power supply | 100 – 240 V, 50 – 60 Hz, 75 W |
| Dimensions | 361 x 158 x 523 mm (W x H x D) |
| Weight | 13.8 kg |

Other


| | |
|------|--|
| Note | Flow cells are not included and must be ordered separately (see Accessories / Spare parts) |
|------|--|

Ordering details:

Device

| | |
|-------|---|
| ADC11 | AZURA® Detector DAD 6.1L Diode array detector DAD 6.1L without flow cell 190 - 1000 nm, incl. test cell |
|-------|---|

Accessories

| | |
|---|--|
| AMC19XA | 10 mm path length, 2µl, 1/16", 50 bar, LightGuide Flow cell cartridge for AZURA® Detector DAD/MWD |
| AMD59XA | 50 mm path length, 6µl, 1/16", 50 bar, High Sensitivity LightGuide Flow cell cartridge for AZURA® Detector DAD/MWD |
| AMC38 | 10 mm path length, 10µl, 1/16", 300 bar, PressureProof Flow cell cartridge for AZURA® Detector DAD/MWD |
| AMB18 | 3 mm path length, 2µl, 1/16", 300 bar, PressureProof Flow cell cartridge for AZURA® Detector DAD/MWD |
| AZL01 | HBST deuterium lamp for AZURA® Detector DAD 6.1L |
| AZL02 | Halogen lamp for AZURA® Detector DAD 6.1L |
| AMKX8KIT | Fiber optics adapter kit for AZURA® Detector DAD/MWD, with fiber optic cables (1x 400 mm and 1x 750 mm) and mounting bracket |
| AMLX8 | Test cell for AZURA® Detector DAD/MWD |
|  AZZ00OC | AZURA® Heat Exchanger for analytical PressureProof flow cells, 14 µl |

AZURA® Detector DAD 2.1L & MWD 2.1L

The AZURA® DAD 2.1L is a highly competitive diode array detector which combines high performance with easy handling at an affordable price.

A wide range of easily exchangeable cartridge flow cells make this device the right choice for fast, standard analytical, semi-preparative and preparative separations with bio-inert or stainless steel wetted parts.

State-of-the-art total reflection flow cells (LightGuide technology) are available for this device providing maximum light throughput (due to total internal reflection) with minimal peak dispersion (due to the small cell volume) to guarantee an optimized S/N ratio.

An optional fiber optics adapter offers the possibility to separate the flow cell spatially from the device and thus provides enhanced security for hazardous, explosive or toxic work processes, as well as protecting the device from leakages at high flow rates.

The newly developed optical unit and intelligent temperature management ensure maximum sensitivity combined with minimal baseline drift.

Furthermore, easy frontal access and improved safety features enable effortless lamp replacement. This eases maintenance and guarantees short downtimes.

The DAD 2.1L comes installed with a deuterium lamp which covers a wavelength range from 190 to 700 nm.



Key features

- Wide application range
- Large choice of flow cells
- Fiber optics adapter available
- Leak management
- Made in Germany



KNAUER offers various software control options:

www.knauer.net/softwarecontrol



For detector accessories see p. 44

Specifications

| Detection | |
|-------------------------|---|
| Detector type | Diode array detector |
| Number of diodes | 256 |
| Pixel pitch | 2 nm/diode |
| Detection channels | 8 (Digital)/4 (Analog) |
| Light source | Deuterium (D ²) lamp with integrated GLP chip |
| Wavelength range | 190 - 700 nm |
| Spectral bandwidth | <10 nm at H _α line (FWHM) /Note: digital bandwidth 1 - 32 nm |
| Slit width | 70 μm |
| Wavelength accuracy | ± 1 nm |
| Noise | ± 5 μAU at 254 nm (ASTM E1657-98) |
| Drift | 400 μAU/h at 254 nm (ASTM E1657-98) |
| Linearity | > 2.0 AU at 274 nm (ASTM E1657-98) |
| Maximum data rate | 100 Hz (LAN)/12.5 Hz (analog) |
| Flow cell | Not included (see Accessories / Spare parts) |
| Time constants | 0.00 / 0.01 / 0.02 / 0.05 / 0.1 / 0.2 / 0.5 / 1.0 / 2.0 / 5.0 / 10.0 s |
| Integration time | Automatic |
| Wavelength verification | Internal holmium filter and deuterium lines |
| Leak sensor | Yes |
| Communication | |
| Inputs | Error (IN), Start (IN), Autozero |
| Outputs | Events 1 - 2 (Relay and TTL compatible, respectively), Error (OUT), + 5 V, Valve + 24 V, Valve (OUT) |
| Analog outputs | 4 x 0 - 5 V, 20 bit, offset adjustable |
| Control | Mobile Control, software, event control, analog, terminal protocol |
| Interfaces | LAN (RJ-45), USB (service only), multi-pin connector, analog (RCA cinch connector) |
| Technical parameters | |
| GLP | Detailed report including lamp recognition, operating hours, lamp operating hours, number of lamp ignitions |
| Display | Mobile Control (optional) |
| Ambient conditions | Temperature range: 4 - 40 °C, 39.2 - 104 °F, Humidity: below 90 % non condensing |
| General | |
| Power supply | 100 - 240 V, 50 - 60 Hz, 75 W |
| Dimensions | 361 x 158 x 523 mm (W x H x D) |
| Weight | 12.2 kg |
| Note | Flow cells are not included and must be ordered separately (see Accessories / Spare parts) |

Ordering details:

| Device | |
|-------------|--|
| ADC01 | AZURA® Detector DAD 2.1L Diode array detector DAD 2.1L without flow cell 190 - 700 nm, incl. test cell |
| ADB01 | AZURA® Detector MWD 2.1L Multiwavelength detector MWD 2.1L, without flow cell 190 - 700 nm, incl. test cell |
| Accessories | |
| AMC19XA | 10 mm path length, 2μl, 1/16", 50 bar, LightGuide Flow cell cartridge for AZURA® Detector DAD/MWD |
| AMD59XA | 50 mm path length, 6μl, 1/16", 50 bar, High Sensitivity LightGuide Flow cell cartridge for AZURA® Detector DAD/MWD |
| AMC38 | 10 mm path length, 10μl, 1/16", 300 bar, PressureProof Flow cell cartridge for AZURA® Detector DAD/MWD |
| AMB18 | 3 mm path length, 2μl, 1/16", 300 bar, PressureProof Flow cell cartridge for AZURA® Detector DAD/MWD |
| A5193 | Deuterium lamp, replacement, for S2520, 10D, 40D, UVD 2.1S, UVD 2.1L, DAD 2.1L, MWD 2.1L |
| AMKX8KIT | Fiber optics adapter kit for AZURA® Detector DAD/MWD, with fiber optic cables (1x 400 mm and 1x 750 mm) and mounting bracket |
| AMLX8 | Test cell for AZURA® Detector DAD/MWD |
| AZZ00OC | AZURA® Heat Exchanger for analytical PressureProof flow cells, 14 μl |

AZURA® Detector UVD 2.1L

The AZURA® UV/VIS Detector UVD 2.1L is a competitively priced HPLC spectrophotometer for routine HPLC applications including fast LC methods. Besides offering excellent technical specifications, this robust detector features a highly flexible and compact design. The UVD 2.1L comes with an installed deuterium lamp which covers a wavelength range from 190 to 750 nm.

Due to a smart design the flow cell is easily accessible and can be changed very quickly. You can choose between a wide range of flow cells for analytical or preparative LC applications with flow rates from 10 µl/min up to 10 l/min.



Specifications

| Detection | |
|----------------------|---|
| Detector type | Variable single wavelength UV detector |
| Detection channels | 1 |
| Light source | Deuterium (D ²) lamp with integrated GLP chip |
| Wavelength range | 190 - 750 nm |
| Spectral bandwidth | 11 nm at H _α line (FWHM) |
| Wavelength accuracy | ± 2.5 nm |
| Wavelength precision | 0.3 nm (ASTM E275-93) |
| Noise | ± 15 µAU at 254 nm (ASTM E1657-98) |
| Drift | 300 µAU/h at 254 nm (ASTM E1657-98) |
| Linearity | > 2.0 AU at 274 nm (ASTM E1657-98) |
| Maximum data rate | 50 Hz (LAN)/20 Hz (Analog) |
| Flow cell | Not included (see Accessories / Spare parts) |
| Time constants | 0.0 / 0.1 / 0.2 / 0.5 / 1.0 / 2.0 / 5.0 / 10.0 s |
| Integration time | Automatic |
| Leak sensor | Yes |
| Communication | |
| Inputs | Error (IN), Start (IN), Autozero, 0 - 10 V Analog (IN) |
| Outputs | Events 1 - 3, + 5 V, 24 V Valve |
| Analog outputs | 1 x 0 - 5 V scalable, 20 bit, offset adjustable |
| Control | Digital: LAN, remote connector/Analog: wavelength control/Manual: Mobile Control (optional) |
| Programming | Timed: wavelength, events, fraction valve, links, wake up (program, link); 9 programs, 50 program lines |
| Technical parameters | |
| GLP | Detailed report incl. lamp recognition, operating hours, lamp operating hours, number of lamp ignitions |
| Display | Mobile Control (optional) |
| Ambient conditions | Temperature range: 4 - 40 °C, 39.2 - 104 °F, Humidity: below 90 % non condensing |
| General | |
| Power supply | 100 - 240 V, 50 - 60 Hz, 65 W |
| Dimensions | 361 x 158 x 523 mm (W x H x D) |
| Weight | 5.9 kg |
| Note | Flow cells are not included and must be ordered separately (see Accessories / Spare parts) |

Key features

- Large choice of flow cells
- Leak management
- 55 years experience
- Made in Germany



KNAUER offers various software control options:
www.knauer.net/softwarecontrol



For detector accessories see p. 44

Ordering details:

| Device | |
|-------------|---|
| ADA01XA | AZURA® Detector UVD 2.1L with deuterium lamp without flow cell, incl. test cell |
| ADA04XA | AZURA® Detector UVD 2.1L Fiber Optics Version with deuterium lamp without flow cell |
| Accessories | |
| A4061XB | 10 mm path length, 10 µl, 1/16", 300 bar, stainless steel, with heat exchanger one sided inlet and outlet, classical KNAUER flow cell |
| A4042 | 3 mm path length, 2 µl, 1/16", stainless steel, classical KNAUER flow cell |
| A5193 | Deuterium lamp, replacement, for S2520, 10D, 40D, UVD 2.1S, UVD 2.1L, DAD 2.1L, MWD 2.1L |
| A4126 | Test cell Holmium Oxid Filter |
| A4146 | Test cell, WG 280 filter stray light |
| A4123 | Test cell |

AZURA® Detector UVD 2.1S

The AZURA® UVD 2.1S is a highly competitive single variable wavelength UV detector for HPLC. It offers excellent technical specifications for routine laboratory work. With its small footprint, it is one of the smallest detectors for HPLC on the market.

The UVD 2.1S comes in the novel small AZURA® housing. The installed deuterium lamp covers a wavelength range from 190 to 500 nm. The UV detector can be controlled with OpenLab EZChrom Edition®, ChromGate®, PurityChrom Bio and ClarityChrom® software, as well as from the front panel (stand-alone operation), via LAN, via RS-232, or through analog input/output; allowing it to be integrated into almost any LC system.

Due to a smart design, the flow cell is easily accessible and can be changed very quickly. Choose between a wide range of flow cells for analytical or preparative LC applications with flow rates from 10 µl/min up to 10 l/min. Also available as a module for AZURA® Assistant ASM 2.2L.



Specifications

| Detection | |
|----------------------|--|
| Detector type | Variable single wavelength UV detector |
| Detection channels | 1 |
| Light source | Deuterium (D ²) lamp with integrated GLP chip |
| Wavelength range | 190 - 500 nm |
| Spectral bandwidth | 13 nm at H _α line (FWHM) |
| Wavelength accuracy | ± 3 nm |
| Wavelength precision | 0.7 nm (ASTM E275-93) |
| Noise | ± 20 µAU at 254 nm (ASTM E1657-98) |
| Drift | 300 µAU/h at 254 nm (ASTM E1657-98) |
| Linearity | > 2.0 AU at 274 nm (ASTM E1657-98) |
| Maximum data rate | 50 Hz (LAN)/20 Hz (Analog)/10 Hz (RS-232) |
| Flow cell | Not included (see Accessories / Spare parts) |
| Time constants | 0.00 / 0.02 / 0.05 / 0.1 / 0.2 / 0.5 / 1.0 / 2.0 s |
| Integration time | Automatic |
| Communication | |
| Inputs | Autozero, Start (IN), Error (either IN or OUT) |
| Outputs | Error (either OUT or IN) |
| Analog inputs | Wavelength 0 - 10 V |
| Analog outputs | 1 x ± 2.5 V scalable, 20 bit |
| Control | Front panel, Mobile Control, software, event control, analog, terminal protocol |
| Interfaces | LAN (RJ-45), RS-232 (SUB-D 9), multi-pin connector, analog (RCA cinch connector) |
| Technical parameters | |
| GLP | Lamp operating hours |
| Display | LED |
| Ambient conditions | Temperature range: 4 - 40 °C, 39.2 - 104 °F, Humidity: below 90 % non condensing |
| General | |
| Power supply | External: input 100 - 240 V, output 24 V DC, 60 W |
| Dimensions | 121 x 129 x 187mm (W x H x D) |
| Weight | 1.5 kg |
| Note | Flow cells are not included and must be ordered separately (see Accessories / Spare parts) |

Key features

- Compact
- Large choice of flow cells
- 55 years experience
- Made in Germany



KNAUER offers various software control options:
www.knauer.net/softwarecontrol



For detector accessories see p. 44

Ordering details:

| Device | |
|-------------|---|
| ADA00 | AZURA® Detector UVD 2.1S with deuterium lamp without flow cell, incl. test cell |
| ADA05 | AZURA® Detector UVD 2.1S Fiber Optics Version with deuterium lamp without flow cell |
| Accessories | |
| A4061XB | 10 mm path length, 10 µl, 1/16", 300 bar, stainless steel, with heat exchanger one sided inlet and outlet, classical KNAUER flow cell |
| A4042 | 3 mm path length, 2 µl, 1/16", stainless steel, classical KNAUER flow cell |
| A4045 | 3 mm path length, 2 µl, 1/16", 30 bar, biocompatible, classical KNAUER flow cell |
| A5193 | Deuterium lamp, replacement, for S2520, 10D, 40D, UVD 2.1S, UVD 2.1L, DAD 2.1L, MWD 2.1L |

AZURA® Detector RID 2.1L

The AZURA® RID 2.1L is a sensitive and competitively priced differential refractometer. It is suitable for detecting compounds with little or no UV activity such as alcohols, sugars, lipids or polymers. This instrument is designed for use in analytical HPLC (high performance liquid chromatography) as well as under certain conditions for GPC (gel permeation chromatography) applications.

The intelligently designed optical unit with advanced temperature control ensures high sensitivity, fast baseline stabilization, and excellent reproducibility. Furthermore, the long-life LED, highly pressure resistant flow cell, improved safety features and enhanced diagnostics functions guarantee easy handling and minimal maintenance. The wide linear dynamic range and 10 ml/min maximum flow rate make the AZURA® RID 2.1L the perfect choice for most laboratory tasks.



Key features

- Temperature controlled optical unit
- Long-life LED
- Pressure resistant flow cell
- 55 years experience
- Made in Germany

Specifications

| Detection | |
|------------------------|--|
| Detector type | Refractive index detector |
| Version | analytical |
| Light source | Long-life LED |
| Detection channels | 1 |
| Refractive index range | 1.00 - 1.75 RIU |
| Noise | ± 2.5 nRIU |
| Drift | 200 nRIU/h |
| Linearity | > 1000 µRIU |
| Flow cell | 5 bar back pressure resistance Flow cell included |
| Max. flow rate | 10 ml/min (pure water) |
| Flow cell volume | 15 µl (43 µl dispersion volume) |
| Wetted materials | Stainless steel / quartz / PTFE |
| Temperature control | OFF, 30 - 55 °C (1 °C increment) |
| Time constants | 0.00 / 0.01 / 0.02 / 0.05 / 0.1 / 0.2 / 0.5 / 1.0 / 2.0 / 5.0 / 10.0 s |
| Maximum data rate | 100 Hz (LAN)/20 Hz (Analog) |
| Autozero | Full range |
| Leak sensor | Yes (internal and external leak management) |
| Communication | |
| Inputs | Error (IN), Start (IN), Autozero, Flush (IN) |
| Outputs | Event 1, Start (OUT), Error (OUT), + 5 V, 24 V Valve |
| Analog outputs | 1 x 0- 2.5 V scalable, 20 bit, offset adjustable |
| Control | Mobile Control, software, event control, analog, terminal protocol |
| Interfaces | 2 x LAN (RJ-45, dual IP-stack), USB (service only), multi-pin connector, analog (cinch connector) |
| Technical parameters | |
| GLP | Detailed report including operating hours, light source operating hours |
| Display | Mobile Control (optional) |
| Ambient conditions | Temperature range: 4 - 40 °C, 39.2 - 104 °F, Humidity: below 90 % non condensing |
| General | |
| Power supply | 100 - 240 V, 50 - 60 Hz, 65 W |
| Dimensions | 361 x 158 x 523 mm (W x H x D) |
| Weight | 10.8 kg |
| Optional accessories | Mobile Control |
| Ordering details: | |
| Device | |
| ADD31 | AZURA® Detector RID 2.1L analytical refractive index detector with flow cell |
| ADD38 | AZURA® Detector RID 2.1L HighFlow preparative refractive index detector with flow cell and external pressure release valve |

KNAUER offers various software control options:
www.knauer.net/softwarecontrol

For detector accessories
 see p. 44

Fluorescence Detector RF-20A/Axs

The fluorescence detector RF-20A provides world-class sensitivity, excellent maintainability and diverse validation / support functions. It supports a wide range of applications in the wavelength range of 200 to 650 nm from conventional analysis to high-performance analysis. With a signal-to-noise ratio of 1200 for the water-Raman band, the fluorescence detector is well suited for trace analysis. The xenon lamp and flow cell are directly accessible on the device, thus allowing a quick and easy handling and maintenance of the device by the user, thereby minimizing downtime. The lamp life is 2000 hours. When replacing the xenon lamp, no adjustment is required.



Specifications

| Detection | |
|----------------------|---|
| Detector type | Fluorescence detector |
| Detection channels | 1 |
| Number of signals | 1 |
| Light source | xenon lamp |
| Wavelength range | 200 - 650 nm |
| Spectral bandwidth | 20 nm |
| Wavelength accuracy | ± 2 nm |
| Wavelength precision | ± 0.2 nm Indicates the precision performance when the power is turned ON in the single wavelength mode and the wavelength is changed. |
| Sensitivity | can be set at three levels: HIGH (x 1), MED (x 32), LOW (x 1024) |
| Wetted materials | SUS316L, PTFE (fluorocarbon polymers), quartz |
| Flow cell volume | 12 µl |
| Time constants | 11 levels can be selected, equivalent to "no filter", 0.05, 0.1, 0.5, 1.0, 1.5, 2.0, 3.0, 6.0, 8.0 and 10.0 seconds |
| Autozero | auto zero function, baseline shift function |
| Communication | |
| Gain | can be set at three levels: x 1, x 4, x 16 |
| Technical parameters | |
| Ambient conditions | operating temperature: 4 to 35 °C, relative humidity: 20 to 85 % (non condensing) |
| General | |
| Power supply | AC220-240 V, 400 VA, 50/60 Hz |
| Dimensions | 260 x 210 x 420 mm (W x H x D) |
| Weight | 16 kg |

Key features

- Pressure resistant flow cell



KNAUER offers various software control options:

www.knauer.net/softwarecontrol



For detector accessories see p. 44

Ordering details:

| Device | |
|-------------|--|
| A59200 | Fluorescence detector RF-20 A 200 - 650 nm incl. accessories and flow cell |
| A59201 | Fluorescence detector RF-20 Axs 200 - 750 nm incl. accessories and flow cell |
| A59203 | Fluorescence detector RF-20 A with photomultiplier from 200 - 900 nm incl. accessories and flow cell |
| A59204 | Fluorescence detector RF-20 Axs with photomultiplier from 200 - 900 nm incl. accessories and flow cell |
| Accessories | |
| A59210 | Xenon lamp for RF-20A/Axs fluorescence detector |
| A59211 | Flow cell for Semi micro LC cell capacity 3 µl, supports temperature control (RF-20Axs only) |
| A59212 | Inert flow cell for RF-20A/Xs, cell capacity 12 µl, contact materials: PEEK, PTFE |

Light scattering detector Sedex LC/85LT/90LT/100LT

Evaporative Light Scattering Detection (ELSD) is a universal modern technology with which every analyte that is less volatile than the mobile phase can be detected. Using the unique Low Temperature technology, this detector allows the achievement very high sensitivity. The technology is gradient compatible and is independent of the absorption characteristics of the eluents. Compounds can be universally measured with this detector (carbohydrates, proteins, peptides, polymers, lipids, steroids, etc.), regardless of their fluorescence, absorption or refractive-index characteristics. Comprehensive SOP protocols for GLP conformity and validation procedures are available.



Note: This product is only available in Germany, Austria and Switzerland. For other countries please send us a request.

Key features

- Long-life LED
- Attractively priced
- Wide application range
- Large choice of nebulizers

Specifications (for Sedex 85 LT)

Detection

| | |
|--------------------|--|
| Detector type | Light scattering detector |
| Detection channels | 1 |
| Light source | selected high efficiency blue LED (470 nm), elapsed-time counter |
| Sensitivity | < 1 ng caffeine (LOD) |
| Maximum data rate | Digital: 100 Hz/Analog: 30 Hz |

Gas requirements

| | |
|--------------------|--|
| Gas | nitrogen preferred |
| Gas flow rate | < 3 l/min |
| Gas inlet pressure | 3.5 bar |
| HPLC flow rate | standard HPLC with 4 nebulizers: 0.2 - 2.5 ml/min ultra high performance LC with 1 nebulizer |
| Maintenance | easily accessible from the front for cleaning |

Heated zone

| | |
|-------------------|-------------------|
| Temperature range | ambient to 100 °C |
|-------------------|-------------------|

Communication

| | |
|----------------------|--|
| Gain | 1 to 12 - factor 2 ¹¹ (2048) |
| Filter | moving average (0 - 10 s) |
| Analog outputs | 0 - 1 V |
| Analog control input | contact closure, TTL for ready, autozero, power down |
| Control | RS-232 |
| Power-down methods | shut-off: gas, LED, heating and/or PMT cleaning mode |

Technical parameters

| | |
|---------|----------------|
| Display | LCD and keypad |
|---------|----------------|

General

| | |
|--------------|---|
| Power supply | 230 V/50 Hz, 1.7 A - 115 V/60 Hz, 1.8 A |
| Dimensions | 250 mm x 480 mm x 550 mm (W x H x D) |
| Weight | 16 kg |

Ordering details:

Device

| | |
|---------|--|
| A0754-1 | Sensitive Light scattering detector ELSD 85LT for univ. detection 0.2 - 2.5 ml/min, 100 Hz including accessories |
| A0754-3 | High sensitive ELSD 90LT for univ. detection for HPLC and ultrafast HPLC, low temperature technology, supports high data rates |
| A0754-5 | Light scattering detector ELSD SEDEX LC for univ. detection 200 µl/min - 2 ml/min |
| A0754-6 | Ultra high sensitive light scattering detector ELSD SEDEX 100LT for univ. detection 200 µl/min - 2 ml/min 100 Hz including accessories, SAGA |

Accessories

| | |
|----------|---|
| A2618-01 | OpenLAB® CDS EZChrom Edition drivers for 80LT, 85LT, 90LT, 100LT and LC from Sedere |
| A1783-4 | Sedex Driver for Chromeleon 7.2; For Sedex 85LT / 90LT; Instrument Controller Class 3 necessary |
| A1783-5 | Sedex Driver for Chromeleon 7.2; For Sedex FP / LC / 100LT; Instrument Controller Class 3 necessary |



KNAUER offers various software control options:
www.knauer.net/softwarecontrol



For detector accessories see p. 44

Electrochemical Detector AZURA® ECD 2.1

With its measurement frequency up to 100 Hz, this electrochemical detector is specially designed for super-fast highly sensitive and selective measurement of oxidisable and reducible substances in (U)HPLC. The AZURA® ECD 2.1 comprises of a thermostat-controlled Faraday's cage, accommodating column and flow cell.

The AZURA® ECD 2.1 unites the three operating modes DC, Pulse and Scan in one instrument. The DC mode covers about 90% of all applications. The pulse mode is important for PAD (Pulsed Amperometric Detection) of e.g. carbohydrates. The scan mode is used to obtain a voltammogram in method optimization.

A digital low-pass filter provides an excellent signal-to-noise ratio. For highest sensitivity, the SenCell flow cell is recommended. The correct flow cell can be chosen from a broad variety of flow cells after our advice.

The AZURA® ECD 2.1 is controlled by software packages ClarityChrom® and Chromeleon®.



Key features

- Compact
- Thermostat included

Specifications

Puls mode

| | |
|------------------|---|
| Measuring range | 10 pA-200 µA in steps of 1, 2 and 5 |
| Filter (cut off) | Advanced Digital Filter, 0.5 -0.001 Hz in steps of 1, 2 and 5 |
| Pulse times | t1: 100-2000 ms; t2: 0-2000 ms; t3: 0-2000 ms in steps of 10 ms |
| Data recording | 20, 40, 60, 80 and 100 ms |

DC mode

| | |
|------------------|---|
| Measuring range | 10 pA-200 µA in steps of 1, 2 and 5 |
| Filter (cut off) | Advanced Digital Filter, 0.5 -0.001 Hz in steps of 1, 2 and 5 |
| Noise | < 2 pA with dummy cell |

Scan mode

| | |
|-----------------|-------------------------------------|
| Measuring range | 10 pA-200 µA in steps of 1, 2 and 5 |
| Scan speed | 1-50 mV/s in steps of 1, 2 and 5 |
| Scan cycle | half, complete, continuous |

Detection

| | |
|--------------------|---------------------------------|
| Detector type | Electrochemical detector |
| Version | DC, pulse, scan operating modes |
| Detection channels | 1 |
| Working potential | -2.0 V to +2.0 V |
| Maximum data rate | 10 Hz |
| Autozero | triggered via key, TTL, RS-232 |

Heated zone

| | |
|-------------------|--|
| Thermostating | oven included |
| Temperature range | from 7°C above ambient temperature to 45°C |

Communication

| | |
|---------|---|
| Control | Parametric control and data-acquisition via LAN port (USB service port) |
|---------|---|


Ordering details:

Device

| | |
|-------|---|
| A1651 | Electrochemical detector AZURA® ECD 2.1 without flow cell |
|-------|---|

Accessories (Flow cells)

| Art. no. | Flow cell type | Typical applications |
|----------|----------------------------------|--|
| A1652 | Flow cell SenCell GC Salt-Bridge | DC mode for Phenols, Polyphenols, Aminophenols |
| A1652-1 | Flow cell FlexCell Au HyREF | 3-Step PAD for Carbohydrates |
| A1652-2 | Flow cell SenCell GC HyREF | DC mode for Phenols, Polyphenols, Bisphenol A |
| A1652-3 | Flow cell SenCell Au HyREF | 4-Step PAD for Carbohydrates |

 KNAUER offers various software control options:
www.knauer.net/softwarecontrol

 For detector accessories see p. 44

CDD-10-AVP

The CDD-10-AVP is a highly sensitive conductivity detector applicable to ion chromatography or organic acid analysis. Low noise, low drift and wide dynamic range assure proven performance of the CDD-10-AVP detector. A special features is the VP key for validation. Flow cell 0.25 µl included.



Specifications

Detection

| | |
|---------------------------|--|
| Detector type | Conductivity detector |
| Detection channels | 1 |
| Measurement range | 0.01 - 52000 µS/cm |
| Noise | < 4 nS/cm |
| Drift | < 25 nS/cm per hour |
| Flow cell volume | 0.25 µl |
| Time constants | 0.05, 0.1, 0.5, 1.0, 1.5, 2.0, 3.0, 6.0, 8.0, 10.0 s |



KNAUER offers various software control options:
www.knauer.net/softwarecontrol



For detector accessories
 see p. 44

Communication

| | |
|----------------|-------------------------------------|
| Outputs | 10 mV recorder terminal, integrator |
|----------------|-------------------------------------|

Ordering details:

Device

| | |
|---------|---|
| A1252-1 | Conductivity detector CDD-10 Avp with flow cell 0.25 µl |
|---------|---|

Accessories

| | |
|---------|---|
| AZB00XA | AZURA® Interface Box IFU 2.1 LAN, A/D converter, 4 channels |
|---------|---|

AZURA® Conductivity Monitor CM 2.1S

The AZURA® CM 2.1S is a reliable conductivity monitor which is usually utilized in FPLC to follow salt gradients. By adding a pH sensor also pH values can be measured.

The contactless measurements of conductivity reduces the risk of carryover to a minimum and makes the device easy to maintain. The fully biocompatible flow cells cover a flow rate of up to 100 ml/min. Choose between ADG30GC and ADG30GD for a CM 2.1S with ready to measure flow cells for either 10 or 100 ml/min maximum flow.



Specifications

Flow cell

| | | |
|--------------------------------|---|--|
| Flow cell type | Contactless conductivity flow cell | |
| Biocompatible | Yes | |
| Connection of flow cell | Female 10-32" UNF or M8x1 thread (PEEK) - both included in shipment | |
| Capillary connection | 1/16" or 1/8" - both included in shipment | |
| Wetted materials | PEEK | |



KNAUER offers various software control options:
www.knauer.net/softwarecontrol

Flow cell features by device

| Device order number | ADG30GC | ADG30GD |
|-------------------------|-----------|------------|
| Flow cell volume | 30 µl | 300 µl |
| Max. flow rate | 10 ml/min | 100 ml/min |
| Maximum pressure | 160 bar | 100 bar |


| Detection | |
|-----------------------|---|
| Measurement accuracy | Conductivity: < 5% full scale end value Temperature: $\pm 1.0^{\circ}\text{C}$ pH: ± 0.5 pH (within 4 - 25°C) |
| Measurement precision | Conductivity: < 2 % of end value or ≤ 5 mS/cm of higher values (measured within 0.1 - 300 mS/cm; pH: ± 0.2 pH (within 4 - 25°C) |
| Measurement range | 0.1 - 999 mS/cm |
| pH measurement | 2 - 12 |
| Maximum data rate | 5 Hz |
| Supported electrodes | All pH electrodes with BNC connector and compatible flow cell |
| Communication | |
| Analog outputs | 2 channels (conductivity and pH value - not active if remote controlled by software) DAC 18 bit |
| Digital outputs | LAN; RS-232 |
| Technical parameters | |
| GLP | Electronic serial number |
| Display | LCD, 2 x 8 characters |
| Ambient conditions | Operating temperature: 4 - 40°C, 39,2 - 104 °F, relative humidity: below 90 %, non condensing |
| General | |
| Power supply | 100 - 240 V, 50 - 60 Hz, max. 20 W |
| Dimensions | 121 x 129 x 187 mm (W x H x D) |
| Weight | 3.2 kg |

Ordering details:




Device

| | |
|---------|--|
| ADG30GC | AZURA® CM 2.1S with flow cell - up to 10 ml/min - conductivity monitor with optional pH measurement |
| ADG30GD | AZURA® CM 2.1S with flow cell - up to 100 ml/min - conductivity monitor with optional pH measurement |

Accessories

| | |
|--|--|
| A4156 | Flow cell CM 2.1S for flow rates up to 10 ml/min |
| A4157 | Flow cell CM 2.1S for flow rates up to 100 ml/min |
|  A70091-2 | pH measuring kit for conductivity monitor CM 2.1S for flow rates up to 100 ml/min, delay volume ~ 80 μl , max. pressure 5 bar |
| A5813 | Flow splitter for CM 2.1S when used with flowrates over 100 ml/min |
| A9854-3 | Mounting bracket AZURA® L for AZURA® UVD 2.1S and AZURA® CM 2.1S (left-sided on AZURA® L) |

pH Sensors

| | |
|---|--|
|  A1943 | AZURA® pH flow cell for AZURA pH sensor |
|  A1942-1 | AZURA® pH dummy electrode for AZURA® pH sensor |
|  A1933-1 | pH electrode for AZURA® pH sensor & CM 2.1S respectively |

Fraction collector Foxy® R1 / R2

The Foxy® R1 fraction collector can be adapted to a broad spectrum of applications. Fractions can be collected into 96 well microplates, standard tube sizes, and bottles. For essentially unlimited volumes, funnel racks can direct fluids to any collection vessel or downstream process.



Specifications

Fraction collection

| | |
|----------------------------|--|
| Brand | Foxy R1 |
| Fractionation modes | drop counting, time intervals, volume intervals, level |
| Max. flow rate | 25 ml/min or 125 ml/min |
| Fraction capacity | consider list of racks in accessories below |
| Diverter valve | drop former (NC): 110 µl waste (NO): 130 µl |
| Wetted materials | valve: PEEK and perfluoroelastomer (FFKM), supplied ferrules: ETFE, supplied valve tubing: PTFE, supplied drain tubing: viny |

| | |
|---------------------------------|--|
| Fractionation control | operator: front panel control via touch screen LCDintegrated systems: direct communication via Ethernet (TCP/IP) and RS-232 serial communications |
| Maximum test tube height | 160 mm |
| RFID rack recognition | no |
| Number of racks | 1 |
| Capillary connection | 1/16" : 25 ml/min 1/8" : 125 ml/min 1/4" : 1000 ml/min |

Communication

| | |
|----------------|-------------|
| Control | LAN, RS-232 |
|----------------|-------------|

Technical parameters

| | |
|---------------------------|---------------------------|
| Conformity | CE, CSA |
| Display | touch screen LCD displays |
| Ambient conditions | 0–40 °C, 32–104 °F |

General

| | |
|---------------------|--|
| Power supply | 100–240 VAC, 50–60 Hz, max. 1 A |
| Dimensions | R1: 311 x 330 x 355 mm (W x D x H) R2 1/8": 311 x 533 x 378 mm (W x D x H) R2 1/4": 311 x 533 x 394 mm (W x D x H) |
| Weight | R1: 7.1 kg R2 1/8": 10.3 kg R2 1/4": 10.4 kg |

Ordering details:

Device

| | |
|---------|--|
| A59100 | Fraction collector Foxy® R1 for 1/16" or 1/8" tubing |
| A59102 | Fraction collector Foxy® R2 for 1/16" or 1/8" tubing |
| A591021 | Fraction collector Foxy® R2 for 1/4" tubing |



KNAUER offers various software control options:

www.knauer.net/softwarecontrol



For purification accessories see p. 53


Accessories

| | |
|---------|--|
| A59122 | Cooling option for Foxy® R1 with cooling hood, cooling plate and accessories |
| A59117 | Cooling rack for 144 tubes 1.5 ml for Foxy® R1 * |
| A59118 | Cooling rack for 72 Falcons 15 ml for Foxy® R1 * |
| A59119 | Cooling rack for 96-Well Microplates for Foxy® R1 * |
| A59105 | Rack for 100 vials 16 mm/max. 20 ml for Foxy® R1/R2 |
| A59104 | Rack for 144 vials 13 mm/max. 9 ml for Foxy® R1/R2 |
| A59111 | Rack for 2 microwell plates 96 well for Foxy® R1/R2 |
| A59114 | Rack for 2 x 9 bottles 480 ml for Foxy® R2 (not suitable for Foxy® R1, bottles too tall) |
| A59110 | Rack for 36 Falcon 50 ml for Foxy® R1/R2 |
| A59108 | Rack for 36 vials 25 mm/max. 70 ml for Foxy® R1/R2 |
| A59107 | Rack for 60 tubes 1.5 ml for Foxy® R1/R2 |
| A59106 | Rack for 72 Falcons 15 ml for Foxy® R1/R2 |
| A591091 | Rack with 26 funnels with vinyl tubing for Foxy R2, up to 1000 ml/min |
| A59109 | Rack with 36 funnels with vinyl tubing for Foxy® R1/R2 |
| A591092 | Scintillation rack for 36 vials 28 mm for Foxy® R1/R2 |
| A70055 | Thermostatting unit -20° to 40°C |
| A70050 | Thermostatting unit -40° to 200°C |

* for Foxy R1 with cooling option

Fraction collector LABOCOL Vario-4000


The LABOCOL Vario-4000 fraction collectors are characterized by their high robustness and optimal ratio of dimensions/benefit. The user is not limited to given rack types. The rack layout can be designed according to individual needs. Freely define the number of fraction vessels and their position. The wide application area make the Vario-4000 series ideal for use in research and development as well as in production. The Vario-4000 models differ in the base area and the flow rate range.




Specifications

Fraction collection

| | | | |
|-----------------------------|---|--|--|
| Brand | LABOCOL Vario-4000 | | |
| Max. flow rate | 100 ml/min for 1/16"; 500 ml/min for 1/8" | | |
| Fraction capacity | consider list of racks in accessories below | | |
| Wetted materials | Stainless steel, PEEK and PTFE | | |
| Number of racks | 3 (Vario-4000) / 5 (Vario-4000 Plus) | | |
| Capillary connection | 1/16" : 100 ml/min 1/8" : 500 ml/min 1/4" : 1000 ml/min | | |

 KNAUER offers various software control options:
www.knauer.net/softwarecontrol

 For purification accessories see p. 53

Communication

| | | | |
|----------------|-------------|--|--|
| Control | LAN, RS-232 | | |
|----------------|-------------|--|--|

Technical parameters

| | | | |
|---------------------------|--------------------|--|--|
| Ambient conditions | 0-40 °C, 32-104 °F | | |
|---------------------------|--------------------|--|--|

General

| | | | | |
|---------------------|---|------------------|------------------|------------------|
| Power supply | 100-240 VAC, 50-60 Hz, max. 2.5 A | | | |
| Dimensions | Vario-4000 | 30 x 50 cm (WxD) | Max. floor space | |
| | Vario-4000 Plus | 46 x 50 cm (WxD) | Vario-4000 | 24 x 41 cm (WxD) |
| | | min. H *: 52 cm | Vario-4000 Plus | 40 x 41 cm (WxD) |
| | | max. H *: 67 cm | | |
| Weight | 8 kg (Vario-4000) / 10 kg (Vario-4000 Plus) | | | |



* with touchpanel

Ordering details:

Device

| | |
|---------|--|
| A591022 | Fraction collector LABOCOL Vario-4000, for 1/16" or 1/8" tubing |
| A591024 | Fraction collector LABOCOL Vario-4000, for 1/4" tubing |
| A591023 | Fraction collector LABOCOL Vario-4000 Plus, for 1/16" or 1/8" tubing |
| A591026 | Fraction collector LABOCOL Vario-4000 Plus, for 1/4" tubing |

Accessories

| | |
|--|---|
| A591029 | Touchpanel for LABOCOL Vario-4000/Vario-4000 Plus |
| A59130 | Rack standard for 80 tubes 18 mm/max. 36 ml/ 15 ml Falcons for LABOCOL Vario-4000/Vario-4000 Plus |
| A59131 | Rack micro for 125 tubes 10.5 mm/max. 9 ml for LABOCOL Vario-4000/Vario-4000 Plus |
| A59132 | Rack prep for 20 tubes 36 mm/max. 140 or 240 ml for LABOCOL Vario-4000/Vario-4000 Plus |
| A59133 | Rack semiprep for 39 tubes 26 mm/max 80 ml for LABOCOL Vario-4000/Vario-4000 Plus |
| A59134 | Rack for 24 Falcon® tubes of 50 ml for LABOCOL Vario-4000/Vario-4000 Plus |
|  A20521 | Micro test tubes, 9 ml, 100 pcs, L 150 mm, OD 10,5 mm for rack A59131 |
|  A20522 | Preparative tubes, 25 pcs, L 284 mm x OD 36 mm, V 240 ml for rack A59132 |

AZURA® Degasser DG 2.1S

Dissolved gases in the solvent can cause bubbles in the flow path of pumps and detectors. Reliable chromatographic separation therefore requires degassing of the solvent. The small analytical 2-channel degasser DG 2.1S is equipped with two degassing chambers and can thus degas two solvents simultaneously.



KNAUER offers various software control options:

www.knauer.net/softwarecontrol

Specifications

Degasser module

| | |
|-------------------------------|--|
| Degasser channels | 2 |
| Max. flow rate/channel | 10 ml/min |
| Recommended flow rate/channel | 2.8 ml/min |
| Degassing method | Gas permeation through a fluoropolymer membrane |
| Degassing chamber volume | 285 µl |
| Solvent applicability | universal, except hydrochloric acid, halogenated hydrocarbons, hexafluoro isopropanol (HFIP) |
| Wetted materials | PTFE, PPS, PEEK, Systec AF™ |
| Pressure decline | 1.37 mm (Hg/ml/min) |
| Maximum pressure stability | 70 psi |

Technical parameters

| | |
|--------------------|---|
| Display | 1 LED |
| Ambient conditions | temperature range: 4 - 40 °C, 39.2 - 104 °F air humidity: below 90 %, non-condensing |

General

| | |
|--------------|--------------------------------|
| Power supply | 85 - 265 V, 50 - 60 Hz, 20 W |
| Dimensions | 121 x 138 x 190 mm (W x H x D) |
| Weight | 2.3 kg |
| Connector | 1/4"-28 UNF female port |

Feature overview

| Order no. | Degasser type | Channels | Max. flow rate | Chamber volume |
|-----------|---------------------|----------|------------------------------------|--------------------|
| AZE02 | analytical | 2 | 10 ml/min (recommended 2 ml/min) | 285 µl per channel |
| AZE03-1 | analytical | 4 | 10 ml/min (recommended 2 ml/min) | 285 µl per channel |
| A5335 | analytical, for GPC | 2 | 10 ml/min (recommended 3 ml/min) | 480 µl per channel |
| A5328 | semi-preparative | 2 | 30 ml/min (recommended 15 ml/min) | 5.3 µl per channel |
| AZE02-1 | preparative | 2 | 200 ml/min (recommended 75 ml/min) | 23 ml per channel |
| AZE03 | preparative | 4 | 200 ml/min (recommended 75 ml/min) | 23 ml per channel |

Ordering details:

Device

| | |
|---------|---|
| AZE02 | Biocompatible 2 channel degasser |
| AZE03-1 | Analytical 4 channel degasser, biocompatible |
| A5335 | Analytical 2 channel GPC degasser |
| A5328 | Semi-preparative 2 channel degasser |
| AZE02-1 | Preparative 2 channel degasser, biocompatible |
| AZE03 | Preparative 4 channel degasser, biocompatible |

AZURA® Valve unifier VU 4.1

The valve drive AZURA® Valve Unifier VU 4.1 enables automatic valve switching. Due to its low switching time, the flow path is blocked only for a very short time, and pressure peaks are reduced to a minimum. Valves are identified via RFID technology, which guarantees an easy valve exchange of KNAUER valves. An additional feature is the easy monitoring of GLP data, which simplifies maintenance such as the exchange of a rotor seal. The display enables user-friendly standalone operation. In addition, the valve drive can be operated with software as well with an optional touch display (Mobile Control), via LAN or analog input/output, by which it can be integrated into nearly every LC system.



Specifications

Communication

| | |
|-------------------|---|
| Interfaces | LAN, display, terminal strip |
| Control | Display, software, event control |
| Inputs | Binary control; Home, Backward/Inject, Forward/Load, Start (IN) |
| Outputs | Trigger out, Event |

General

| | |
|---------------------------|---|
| Power supply | External DC 24V, 65 W |
| Dimensions | 80 x 123 x 192 mm (W x H x D) |
| Weight | 1.9 kg |
| Ambient conditions | Temperature range: 4-40 °C; 39.2-104 °F below 90 % humidity (non condensing) |

Key features

- One valve drive for all valves
- Ultra fast switching cycle
- Easy maintenance
- Compact
- Multiple interfaces and drivers available




KNAUER offers various software control options:
www.knauer.net/softwarecontrol



For valve accessories see p. 48



Valve drive VU 4.1 (AWA01XA) with 6 port 2-position valve (AVC28AC)



AVS34CE

AVN94CE

AVR38AC

AVC38AC

AVS62CE

Ordering details:

Device


AWA01XA VU 4.1 valve drive for V 4.1 valves

Accessories

A9854-3 Mounting bracket AZURA® L for AZURA® Valve Unifier VU 4.1 (both-sided) or AZURA® UVD 2.1S and AZURA® CM 2.1S (left-sided on AZURA® L)




Valves for valve drive VU 4.1

Manual valves*

| Ports | Stator material | Rotor material | Max. pressure [bar] | Bore size [mm] | Capillary connection | Order number |
|-------|----------------------|----------------|---------------------|----------------|--------------------------|---|
| 6 | SST DLC ¹ | POM | 100 | 0.75 | 1/16 (UNF 10-32) | AVJ23AF  |
| 6 | SST DLC ¹ | PEEK | 500 | 0.75 | 1/16 (UNF 10-32) | AVJ26AE |
| 6 | SST DLC ¹ | Vespel | 1200 | 0.3 | 1/16 (UNF 10-32) | AVI28AC |
| 6 | PEEK | PEEK | 240 | 0.75 | 1/16 (UNF 10-32) | AVG24CE |
| 6 | SST DLC ¹ | PEEK | 300 | 1.5 | 1/8" (UNF 1/4-28, coned) | AVK25AE |
| 6 | PEEK | PEEK | 100 | 2.0 | 1/8" (UNF 1/4-28, coned) | AVL23CE |
| 8 | SST DLC | PEEK | 500 | 0.75 | 1/16 (UNF 10-32) | AVJ36AE |
| 8 | SST DLC ¹ | Vespel | 1200 | 0.3 | 1/16 (UNF 10-32) | AVI38AC |






* The mounting bracket A9853 is required to mount the manual valves to an AZURA® L device.

2-position valves

| Ports | Stator material | Rotor material | Max. pressure [bar] | Bore size [mm] | Capillary connection | Order number |
|-------|----------------------|----------------|---------------------|----------------|--------------------------|---|
| 6 | SST DLC ¹ | POM | 100 | 0.75 | 1/16 (UNF 10-32) | AVD23AF  |
| 6 | SST DLC | PEEK | 500 | 0.75 | 1/16 (UNF 10-32) | AVD26AE |
| 6 | SST DLC ¹ | Vespel | 1200 | 0.3 | 1/16 (UNF 10-32) | AVC28AC |
| 6 | PEEK | PEEK | 240 | 0.75 | 1/16 (UNF 10-32) | AVD24CE |
| 6 | SST DLC ¹ | PEEK | 300 | 1.5 | 1/8" (UNF 1/4-28, coned) | AVE25AE |
| 6 | SST DLC ¹ | PEEK | 300 | 1.5 | 1/8" (UNF 1/4-28, coned) | AVE25AI*  |
| 6 | PEEK | PEEK | 100 | 2.0 | 1/8" (UNF 1/4-28, coned) | AVF23CE |
| 8 | SST DLC ¹ | PEEK | 500 | 0.75 | 1/16 (UNF 10-32) | AVD36AE |
| 8 | SST DLC ¹ | Vespel | 1200 | 0.3 | 1/16 (UNF 10-32) | AVC38AC |
| 8 | PEEK | PEEK | 50 | 2.0 | 1/8" (UNF 1/4-28, coned) | AVF32CE  |

* Special version of AVE25AE with 2-channel rotor seal instead of 3 channels.

Multiposition valves

| Ports | Stator material | Rotor material | Max. pressure [bar] | Bore size [mm] | Capillary connection | Order number |
|-------|----------------------|----------------|---------------------|----------------|--------------------------------|---|
| 6 | SST DLC ¹ | POM | 100 | 0.75 | 1/16 (UNF 10-32) | AVS23AF  |
| 6 | SST DLC ¹ | PEEK | 500 | 0.75 | 1/16 (UNF 10-32) | AVS26AE |
| 6 | SST DLC ¹ | Vespel | 1200 | 0.3 | 1/16 (UNF 10-32) | AVR28AC |
| 6 | SST DLC ¹ | PEEK | 300 | 1.5 | 1/8" (UNF 1/4-28 coned) | AVT25AE |
| 8 | SST DLC ¹ | PEEK | 300 | 0.75 | 1/16 (UNF 10-32) | AVS35AE |
| 8 | SST DLC ¹ | PEEK | 500 | 0.75 | 1/16 (UNF 10-32) | AVS36AE |
| 8 | SST DLC ¹ | Vespel | 1200 | 0.3 | 1/16 (UNF 10-32) | AVR38AC |
| 8 | PEEK | PEEK | 240 | 0.75 | 1/16 (UNF 10-32) | AVS34CE |
| 8 | SST DLC ¹ | PEEK | 200 | 1.5 | 1/8" (UNF 1/4-28 coned) | AVT34AE |
| 8 | PEEK | PEEK | 50 | 2.0 | 1/8" (UNF 1/4-28, flat-bottom) | AVU32GE  |
| 8 | PEEK | PEEK | 50 | 2.0 | 1/8" (UNF 1/4-28 coned) | AVU32CE |
| 12 | SST DLC ¹ | PEEK | 100 | 1.5 | 1/8" (UNF 1/4-28 coned) | AVT53AE |
| 12 | PEEK | PEEK | 100 | 1.5 | 1/8" (UNF 1/4-28 coned) | AVT53CE  |
| 16 | SST DLC ¹ | POM | 100 | 0.75 | 1/16 (UNF 10-32) | AVQ63AF  |
| 16 | SST DLC ¹ | PEEK | 500 | 0.6 | 1/16 (UNF 10-32) | AVQ66AE |
| 16 | PEEK | PEEK | 50 | 0.75 | 1/16 (UNF 10-32) | AVS62CE |
| 16 | PEEK | PEEK | 150 | 0.75 | 1/16 (UNF 10-32) | AVS63CE  |

¹ stainless steel coated with diamond-like carbon

Special purpose valves*

| Valves | Capillary connection | Max. pressure [bar] | Bore size [mm] | Order number |
|--|----------------------|---------------------|----------------|--------------|
| Column selection valve, biocompatible. Allows switching of up to 5 columns incl. bypass and reverse flow option. | 1/16" (UNF 10-32) | 50 | 0.4 | AVZ52CE |
| Multi-injection valve, biocompatible. Allows manual and automated sample loading as well as direct injection. | 1/16" (UNF 10-32) | 240 | 0.75 | AVN94CE |
| Multi-injection valve, stainless steel. Allows manual and automated sample loading as well as direct injection. | 1/16" (UNF 10-32) | 500 | 0.75 | AVN96AE |

*for detailed information please check our website: www.knauer.net/valves

K-7400S Semi-Micro Osmometer

KNAUER is one of the pioneers in the field of osmometry and known for its reliable and user friendly instruments for many decades.

Our freezing point osmometer K-7400S allows the easy and fast determination of the osmolality of various aqueous solutions. Also, the freezing point depression of the samples can be read. The proven technology of freezing point determination in combination with the robust and intelligent design of the device allows reproducible measurements.

The instrument is equipped with a peltier cooler and an integrated microprocessor controlling the automated measurement. The freezing point osmometer is a standalone device that optionally can be equipped with a printer. Furthermore, the device can be controlled via the EuroOsmo 7400 software. The software automatically plots the temperature curve for each measurement and calibration and allows saving of the measured values. In addition, the data can optionally be exported into various file formats for archival storage.

The measurement specifications of the KNAUER Semi-Micro Osmometer K-7400S comply with the European Pharmacopoeia for osmolality (Ph. Eur. 2.2.35, 10/2021) in the pharmaceutically relevant range of 0-400 mOsmol/kg.



Key features

- Made in Germany
- 55 years experience
- Fast measurements

Specifications

Technical parameters

| | |
|-------------------------|---|
| Sample volume | 50 - 150 µl |
| Osmolality range | 0 - 2000 mOsmol/kg |
| Resolution | Osmolality: integer value without decimal part, e.g. 850 mOsmol Temperature: value with three digits, e.g. -1.576 °C |
| Test time | ~ 2 min |
| Precision | SD ≤ 4 mOsmol/kg [0 – 400 mOsmol/kg] RSD ≤ 1 % [400 – 2000 mOsmol/kg] |
| Linearity | ± 1 % [0 - 1500 mOsmol/kg] ± 1.5 % [0 - 2000 mOsmol/kg] |
| Calibration | Two-point calibration (0 mOsmol/kg and one freely selectable osmolality) optional: Three-point calibration (0 Osmol/kg and two freely selectable osmolalities) |

General

| | |
|---------------------------|--|
| Power supply | 100 - 240 V, 50 - 60 Hz, 70 W |
| Dimensions | 160 x 182 x 340 mm (W x H x D) |
| Weight | 5.3 kg |
| Ambient conditions | Temperature range: 10-35 °C Rel. humidity: 20-80 % (non-condensing) |

Communication

| | |
|-------------------|--|
| Interfaces | RS-232 port |
| Control | Keypad (LED display, 2 rows with 24 characters) optional: EuroOsmo7400 Software |

Ordering details:

Device

A0006AC Osmometer for the determination of osmolality or freezing point of aqueous solutions

Accessories

| | |
|---------|--|
| A0840-2 | Measuring head for plastic sample tubes; compatible with the K-7400 and the K-7400S Semi-Micro Osmometer |
| A3705 | EuroOsmo 7400 - software for K-7400 and K-7400S osmometers |
| A3711 | Plain paper printer for freezing point osmometer K-7400 and K-7400S |
| A13270 | Barcode Scanner with USB cable, for EuroOsmo 7400 |
| A0011XB | Upgrade Kit for K-7400 required to use measuring head for sample tubes (A0840-2) |
| A0272 | 500 Pack of plastic sample tubes for Semi-Micro Osmometer K-7400S |

For more osmometry accessories and standards see p. 52

Pump accessories



AZC00



A5325



A2056



A5324

Eluent trays & bottles

Recommended for AZURA systems: Set of 4 eluent bottles 1000 ml, 1 bottle 250 ml for piston backflushing, incl. closed caps, GL45 use with AZURA tubing kit

A5324-3

Set of 4 eluent bottles 1000 ml, 1 bottle 250 ml for piston backflushing, incl. caps for eluent tubing, GL45

A5324-2

Set of 4 eluent bottles 1000 ml, incl. caps for eluent tubing, GL45

A5324

Set of eluent supply bottles, 3 x 2.5 l brown glass bottles (borosilicate glass) with special round bottom for minimal eluent remainder, for preparative HPLC/FPLC, includes screw-type cap

A70037

Eluent supply bottle plastic 2 l incl. cap and tubing for IC and ECD systems

A70038

AZURA® Eluent tray E 2.1L for AZURA® devices with a capacity of 6 x 1 l bottles or 4 x 2.5 l bottles or 2 x 5 l bottles, (delivery without bottles)

AZC00

Eluent bottle 1000 ml, clear glass, incl. cap for eluent tubing, GL45

A5325

Eluent supply bottle 2000 ml, GL45 thread, round, clear glass, without screw cap

A59158-1

250 ml bottle for piston backflushing

A2056

Waste can 2.5 l with GL45 screw top, UN-approved, 153 x 115 x 202 mm

A59173

Waste can, 10 l with GL45 screw top, UN-approved, 192x317x231 mm

A59256



A5390



A5398



A5396

Mass flow controllers*

Mini CORI-Flow (M13) Mass flow controller incl. mounting block, Flow: 1 - 50 ml/min, stainless steel 316

A5390

Mini CORI-Flow (M14) Mass flow controller incl. mounting block, Flow: 2 - 833 ml/min, stainless steel 316

A5391

Mini CORI-Flow (M14) Mass flow controller incl. mounting block, Flow: 2 - 833 ml/min, stainless steel 316, Profibus

A5391P

Mini CORI-Flow (M13) Mass flow controller incl. mounting block, Flow: 1 - 50 ml/min, stainless steel 316, Profibus

A5393

Mini CORI-Flow (M12) Mass flow controller incl. mounting block, Flow: 0,03 - 3,3 ml/min, stainless steel 316

A5394

Mini CORI-Flow (M13) Mass flow controller incl. mounting block, Flow: 1 - 50 ml/min, Hastelloy-C22

A5395


Mini CORI-Flow (MI140) Mass flow controller, Flow: 2 - 833 ml/min, stainless steel 316, Profinet

A5398


HI-TEC Bright display for Mini CORI-Flow mass flow controller (display, setpoint and counter)

A5396

*analog and bus versions on request



AZZ00MB



AZZ00MC



AZZ10ME



A5830

Static mixers

| | |
|--|---------|
| AZURA® HPLC mixer up to 1000 bar, 50 µl mixing volume, stainless steel | AZZ00MB |
| AZURA® HPLC mixer up to 1000 bar, 100 µl mixing volume, stainless steel | AZZ00MC |
| AZURA® HPLC mixer up to 1000 bar, 200 µl mixing volume, stainless steel | AZZ00MD |
| AZURA® HPLC mixer up to 40 MPa, 250 µl mixing volume, PEEK (biocompatible) | AZZ10ME |
| HyperShear Static Mixer, 1.5 ml, 1-40 ml/min, max. 414 bar, stainless steel and PEEK, incl. mounting brackets for AZURA® L devices (A9853-8) | A5830 |
| Mounting bracket AZURA® L for Hypershear mixing chambers | A9853-8 |




A0285

Dynamic mixers

| | |
|--|---------|
| Dynamic mixing chamber (250 V), stainless steel, analytical, 1/16", up to 420 bar, 1740 µl mixing volume | A0285 |
| Dynamic mixing chamber (115 V), stainless steel, analytical, 1/16", up to 420 bar, 1740 µl mixing volume | A02851 |
| Dynamic mixing chamber (250 V), stainless steel, preparative, 1/8", up to 250 bar, 5.9 ml mixing volume | A0581 |
| Dynamic mixing chamber (115 V), stainless steel, preparative, 1/8", up to 250 bar, 5.9 ml mixing volume | A05811 |
| Dynamic mixing chamber (250 V), titanium, analytical, 1/16", up to 420 bar, 1740 µl mixing volume | A0275 |
| Dynamic mixing chamber (115 V), titanium, analytical, 1/16", up to 420 bar, 1740 µl mixing volume | A02751 |
| Dynamic mixing chamber (250 V), titanium, preparative, 1/8", up to 250 bar, 5.9 ml mixing volume | A70581 |
| Dynamic mixing chamber (115 V), titanium, preparative, 1/8", up to 250 bar, 5.9 ml mixing volume | A705811 |
| Extension unit for dynamic mixer A70581/A705811 | A2515 |



A3373



A3374



A3375



A3364

Solvent filters & inlet tubing

Mobile Phase Filter, stainless steel, 2 μm , 1/8" pipe OD, suitable for all analytical HPLC systems, max. flow rate 50 ml/min

A3373

Mobile Phase Filter, stainless steel, 20 μm , for 1/8" OD, compatible with the AZURA® Tubing Kit (A9650), suitable for all analytical and semi preparative HPLC systems, max. flow rate 100 ml/min

A3374

Mobile Phase Filter, stainless steel, 10 μm , for 1/8" OD, compatible with the AZURA® Tubing Kit (A9650), suitable for all analytical HPLC systems, max. flow rate 50 ml/min

A3375

Mobile Phase Filter, Biocompatible PE, 20 μm , 1/8" pipe OD, suitable for all FPLC systems, max. flow rate 500 ml/min

A3364

AZURA® Tubing kit with cap and solvent filter (A3375, stainless steel, 10 μm), suitable for all analytical HPLC systems

A9650

AZURA® Tubing kit bio with cap and insert, solvent filter inlet and fittings, 1set

A96507

Inlet-bushing kit with 1/4"-PTFE Tubing and 20 μm stainless steel solvent filter (up to 250 ml/min)

A58207



AZZ00NA



AZZ00NB



AZZ10NB



FZZ2

Pulse dampers

KNAUER Pulse Damper, Low Volume, 275 μl , stainless steel, 1/16", 1000 bar

AZZ00NA

KNAUER Pulse Damper, High Volume, 290 μl , stainless steel, 1/16", 1000 bar

AZZ00NB

KNAUER Pulse Damper, 290 μl , PEEK, 1/16", 225 bar

AZZ10NB

Mounting Bracket KNAUER Pulse Damper

FZZ2



A9861



A9868



A98611



A58267

Pump head inlet fittings

| | |
|--|---------|
| Pump head inlet for AZURA® Pump P 2.1L, BlueShadow 80P, 1/4" (NPT), stainless steel | A9861 |
| Pump head inlet for AZURA® Pump P 2.1L, Set, 1/2"-20 UNF, PEEK with CTFE (Kel-F) adapter, including tubing 1/4" PTFE | A9868 |
| Inlet bushing for prep pump heads, adapter to 3/8" tube stub | A98611 |
| Inlet bushing for binary LPG prep pump heads, LPG inlet to 3/8" tube stub | A98612 |
| Inlet bushing for LPG prep pump heads, LPG ternary inlet to 3/8" tube stub | A98613 |
| Male connector to connect a 1/4" OD capillary to a 1/4" male NPT adapter (material: stainless steel) for A9861 | A58267 |
| Adapter female 1/4-28 Flat Bottom to 1/2-20 UNF, for pump head inlet A9868 & 1/4" VICI valves, material: PEEK | A142605 |
| Male connector to connect a 4 mm OD capillary to a 1/4" male NPT adapter (material: stainless steel) for A9861 | A58268 |
| Male connector to connect a 1/8" OD capillary to a 1/4" male NPT adapter (material: stainless steel) for A9861 | A58269 |
| Inlet-bushing kit for P 2.1S, P 4.1S, P6.1L, 40P and S1050 pumps for pump heads 10 ml (1/8" capillaries) | A58202 |
| Inlet-bushing kit for P 2.1S, P 4.1S, P6.1L, 40P and S1050 pumps for pump heads 10 ml (1/16" capillaries) | A58203 |
| Inlet-bushing kit for P 2.1S, P 4.1S, P6.1L, 40P and S1050 pumps for pump heads 50 ml (1/8" capillaries) | A58204 |
| Inlet-bushing kit for P 2.1S, P 4.1S, P6.1L, 40P and S1050 pumps for pump heads 50 ml (1/16" capillaries) | A58205 |
| Inlet-bushing kit with 1/4"-PTFE Tubing and 20 µm stainless steel solvent filter (up to 250 ml/min) | A58207 |



A5822



A7200

Pump head outlet fittings

| | |
|---|-------|
| Outlet-bushing kit 1/8" tube stub for S1800, 80P and P 2.1L pumps | A5822 |
| Adapter to connect a capillary with 1/16" OD (thread: 10-32 UNF) to AZURA® Pump P 2.1L or BlueShadow Pump 80P outlet (1/8", M8x1 thread), material: stainless steel, 2 pcs. | A7200 |



AHB40XA



AHB32



AHC20



AHB40CA

Replacement pump heads for analytical AZURA® pumps

| | |
|---|---------|
| Pump head 10 ml, stainless steel, 862 bar | AHB40XA |
| Pump head 10 ml, stainless steel, 400-700 bar | AHB40 |
| Pump head 10 ml, ceramic with PEEK bushings, 400 bar | AHB32 |
| Pump head 10 ml, ceramic with Ti-bushings, 400 bar | AHB32DA |
| Pump head 10 ml, Hastelloy-C, 400 bar, for corrosive chemicals | AHB43 |
| Pump head 50 ml, stainless steel, 300 bar | AHC20 |
| Pump head 50 ml, ceramic, 200 bar | AHC22 |
| Pump head 50 ml, Hastelloy-C, 300 bar, for corrosive chemicals | AHC23 |
| Pump head 10 ml, stainless steel, 700 bar, for aqueous solutions | AHB40FA |
| Pump head 10 ml, ceramic with Titanium bushings, 400 bar, for aqueous solutions | AHB32GA |
| Pump head 50 ml, stainless steel, 300 bar, for aqueous solutions | AHC20FA |
| Pump head 50 ml, ceramic, 200 bar, for aqueous solutions | AHC22FA |
| Pump head 5 ml, stainless steel, 1000 bar | AHA60 |
| Pump head 10 ml, stainless steel, for normal phase applications | AHB40BA |
| Pump head 10 ml, stainless steel, 700 bar, for high-temperature applications | AHB40CA |
| Pump head 50 ml, stainless steel, for normal phase applications | AHC20BA |
| Pump head 50 ml, stainless steel, 300 bar, for high-temperature applications | AHC20CA |



A4029-1



A4029V2



A4021-1



A4021V2

Replacement pump heads for preparative AZURA® pumps

| | |
|---|---------|
| Pump head 100 ml, stainless steel, 400 bar | A4029-1 |
| Pump head 100 ml, titanium, 400 bar | A4029V2 |
| Pump head 250 ml, stainless steel 200 bar | A4021-1 |
| Pump head 250 ml, titanium, 200 bar | A4021V2 |
| Pump head 500 ml, stainless steel , 100 bar | A4038-1 |
| Pump head 500 ml, titanium, 100 bar | A4038V2 |
| Pump head 1000 ml, stainless steel, 50 bar | A4022-1 |
| Pump head 1000 ml, titanium, 50 bar | A4022V2 |



A06840



A06841



A068411



A1122

Check valves

| | |
|--|---------|
| Check valve unit for 10/50 ml pump heads, for dosing applications, Bore: Ø1.4 mm, Ball: Ø1.75 mm | A06840 |
| Check valve unit for 10 ml pump heads, for HPLC applications, Bore: Ø0.7 mm, Ball: Ø1.75 mm | A06841 |
| Spring-loaded check valve unit for 10 ml/50 ml pump heads, for normal phase applications, Bore: Ø1.4 mm, Ball: Ø1.75mm | A068411 |
| Check valve unit (KEL-F) for 10 ml pump head, for aggressive substances, Bore: Ø0.7 mm, Ball: Ø1.75mm | A068412 |
| Check valve unit for 50 ml pump heads, for HPLC applications, Bore: Ø1.2 mm, Ball: Ø1.75 mm | A06842 |
| Check valve unit (KEL-F) for 50 ml pump head, for aggressive substances, Bore: Ø1.2 mm, Ball: Ø1.75 mm | A068422 |
| Check valve unit stainless steel/PEEK for 500 ml and 1000 ml pump heads | A1080 |
| Check valve unit titanium/PEEK for 500 ml and 1000 ml pump heads | A1080V1 |
| Check valve unit titanium/KEL-F for 500 ml and 1000 ml pump heads | A1080V2 |
| Check valve unit stainless steel/KEL-F for 500 ml and 1000 ml pump heads | A1080V3 |
| Check valve unit stainless steel/PEEK for 100 ml and 250 ml pump heads | A1122 |
| Check valve unit titanium/PEEK for 100 ml and 250 ml pump heads | A1122-1 |
| Check valve unit titanium/KEL-F for 100 ml and 250 ml pump heads | A1122-2 |
| Check valve unit stainless steel/KEL-F for 100 ml and 250 ml pump heads | A1122-3 |



AZZ00AA




AZZ00AB




AZZ10AB

LPG modules

| | |
|---|---------|
| LPG module for Pump P 2.1L binary up to 800 ml/min (stainless steel) | AZZ00AA |
| LPG module for Pump P 2.1L ternary up to 220 ml/min (stainless steel) | AZZ00AB |
| LPG module for Pump P 2.1L ternary up to 220 ml/min (PEEK) | AZZ10AB |



A2034-1



A2035-1



A57024



A57036-1

Temperature control

| | |
|--|----------|
| Pump head cooling and heating device for 100/250/500/1000 ml/min pump heads | A2034-1 |
| Pump head cooling and heating device for 10 and 50 ml/min pump heads | A2035-1 |
| Temperature controller for column heating sleeve | A57024 |
| Heating solution for 10 and 50 ml/min pump heads, incl. heating plate and insulation sleeve | A57036-1 |
| Heating solution for 10 and 50 ml/min pump heads, incl. heating plate and insulation sleeve (without temperature controller) | A57037-1 |



A96423



A96424



A96425




A9670

Maintenance kits for AZURA® Pumps

| | |
|---|---------|
| Maintenance kit for AZURA® Pump P 2.1S / P 4.1S, AZURA® Pump P 6.1L, BlueShadow 40P, 10 ml pump head, including 1 set of gaskets, 2 piston rods, 2 sapphire backing rings, 2 O-rings, for organic solvents | A96423 |
| Maintenance kit for AZURA® Pump P 4.1S & P 2.1S, AZURA® Pump P 6.1L, BlueShadow 40P, 50 ml pump head, including 1 set of gaskets, 2 piston rods, 2 sapphire backing rings, 2 O-rings | A96424 |
| Maintenance kit for AZURA® Pump P 4.1S & P 2.1S, AZURA® Pump P 6.1L, BlueShadow 40P, 50 ml pump head, including 1 set of gaskets, 2 piston rods, 2 sapphire backing rings, 2 O-rings, for aqueous solutions | A964242 |
| Maintenance kit 100 ml for AZURA® Pump P 2.1L and BlueShadow 80P, including 2 sets of gaskets, 2 piston rods, 2 springs | A96425 |
| Maintenance kit 250 ml for AZURA® Pump P 2.1L and BlueShadow 80P, including 2 sets of gaskets, 2 piston rods, 2 springs | A96426 |
| Maintenance kit 500 ml for AZURA® Pump P 2.1L and BlueShadow 80P, including 2 sets of gaskets, 2 piston rods, 2 springs | A96427 |
| Maintenance kit 1000 ml for AZURA® Pump P 2.1L and BlueShadow 80P, including 2 sets of gaskets, 2 piston rods, 2 springs | A96428 |
| Rebuild-Kit for one pump head AZURA® Pump P 2.1L and BlueShadow 80P (100 ml/250 ml), Venting screw KEL-F, Check valve unit KEL-F, O-ring | A58211 |
| Rebuild-Kit for one pump head AZURA® Pump P 2.1L and BlueShadow 80P (500ml/1000 ml), Venting screw KEL-F, Check valve unit KEL-F, O-ring | A58212 |
| Rebuild-Kit Kel-F for AZURA® Pump P 4.1S/P 2.1S/P 6.1L/40P, 10 ml/min pump head | A5821-1 |
| Rebuild-Kit Kel-F for AZURA® Pump P 4.1S/P 2.1S/P 6.1L/40P, 50 ml/min pump head | A5821-2 |
| Rebuild-Kit for aqueous eluents (for P2.1S/P4.1S/P 6.1L/40P with 10 ml pump head) | A5823 |
| Rebuilt-Kit for aqueous eluents (for AZURA® P 4.1S, P 2.1S, P 6.1L and BlueShadow 40P (50 ml pump head)) | A5823-1 |
| Maintenance tool kit for 10 ml pump heads | A9670 |
| Maintenance tool kit for 50 ml pump heads | A9671 |
| Maintenance tool kit for 100 ml pump heads | A9672 |
| Maintenance tool kit for 250 ml pump heads | A9673 |
| Maintenance tool kit for 500 ml pump heads | A9674 |
| Maintenance tool kit 1000 ml pump heads | A9675 |
| Filter Cartridge for pump P 6.1L/40P, Titanium frit, 2 µm pore size, 50 ml/min maximum flow, High capacity filter, 60 µl volume, 3 pcs. | A9661 |
| Filter Cartridge for pump P 6.1L/40P, Stainless steel frit, 2 µm pore size, 10 ml/min maximum flow, Volume optimized filter, 20 µl volume, 3 pcs. | A96601 |
| Replacement filters for pump P 6.1L/40P ceramic for serial no 2109 and newer, PEEK frit, 2 µm pore size, 3 pcs. | A9663 |

Autosampler accessories - spare parts



A500526

Buffer tubings kits

Buffer tubing for AS 6.1L, 500 µl; Spare part for AAA00AA, AAA01AA, AAA10AA, AAA11AA, AAA50AA & AAA51AA A500525

Buffer tubing for AS 6.1L, 1000 µl; Spare part for AAA20AA & AAA21AA A500526

Buffer tubing for AS 6.1L, 2000 µl; Spare part for AAA31AA, AAA40AA & AAA41AA A500527



A500519

Rotor seals

Rotor seal for AS 6.1L, 700 bar, Vespel; Spare part for AAA00AA & AAA01AA A500519

Rotor seal for AS 6.1L, 1000 bar, Vespel; Spare part for AAA50AA & AAA51AA A500520

Rotor seal for AS 6.1L, 1240 bar, Vespel; Spare part for AAA10AA & AAA11AA A500521

Rotor seal for AS 6.1L, 345 bar, PEEK; Spare part for AAA20AA & AAA21AA A500522

Rotor seal for AS 6.1L, 200 bar, ValconH; Spare part for AAA40AA & AAA41AA A500523

Rotor seal for AS 6.1L, 200 bar, ValconE; Spare part for AAA31AA A500524



A500515




A500512

Syringes

100 µl Syringe for AS 6.1L; alternative to standard configuration A500515

250 µl Syringe for AS 6.1L; Spare part for AAA00AA, AAA01AA, AAA10AA, AAA11AA, AAA50AA & AAA51AA A500512

500 µl Syringe for AS 6.1L; Spare part for AAA20AA & AAA21AA A500513

 1000 µl Syringe for AS 6.1L; alternative to standard configuration A500539

2500 µl Syringe for AS 6.1L; Spare part for AAA31AA, AAA40AA & AAA41AA A500514



A50078



A50077

Sample loops

| | |
|---|---------|
| Sample loop incl. fittings, 10 µl, stainless steel; Spare part for AAA10AA & AAA11AA | A50078 |
| Sample loop incl. fittings, 10 ml, stainless steel; Spare part for AAA40AA & AAA41AA | A500509 |
| Sample loop incl. fittings, 10 ml, PEEK; Spare part for AAA31AA | A500511 |
| Sample loop incl. fittings, 100 µl, stainless steel; Spare part for AAA00AA, AAA01AA, AAA50AA & AAA51AA | A50077 |
| Sample loop incl. fittings, 100 µl, PEEK; Spare part for AAA20AA & AAA21AA | A500510 |
| Sample loop incl. fittings, 250µl, SSt for AZURA AS 6.1L & 3950 | A500528 |

Sample needles

| | |
|--|---------|
| Sample needle for preparative autosampler AS 6.1L; Spare part AAA40AA & AAA41AA | A500516 |
| Sample needle for biocompatible autosampler AS 6.1L; Spare part AAA20AA & AAA21AA | A500517 |
| Sample needle for bio-preparative autosampler AS 6.1L; Spare part AAA31AA | A500518 |
| Sample needle for analytical autosampler AS 6.1L; Spare part for AAA00AA, AAA01AA, AAA10AA, AAA11AA, AAA50AA & AAA51AA | A64700 |



A50058

Air needles

| | |
|--|---------|
| Air needle for AS 6.1L; 50 mm protrusion length - standard for AAA40AA & AAA41AA | A500529 |
| Air needle for AS 6.1L; 56 mm protrusion length | A500530 |
| Air needle for AS 6.1L; 62 mm protrusion length - standard for all autosampler versions except AAA40AA & AAA41AA | A50058 |
| Air needle for AS 6.1L; 68 mm protrusion length | A500531 |
| Air needle for AS 6.1L; 74 mm protrusion length | A500532 |
| Air needle for AS 6.1L; 80 mm protrusion length | A500533 |
| Set of air needles for autosampler AS 6.1L, 1 pc. of each length | A50059 |

Optional accessories

| | |
|--|---------|
| Fuse (2.5 A) for AS 6.1L, 2 pcs. | A500534 |
| Rectangular bottle (250 ml, PE) for wash or transport solution | A500535 |
| Waste tube for AS 6.1L, silicone, 1 m | A500536 |
| Waste tube for AS 6.1L, PTFE, 1 m | A500537 |



A0638-6



A0638-7



A18201-3



A15854

Vial kits for analytical & preparative HPLC

Vials: Screw neck vials N9 (ø 11.6 mm), 1.5 ml, clear glass, flat bottom, wide opening, with label & scale Caps: Screw caps blue with septum silicone beige/PTFE white 100 pcs each

A0638-6

Vials: Screw neck vials N9 (ø 11.6 mm), 1.5 ml, amber glass, flat bottom, wide opening, with label & scale Caps: Screw caps blue with septum silicone beige/PTFE white 100 pcs each

A0638-7

Vials: Screw neck vials N9 (ø 11.6 mm), 1.5 ml, clear glass, flat bottom, wide opening, with label & scale Caps: Screw caps blue with septum rubber red/TEF colourless 100 pcs each

A0638-8

Vials: Screw neck vials N9 (ø 11.6 mm), 1.5 ml, amber glass, flat bottom, wide opening, with label & scale Caps: Screw caps blue with septum rubber red/TEF colourless 100 pcs each

A0638-9

Microinserts 0.1 ml for Screw neck vials N9, 1.5 ml, 100 pcs.

A18201-3

Vials: Screw neck vials N18 (ø 22.5 mm), 10 ml, clear glass, round bottom Caps: Screw caps, magnetic, with septum rubber red/TEF colourless 100 pcs each

A15854

Maintenance kits

Maintenance kits include the following device specific items: sample needle, air needle, rotor seal, buffer tubing, syringe, syringe valve and wash position.

Preventive maintenance kit for AS 6.1L/3950 (1240 bar/1000 bar) AAA10AA/AAA11AA/A50060/A500601/A50070/A500701

A5009-1

Preventive maintenance kit for AS 6.1L/3950 (700 bar) AAA00AA/AAA01AA/A50080/A50081

A5009-2

Preventive maintenance kit for AS 6.1L/3950 (345 bar, Bio) AAA20AA/AAA21AA/A50052-1/A50053/A50053-1

A5009-3

Preventive maintenance kit for AS 6.1L/3950 (200 bar, Prep) AAA40AA/AAA41AA/A50054-1/A50056-1

A5009-4

Preventive maintenance kit for AS 6.1L/3950 (200 bar, Prep-Bio) AAA31AA/A50055-2


A5009-5

Preventive maintenance kit for AS-1 (1000 bar) A63500/A63501/A63502

A5009-6

Preventive maintenance kit for AS 6.1L (862 bar) AAA50AA/AAA51AA


A5009-7



A50050



A500502



A500505



A500507

Vial plates

Vial plate for 48x 1.5 ml vials for autosampler 3950 and AZURA® AS 6.1L, 1pcs.

A50050

Vial plate for 84x1.5 ml and 3x10ml vials for autosampler 3950 and AZURA® AS 6.1L, 1pcs.

A500501

Prep vial plate for 12 x 10 ml for autosampler 3950 and AZURA® AS 6.1L, 1pcs.

A500502

Vial plate for 108x 1.5 ml vials for autosampler 3950 and AZURA® AS 6.1L, 1pcs.

A500505


Prep vial plate for 30 x 10ml for autosampler 3950 and AZURA® AS 6.1L, 1pcs.

A500507

Detector accessories



AMC19XA



A4045




A4061V2




A4061XB

Flow cells 1/16"

| | |
|--|---------|
| 0.5 mm path length, 3 μ l, 1/16", 200 bar, stainless steel, classical KNAUER flow cell | A4069 |
| 0.5 mm path length, 3 μ l, 1/16", 100 bar, biocompatible, classical KNAUER flow cell | A4095 |
| 3 mm path length, 2 μ l, 1/16", 300 bar, PressureProof Flow cell cartridge for AZURA® Detector DAD/MWD | AMB18 |
| 3 mm path length, 2 μ l, 1/16", stainless steel, classical KNAUER flow cell | A4042 |
| 3 mm path length, 2 μ l, 1/16", 30 bar, biocompatible, classical KNAUER flow cell | A4045 |
| 10 mm path length, 10 μ l, 1/16", 300 bar, stainless steel, for 50D, S2550 and MW-1, classical KNAUER flow cell | A4061V2 |
| 10 mm path length, 10 μ l, 1/16", 300 bar, stainless steel, with heat exchanger one sided inlet and outlet, classical KNAUER flow cell | A4061XB |
| 10 mm path length, 2 μ l, 1/16", 50 bar, LightGuide Flow cell cartridge for AZURA® Detector DAD/MWD | AMC19XA |
| 10 mm path length, 10 μ l, 1/16", 300 bar, PressureProof Flow cell cartridge for AZURA® Detector DAD/MWD | AMC38 |
| 10 mm path length, 2.4 μ l, 1/16", 100 bar, biocompatible, fiber optic connectors, for PDA-1 | A64150 |
| 50 mm path length, 6 μ l, 1/16", 50 bar, High Sensitivity LightGuide Flow cell cartridge for AZURA® Detector DAD/MWD | AMD59XA |
| 50 mm path length, 10 μ l, 1/16", 100 bar, biocompatible, fiber optics connectors, for PDA-1 | A64151 |




A4066




A4067

Flow cells 1/8"

| | |
|---|-------|
| 2 mm path length, 1/8", 200 bar, stainless steel, changeable to 0.5/1.25 mm, classical KNAUER flow cell | A4066 |
| 2 mm path length, 1/8", 100 bar, biocompatible, changeable to 0.5/1.25 mm, classical KNAUER flow cell | A4067 |



A4068



A4068-2

Flow cells 1/4"

| | |
|--|---------|
| 2 mm path length, 1/4" angular connections, 200 bar, stainless steel, changeable to 0.5/1.25 mm, without fittings, classical KNAUER flow cell | A4068 |
| 2 mm path length, 1/4" straight connections, 200 bar, stainless steel, changeable to 0.5/1.25 mm, without fittings, classical KNAUER flow cell | A4068-2 |



A4044



A4044HT



AMKX8KIT



A4047

Flow cells 1/16" fiber optics

0.5 mm path length, 3 μ l, 1/16", 200 bar, stainless steel, fiber optic connectors, classical KNAUER flow cell A4089

0.5 mm path length, 3 μ l, 1/16", 100 bar, biocompatible, fiber optic connectors, classical KNAUER flow cell A4096

3 mm path length, 2 μ l, 1/16", 30 bar, biocompatible, fiber optic connectors, classical KNAUER flow cell A4047

3 mm path length, 2 μ l, 1/16", 300 bar, 85 °C, stainless steel, fiber optic connectors, classical KNAUER flow cell A4044HT

3 mm path length, 2 μ l, 1/16", 300 bar, stainless steel, fiber optic connectors, classical KNAUER flow cell A4044

10 mm path length, 10 μ l, 1/16", 300 bar, stainless steel, fiber optic connectors, classical KNAUER flow cell A4074

Fiber optics adapter kit for AZURA® Detector DAD/MWD, with fiber optic cables (1x 400 mm and 1x 750 mm) and mounting bracket AMKX8KIT



A4078



A4079

Flow cells 1/8" fiber optics

2 mm path length, 1/8", 200 bar, stainless steel, fiber optic connectors, changeable to 0.5/1.25 mm, classical KNAUER flow cell A4078


2 mm path length, 1/8", 100 bar, biocompatible, fiber optic connectors, changeable to 0.5/1.25 mm, classical KNAUER flow cell A4079



A4104

Nano flow cell

3 mm path length, 6 nl, 375 μ m OD, 50 μ m ID, 300 bar, fused silica, fiber optic connectors A4104



A4081



A4153



A4152



A4152-1

Flow cells larger than 1/8" fiber optics

2 mm path length, 1/4" angular connections, 200 bar, stainless steel, fiber optic connectors, changeable to 0.5/1.25 mm A4081

2 mm path length, 1/4" TRI-Clamp connections, 80 bar, stainless steel, fiber optic connectors, changeable to 0.5/1.25 mm A4153

2 mm path length, 3/8" TRI-Clamp connections, 80 bar, stainless steel, fiber optic connectors, changeable to 0.5/1.25 mm A4152

7 mm path length, 3/8" TRI-Clamp connections, 10 bar, biocompatible, fiber optic connectors A4152-1

2 mm path length, 1/2" TRI-Clamp connections, 80 bar, stainless steel, fiber optic connectors, changeable to 0.5/1.25 mm A4154

10 mm path length, 3/8" TRI-Clamp connections, 10 bar, biocompatible, fiber optic connectors A4154-1

2 mm path length, 3/4" TRI-Clamp connections, 80 bar, stainless steel, fiber optic connectors, changeable to 0.5/1.25 mm A4155



A0740



A0743



A0740HT

Fiber optic cables


Fiber optic cables (2 pc), 750 mm, 2x SMA 905 600/660, polycrystalline A0740

Fiber optic cables (2 pc), custom made sizes, 2x SMA 905 600/660, polycrystalline A0743

Fiber optic cables (2 pc), 750 mm, 2x SMA 905 600/660, polycrystalline, up to 85°C A0740HT



A9842



A9843



A9844

Waste tubing kits

Waste tubing kit for AZURA® RID 2.1L, 1/16" A9841

Waste tubing kit for LightGuide flow cells, 1/16" A9842

Waste tubing kit for UV flow cells, 1/16" A9843

Waste tubing kit for UV flow cells, 1/8" A9844



AZL01



A4072



A4448



A5197

Lamps

| | |
|--|---------|
| Deuterium lamp for Smartline S2500 and S2600 detectors | A4071 |
| Deuterium lamp for Smartline PDA detectors K-2800, S2800 and S2850 | A4447V1 |
| HBST deuterium lamp for AZURA® Detector DAD 6.1L | AZL01 |
| HBST deuterium lamp for PLATINblue MW-1 and PDA-1 detectors | A64210 |
| Deuterium lamp, replacement, for S2520, 10D, 40D, UVD 2.1S, UVD 2.1L, DAD 2.1L, MWD 2.1L | A5193 |
| HBST deuterium lamp for Smartline UV and UV/VIS detectors 2550 and BlueShadow 50D | A5194 |
| Halogen lamp for AZURA® Detector DAD 6.1L | AZL02 |
| Halogen lamp for converting Smartline UV 2500 detector into Smartline VIS 2500 detector | A4073 |
| Halogen lamp for converting Smartline UV 2600 detector into Smartline VIS 2600 detector | A4073XA |
| Halogen lamp for PLATINblue MW-1 detector | A64200 |
| Halogen lamp for PLATINblue PDA-1 detector | A64201 |
| Halogen lamp for Smartline 2500 detector, VIS version | A4072 |
| Halogen lamp for Smartline 2600 detector, VIS Version | A4072XA |
| Halogen lamp for Smartline PDA 2800 and 2850 detectors | A4448 |
| Halogen lamp for Smartline UV/VIS detector 2550 | A5195 |
| Halogen lamp for upgrading Smartline UV detector 2550 to Smartline UV/VIS detector 2550 | A5197 |
| LED for Sedex 80LT and Sedex 85LT light scattering detectors | A07541 |
| Xenon lamp for RF-10AXL fluorescence detector | A0753 |
| Xenon lamp for RF-20A/Axs fluorescence detector | A59210 |



AZZ00OC

External heat exchangers

| | |
|--|---------|
| AZURA® Heat Exchanger for analytical PressureProof flow cells, 14 µl | AZZ00OC |
|--|---------|



A5815

Flow splitters

| | |
|---|-------|
| Adjustable preparative post column flow splitter. Calibrated on customer flow rate, split ratio and eluent. | A5815 |
| Passive post column flow splitter for analytical applications. Calibrated with water for split ratios of 15:1 - 250:1 | A1770 |
| Flow splitter for CM 2.1S when used with flowrates over 100 ml/min | A5813 |
| Adjustable semi-preparative post column flow splitter. Calibrated on customer flow rate, split ratio and eluent | A5816 |

Valve accessories - spare parts

Valves, rotor seals and stators

| Description | Valve | Rotor seal | Stator |
|--|---------|------------|---------|
| 6 Port 2-position valve, 1/16", 1200 bar | AVC28AC | A0611.2 | A205118 |
| 8 Port 2-position valve, 1/16", 1200 bar | AVC38AC | A0611.3 | A205120 |
| 6 Port 2-position valve, 1/16", 100 bar | AVD23AF | A205107 | A205140 |
| 6 Port 2-position valve, 1/16", 240 bar, bioinert | AVD24CE | A205101 | A205102 |
| 6 Port 2-position valve, 1/16", 500 bar | AVD26AE | A205145 | A205140 |
| 8 Port 2-position valve, 1/16", 500 bar | AVD36AE | A205144 | A205142 |
| 6 Port 2-position valve, 1/8", 300 bar | AVE25AE | A205147 | A205146 |
| 6 Port 2-position valve, 1/8", 300 bar, 2 channel | AVE25AI | A205162 | A205146 |
| 6 Port 2-position valve, 1/8", 100 bar, bioinert | AVF23CE | A205157 | A205156 |
| 8 Port 2-position valve, 1/8", 50 bar, bioinert | AVF32CE | A205111 | A205130 |
| Manual 6 Port 2-position valve, 1/16", 240 bar, bioinert | AVG24CE | A205101 | A205102 |
| Manual 6 Port 2-position valve, 1/16", 1200 bar | AVI28AC | A0611.2 | A205118 |
| Manual 8 Port 2-position valve, 1/16", 1200 bar | AVI38AC | A0611.3 | A205120 |
| Manual 6 Port 2-position valve, 1/16", 100 bar | AVJ23AF | A205165 | A205140 |
| Manual 6 Port 2-position valve, 1/16", 500 bar | AVJ26AE | A205145 | A205140 |
| Manual 8 Port 2-position valve, 1/16", 500 bar | AVJ36AE | A205144 | A205142 |
| Manual 6 Port 2-position valve, 1/8", 300 bar | AVK25AE | A205147 | A205146 |
| Manual 6 Port 2-position valve, 1/8", 100 bar, bioinert | AVL23CE | A205157 | A205156 |
| 8 Port multiinjection valve, 1/16", 240 bar, bioinert | AVN94CE | A205131 | A205132 |
| 8 Port multiinjection valve, 1/16", 500 bar | AVN96AE | A205160 | A205161 |
| 16 Port multiposition valve, 1/16", 100 bar | AVQ63AF | A205113 | A205152 |
| 16 Port multiposition valve, 1/16", 500 bar | AVQ66AE | A205151 | A205152 |
| 6 Port multiposition valve, 1/16", 1200 bar | AVR28AC | A0880.2 | A205118 |
| 8 Port multiposition valve, 1/16", 1200 bar | AVR38AC | A0880.4 | A205120 |
| 6 Port multiposition valve, 1/16", 100 bar | AVS23AF | A205109 | A205140 |
| 6 Port multiposition valve, 1/16", 500 bar | AVS26AE | A205139 | A205140 |
| 8 Port multiposition valve, 1/16", 240 bar, bioinert | AVS34CE | A205103 | A205104 |
| 8 Port multiposition valve, 1/16", 300 bar | AVS35AE | A205103 | A205142 |
| 8 Port multiposition valve, 1/16", 500 bar | AVS36AE | A205103 | A205142 |
| 16 Port multiposition valve, 1/16", 50 bar, bioinert | AVS62CE | A205105 | A205106 |
| 16 Port multiposition valve, 1/16", 150 bar, bioinert | AVS63CE | A205105 | A205106 |
| 6 Port multiposition valve, 1/8", 300 bar | AVT25AE | A205148 | A205146 |
| 12 Port multiposition valve, 1/8", 100 bar | AVT53AE | A205155 | A205154 |
| 12 Port multiposition valve, 1/8", 100 bar, bioinert | AVT53CE | A205155 | A205164 |
| 8 Port multiposition valve, 1/8", 50 bar, bioinert | AVU32CE | A205129 | A205130 |
| 8 Port multiposition valve, 1/8", 50 bar, bioinert | AVU32GE | A205129 | A205153 |
| 8 Port multiposition valve, 1/8", 200 bar | AVT34AE | A205149 | A205150 |

Fittings for 1/8" valves of V 4.1 valve generation

The ports of V 4.1 valves for 1/8" capillaries are coned with a UNF 1/4-28 thread. Appropriate fittings are provided with each valve. The tables below give an overview on the available fittings.




A7205



A7206



A7207



A7212

NEW Bushings for 1/8" UNF 1/4-28 coned

| | |
|--|-------|
| 1/8" Bushing, short, for UNF 1/4-28 thread, SSt | A7205 |
| 1/8" Bushing, long, for UNF 1/4-28 thread, SSt | A7206 |
| 1/8" Bushing, long, UNF 1/4-28 thread, SSt, for biconical sealing | A7207 |
| 1/8" Blind fitting, for UNF 1/4-28 thread, SSt | A7208 |
| 1/8" Bushing with integrated sealing ring, for UNF 1/4-28 thread, PEEK | A7209 |
| 1/8" Bushing for biconical sealing, UNF 1/4-28 thread, PEEK | A7210 |
| 1/8" Bushing with integrated seal ring, for UNF 1/4-28 thread, PCTFE | A7211 |
| 1/8" Blind plug, for UNF 1/4-28 thread, PEEK | A7212 |



A7213



A7214



A7215



A7217

NEW Ferrules, seal rings and clamp rings for 1/8" UNF 1/4-28 coned

| | |
|---|-------|
| 1/8" Ferrule for wrench-tight fittings, for ports with UNF 1/4-28 thread, SSt | A7213 |
| Split-grooved clamp ring for 1/8" capillary, for ports with UNF 1/4-28 thread, SSt | A7214 |
| Split-grooved clamp ring for 1/8" capillary, for ports with UNF 1/4-28 thread, Titanium | A7215 |
| Biconical seal ring for 1/8" capillary, for ports with UNF 1/4-28 thread, PFTE | A7216 |
| Seal ring for 1/8" capillary, for ports with UNF 1/4-28 thread, PEEK | A7217 |



A7218



A7219



A7220



A7221

NEW Adapters and couplings for 1/8" UNF 1/4-28 coned

| | |
|---|-------|
| Coupling to connect two 1/8" capillaries, 1/8" (UNF 1/4-28) to 1/8" (UNF 1/4-28), SSt | A7218 |
| Coupling to connect two 1/8" capillaries, 1/8" (UNF 1/4-28) to 1/8" (UNF 1/4-28), Titanium | A7219 |
| Coupling to connect two 1/8" capillaries, 1/8" (UNF 1/4-28) to 1/8" (UNF 1/4-28), PEEK | A7220 |
| Coupling to connect 1/16" with 1/8" capillary 1/8" (UNF 1/4-28) to 1/16" (UNF 10-32), SSt | A7221 |
| Coupling to connect 1/16" with 1/8" capillary, 1/8" (UNF 1/4-28) to 1/16" (UNF 10-32), Titanium | A7222 |
| Coupling to connect 1/16" with 1/8" capillary, 1/8" (UNF 1/4-28) to 1/16" (UNF 10-32), PEEK | A7223 |
| Adapter to connect a capillary with 1/16" OD (thread: 10-32 UNF) to 1/8" (thread: 1/4-28 UNF coned), SSt | A7204 |
| Adapter to connect a capillary with 1/16" OD (thread: 10-32 UNF) to 1/8" (thread: 1/4-28 UNF coned), PEEK | A7224 |

Sample loops

Sample loops 1/16" SST incl. fittings




These stainless steel sample loops are designed to be used for 1/16" injection valves with a UNF 10-32 thread.

| | |
|--|---------|
| Sample loop, 1 µl, stainless steel, 0.1 mm ID | A05642 |
| Sample loop, 2 µl, stainless steel, 0.1 mm ID | A05643 |
| Sample loop, 5 µl, stainless steel, 0.25 mm ID | A05644 |
| Sample loop, 10 µl, stainless steel, 0.25 mm ID | A05645 |
| Sample loop, 20 µl, stainless steel, 0.25 mm ID | A05646 |
| Sample loop, 50 µl, stainless steel, 0.45 mm ID | A05647 |
| Sample loop, 100 µl, stainless steel, 0.45 mm ID | A05648 |
| Sample loop, 200 µl, stainless steel, 1 mm ID | A0565 |
| Sample loop, 500 µl, stainless steel, 1 mm ID | A0566 |
| Sample loop, 1000 µl, stainless steel, 1 mm ID | A0567 |
| Sample loop, 2000 µl, stainless steel, 1 mm ID | A0568 |
| Sample loop, 5000 µl, stainless steel, 1.6 mm ID | A0586-2 |

Sample loops 1/8" SST incl. fittings

These stainless steel sample loops are designed to be used for 1/8" injection valves. Based on the port geometry of the valve we offer two variants. For our older valve generations (e.g. V 2.1 valves) please use the one with M8x1 fittings.

For our current V 4.1 valve generation please choose the variant with UNF 1/4-28 fittings. If you're not sure which sample loop to select, you can check the thread specification for your individual valve on our website.



| | |
|---|---------|
| 1 ml sample loop, stainless steel, 2.2 mm ID, incl. M8x1 fittings | A1043 |
|  1 ml sample loop, stainless steel, 2.2 mm ID, incl. UNF 1/4-28 fittings | A142609 |
| 2 ml sample loop, stainless steel, 1.6 mm ID, incl. M8x1 fittings | A1044 |
|  2 ml sample loop, stainless steel, 1.6 mm ID, incl. UNF 1/4-28 fittings | A142610 |
| 10 ml sample loop, stainless steel, 1.76 mm ID, incl. M8x1 fittings | A0843 |
|  10 ml sample loop, stainless steel, 1.76 mm ID, incl. UNF 1/4-28 fittings | A142611 |

Sample loops 1/16" PEEK incl. fittings

These PEEK sample loops are designed to be used for 1/16" injection valves with a UNF 10-32 thread.

| | |
|--|---------|
| Sample loop, 10 µl, PEEK, 345 bar, 0.25 mm ID | A1058 |
| Sample loop, 20 µl, PEEK, 345 bar, 0.25 mm ID | A1059-1 |
| Sample loop, 20 µl, PEEK, 345 bar, 0.5 mm ID | A1059 |
| Sample loop, 50 µl, PEEK, 240 bar, 0.75 mm ID | A1060 |
| Sample loop, 100 µl, PEEK, 240 bar, 0.75 mm ID | A0508 |
| Sample loop, 200 µl, PEEK, 240 bar, 0.75 mm ID | A1061 |
| Sample loop, 500 µl, PEEK, 240 bar, 0.75 mm ID | A1057 |
| Sample loop 1000 µl, PEEK, 240 bar, 0.75 mm ID | A0423 |
| Sample loop 2000 µl, PEEK, 240 bar, 0.75 mm ID | A0785 |

Sample loops 1/8" PEEK incl. fittings

| | |
|--|---------|
| 5 ml sample loop, PEEK, 50 bar, 1/16" ID, incl. M8x1 fittings | A78980 |
|  5 ml sample loop, PEEK, 50 bar, 1/16" ID, incl. UNF 1/4-28 fittings | A142612 |
| 10 ml sample loop, PEEK, 50 bar, 1/16" ID, incl. M8x1 fittings | A78985 |
|  10 ml sample loop, PEEK, 50 bar, 1/16" ID, incl. UNF 1/4-28 fittings | A142613 |



A1159-2



A1160-1

Variloop for sample loading

| | |
|---|---------|
| KNAUER VariLoop S, 10 ml, 1/16", stainless steel, variable injection volume and multiple injections | A1054-2 |
| KNAUER VariLoop S, 10 ml, 1/8", stainless steel, variable injection volume and multiple injections | A1159-2 |
| KNAUER VariLoop L, 40 ml, 1/16", stainless steel, variable injection volume and multiple injections | A1055-1 |
| KNAUER VariLoop L, 40 ml, 1/8", stainless steel, variable injection volume and multiple injections | A1160-1 |
| Superloop, 150 ml, variable injection volumes and multiple injections, movable piston protects against dilution, 1/16", Glass, 20 bar | A1928 |



A0723

TECH TIP

For full-loop injections an overfilling with sample of two to five loop volumes is recommended to ensure precise and reproducible results. Therefore, choose a syringe that exceeds the loop volume by the mentioned factor.

Injection syringes for 1/16" injection port

| | |
|---------------------------|-------|
| Injection syringe 10 µl | A0723 |
| Injection syringe 25 µl | A0724 |
| Injection syringe 50 µl | A0725 |
| Injection syringe 100 µl | A0726 |
| Injection syringe 250 µl | A0727 |
| Injection syringe 500 µl | A0728 |
| Injection syringe 1000 µl | A0729 |
| Injection syringe 2500 µl | A0730 |



A0653

Luer-Lock glass syringes for 1/8" injection port

| | |
|--------------------------------|-------|
| Luer-Lock glass syringe, 10 ml | A0573 |
| Luer-Lock glass syringe, 20 ml | A0653 |



A0555



A0328



A03281



A0505

Loop filling ports

Loop filling port for the injection port of manual injection valves to make the insertion of different sized injection needles safer.

A0555

Injection Port, stainless steel, 1/16"

A0328

Injection Port, PEEK, 1/16"

A03281

Injection Port, stainless steel, 1/8"

A0505

Injection Port, PEEK, 1/8"

A05051

Osmometry accessories



A01242



A0272



A7013



A02330

Pack of 12 ampules NaCl calibrating solution, 300 mOsmol/kg

A01240

Pack of 12 ampules NaCl calibrating solution, 400 mOsmol/kg

A01241-1

Pack of 12 ampules NaCl calibrating solution, 850 mOsmol/kg

A01250

Pack of 12 ampules NaCl calibrating solution, 100 mOsmol/kg

A01242

Pack of 12 ampules NaCl calibrating solution, 2000 mOsmol/kg

A01248

100 Pack of plastic sample tubes for Semi-Micro Osmometer K-7400S

A02721

500 Pack of plastic sample tubes for Semi-Micro Osmometer K-7400S

A0272

1000 Pack of plastic sample tubes for the Semi-Micro Osmometer K-7400S

A0720

Printer paper for the plain paper printer A3711 (60 m roll)

A7013

Ribbon cartridge for the plain paper printer A3711 (black)

A7014

Cleaning tissue, lint-free, for thermistor cleaning

A02330

Purification accessories



A70054V3



A70054V4




A57024




A57026

Eluent & column heating

| | |
|---|----------|
| Eluent heating device (1 channel), 1/16", temperature range: ambient to 100°C, 5,7" display, clean room compatible | A70054V3 |
| Eluent heating device (2 channels), 1/16", temperature range: ambient to 100°C, 5,7" display, clean room compatible | A70054V4 |
| Temperature controller for column heating sleeve | A57024 |
| Heating sleeve for HPLC column 150 x 20 mm HM D=25..57*L=193 mm 100 °C, 230 V, 200 W, Pt100 | A57026 |
| Heating sleeve for HPLC column 250 x 20 mm HM D=25..57*L=293 mm 100 °C, 230 V, 200 W, Pt100 | A57027 |
| Heating sleeve for HPLC column 150 x 30 mm HM D=38..70*L=203 mm 100 °C, 230 V, 400 W, Pt100 | A57028 |
| Heating sleeve for HPLC column 250 x 30 mm HM D=38..70*L=303 mm 100 °C, 230 V, 500 W, Pt100 | A57029 |
| Heating sleeve for HPLC column 150 x 50 mm HM D=60..100*L=211 mm 100 °C, 230 V, 500 W, Pt100 | A57030 |
| Heating sleeve for HPLC column 250 x 50 mm HM D=60..100*L=311 mm 100 °C, 230 V, 800 W, Pt100 | A57031 |
| Heating sleeve for HPLC costum made up to 350 x 50 mm | A57032 |
| Heating sleeve for HPLC costum made up to 350 x 50 mm (moisture proof, for clean room use) | A57034 |



AZG10



AZG10-1



AZG10-2



A70083

Purification

| | |
|--|----------|
| Pressure control for delta pressure measurement up to 250 ml/min for 1/16" and 1/8". Includes interface box. | AZG10 |
| External pressure sensor up to 250 ml/min for 1/16" and 1/8". | AZG10-1 |
| External pressure sensor for up to 1000 ml/min for 1/4", 0 - 10 bar, analog output | AZG10-2 |
| External pressure sensor for up to 1000 ml/min for 1/4", 0 - 10 bar, LAN | AZG10-3 |
| External pressure sensor for up to 1000 ml/min for 1/4", 0 - 10 bar, LAN, biocompatible | AZG10-4 |
| External pressure sensor for up to 1000 ml/min for 1/8", 0 - 50 bar, LAN | AZG10-5 |
| Air sensor (1/16") for AZURA® Bio LC with one air sensor and wiring for up to 4 air sensors | A70092 |
| Additional air sensor for AZURA® Bio LC for 1/16" tubing | A70092-1 |
| Air sensor (1/8") for AZURA® Bio LC with one air sensor and wiring for up to 4 air sensors | A70093 |
| Additional air sensor for AZURA® Bio LC for 1/8" tubing | A70093-1 |
| Air sensor (1/4") for AZURA® Bio LC with one air sensor and wiring for up to 4 air sensors | A70083 |
| Additional Air Sensor for AZURA® BIO LC for 1/4" tubing | A70083-1 |




Consumables

Fittings and bushings

KNAUER K-connect fittings

The K-Connect system consists of a bushing, a split-grooved clamping ring, and a polymer sealing ring. The split-grooved clamping ring and polymer sealing are slipped over the capillary “back to back”, while the bushing tightens all parts. K-Connect fingertight fittings can optionally be tightened further using wrenches if a higher backpressure resistance is needed.



| Article number | Description | Material | For capillary OD [inch] | Thread | Included ferrules | Max.back-pressure [bar] | Amount in set | Picture |
|----------------|--------------------------|----------|-------------------------|-----------|--|-------------------------|---------------|---|
| A9646 | Fingertight Fitting,long | PEEK | 1/16" | UNF 10/32 | Biconical sealing rings A1022 | n/a | 2 |  |
| A9646-1 | Fingertight Fitting,long | PEEK | 1/16" | UNF 10/32 | Biconical sealing rings A1070 | n/a | 10 | |
| A9645 | Fingertight Fitting,long | SST | 1/16" | UNF 10/32 | Split-grooved clamping rings A0484 and polymer Sealing rings A0139 | 1200 | 2 |  |
| A9645-1 | Fingertight Fitting,long | SST | 1/16" | UNF 10/32 | Split-grooved clamping rings A0484 and polymer Sealing rings A0139 | 1200 | 10 | |
| A9647 | Standard Fitting | SST | 1/16" | UNF 10/32 | Stainless steel ferrules A0110 | 1200 | 2 |  |
| A9647-1 | Standard Fitting | SST | 1/16" | UNF 10/32 | Stainless steel ferrules A0110 | 1200 | 10 | |

Flat bottom fittings


| Article number | Description | Material | For capillary OD [inch] | Thread | Included ferrules | Amount in set | Picture |
|----------------|--|----------------------|-------------------------|--------|-------------------|---------------|---|
| A5829 | Bushings flat bottom, super flangeless | PEEK | 1/8" | 1/4-28 | without ferrules | 10 |  |
| A58291 | Bushings flat bottom, super flangeless | PEEK | 1/16" | 1/4-28 | without ferrules | 10 |  |
| A58292 | Ferrules for super flangeless fittings, with lock ring | PEEK/Stainless steel | 1/16" | | | 10 |  |
| A58293 | Ferrules for super flangeless fittings, with lock ring | PEEK/Stainless steel | 1/8" | | | 10 |  |
| A58294 | Ferrules for super flangeless fittings, with lock ring | ETFE/Stainless steel | 1/8" | | | 10 |  |

Dynaseal Fittings

The DYNASEAL connecting system connects capillaries made out of stainless steel and PEEK, as well as PTFE and Tefzel tubings with low dead volumes. It allows maintenance-free operation and provides a long life. Suitable for UNF-threads of type 10/32.


The system consists of a bushing ①, a split-grooved clamping ring ③ and a polymer ferrule ④. The split-grooved clamping ring and polymer ferrule are slipped over the capillary ② "back to back", while the bushing tightens all parts. Thus, leak-free operation is made possible.

DYNASEAL connections are pressure stable up to 450 bar. DYNASEAL can be optionally used with double-cone sealings made out of PEEK. In this case, pressure stability is accordingly reduced to 150 bar.










| Article number | Description | Material | For capillary OD [inch] | Thread | Included ferrules | Max.back-pressure [bar] | Amount in set | Picture |
|----------------|-----------------------------------|----------|-------------------------|-----------|--|-------------------------|---------------|---|
| A0108 | Dynaseal bushings, short | SST | 1/16" | UNF 10/32 | Split-grooved clamping rings A0484 & polymer Sealing rings A0139 | 450 | 4 |  |
| A1021 | Dynaseal bushings, short | SST | 1/16" | UNF 10/32 | without ferrules | depends on ferrule | 10 |  |
| A0181 | Dynaseal bushings, long | SST | 1/16" | UNF 10/32 | Split-grooved clamping rings A0484 & polymer Sealing rings A0139 | 450 | 3 |  |
| A1064 | Dynaseal bushings, long | SST | 1/16" | UNF 10/32 | without ferrules | depends on ferrule | 5 |  |
| A1020 | Dynaseal bushings, short | SST | 1/16" | UNF 10/32 | Biconical sealing rings A1022 | 150 | 10 |  |
| A1069 | Dynaseal bushings, long | SST | 1/16" | UNF 10/32 | Biconical sealing rings A1022 | 150 | 5 |  |
| A0736 | Dynaseal bushings, long | SST | 1/8" | M8x1 | Split-grooved clamping rings A1239 & polymer Sealing rings A0232 | n/a | 4 |  |
| A0735 | Dynaseal bushings, long | SST | 1/8" | M8x1 | without ferrules | depends on ferrule | 4 |  |
| A0644 | Dynaseal bushings, short | SST | 1/8" | M8x1 | Split-grooved clamping rings A1239 & polymer Sealing rings A0232 | n/a | 4 |  |
| A1067 | Dynaseal bushings, short | SST | 1/8" | M8x1 | without ferrules | depends on ferrule | 4 |  |
| A1201 | Dynaseal bushings, long, hex head | SST | 1/8" | M8x1 | without ferrules | depends on ferrule | 4 |  |
| A1201-1 | Dynaseal bushings, long, hex head | PEEK | 1/8" | M8x1 | without ferrules | depends on ferrule | 4 |  |

Standard fittings, stainless steel

| Article number | Description | Material | For capillary OD [inch] | Thread | Included ferrules | Max.back-pressure [bar] | Amount in set | Picture |
|----------------|-------------------------------------|----------|-------------------------|-----------|-------------------|-------------------------|---------------|---|
| A0112 | Bushings short, wrench caliber 1/4" | SST | 1/16" | UNF 10/32 | without ferrules | 1200 | 10 |  |
| A0113 | Bushings short, wrench caliber 1/4" | SST | 1/16" | UNF 10/32 | without ferrules | 1200 | 25 | |
| A0115 | Bushings long, wrench caliber 1/4" | SST | 1/16" | UNF 10/32 | without ferrules | 1200 | 3 | |
| A0116 | Bushings long, wrench caliber 1/4" | SST | 1/16" | UNF 10/32 | without ferrules | 1200 | 10 | |
| A0830 | Bushings, wrench caliber 10 | SST | 1/8" | M8x1 | without ferrules | n/a | 10 | |

Standard fittings, PEEK & polymer

| Article number | Description | Material | For capillary OD [inch] | Thread | Included ferrules | Max.back-pressure [bar] | Amount in set | Picture |
|----------------|--------------------------|----------|-------------------------|-----------|-------------------------|-------------------------|---------------|---|
| A0141 | Bushings knurled, short | Polymer | 1/16" | UNF 10/32 | without ferrules | depends on ferrule | 10 |  |
| A0142 | Bushings knurled, short | Polymer | 1/16" | UNF 10/32 | without ferrules | depends on ferrule | 30 | |
| A0144 | Bushings knurled, long | Polymer | 1/16" | UNF 10/32 | without ferrules | depends on ferrule | 10 |  |
| A0145 | Bushings knurled, short | Polymer | 1/16" | UNF 10/32 | integrated sealing cone | n/a | 10 |  |
| A0584 | Bushings short | PEEK | 1/16" | UNF 10/32 | integrated sealing cone | n/a | 10 |  |
| A0733 | Bushings short | Polymer | 1/8" | M8x1 | integrated sealing cone | n/a | 10 |  |
| A2501 | Bushings short, hex-head | PEEK | 1/16" | UNF 10/32 | integrated sealing cone | 350 | 1 |  |
| A25011 | Bushings short, hex-head | PEEK | 1/16" | UNF 10/32 | integrated sealing cone | 350 | 5 | |
| A2502 | Bushings long, hex-head | PEEK | 1/16" | UNF 10/32 | integrated sealing cone | 350 | 1 |  |
| A25021 | Bushings long, hex-head | PEEK | 1/16" | UNF 10/32 | integrated sealing cone | 350 | 5 | |



A0146



A0582



A0734

Blind fittings / Plugs

| | |
|--|-------|
| 10 Blind plugs, 1/16", knurled, UNF 10-32, short, PETP | A0146 |
| 30 Blind plugs, 1/16", knurled, UNF 10-32, short, PETP | A0147 |
| 10 Blind plugs, 1/16", knurled, UNF 10-32, short, PEEK | A0582 |
| 10 Blind plugs, 1/8", knurled, M8x1, short, PETP | A0734 |

Ferrules and clamping rings



A0139



A1062



A0232

Sealing rings

| | |
|---|-------|
| 30 Sealing rings for capillaries with 1/16" OD, PETP | A0139 |
| 100 Sealing rings for capillaries with 1/16" OD, PETP | A0140 |
| 10 Sealing rings for capillaries with 1/16" OD, PEEK | A1062 |
| 10 Sealing rings for capillaries with 1/8" OD, PETP | A0232 |
| 10 Sealing rings for capillaries with 1/8" OD, PEEK | A1063 |



A1070



A1022



A0738

Biconical sealing rings

| | |
|--|-------|
| 10 Biconical sealing rings for 1/16", PEEK | A1070 |
| 10 Biconical sealing rings for 1/16", PETP | A1022 |
| 10 biconical sealing rings for 1/8", PETP | A0738 |



A0110



A0874



A01101

Ferrules for capillaries

| | |
|---|--------|
| 30 Ferrules for capillaries with 1/16" OD, stainless steel | A0110 |
| 100 Ferrules for capillaries with 1/16" OD, stainless steel | A0111 |
| 10 Ferrules for capillaries with 1/8" OD, stainless steel | A0874 |
| 10 Ferrules for capillaries with 1/16" OD, Hastelloy | A01101 |
| 10 Ferrules for capillaries with 1/16" OD, titanium | A01102 |



A0112



A0115

Bushings for capillaries, SST

| | |
|---|-------|
| 10 Bushings for capillaries with 1/16" OD, stainless steel, wrench caliber 1/4", UNF 10-32, short | A0112 |
| 25 Bushings for capillaries with 1/16" OD, stainless steel, wrench caliber 1/4", UNF 10-32, short | A0113 |
| 3 Bushings for capillaries with 1/16" OD, stainless steel, wrench caliber 1/4", UNF 10-32, long | A0115 |
| 10 Bushings for capillaries with 1/16" OD, stainless steel, wrench caliber 1/4", UNF 10-32, long | A0116 |
| 10 Bushings for capillaries with 1/8" OD, M8x1, wrench caliber 10, stainless steel | A0830 |



A0484



A1239

Split-grooved clamping rings

| | |
|--|-------|
| 4 Split-grooved clamping rings for capillaries with 1/16" OD | A0484 |
| 4 Split-grooved clamping rings for capillaries with 1/8" OD | A1239 |
| 100 Split-grooved clamping rings for capillaries with 1/16" OD | A0482 |

Couplings & adapters



A0148



A0233



A0233-1



A1407

Polymer couplings

Coupling to connect 2 capillaries with 1/16" OD (material: PEEK/PETP, thread: UNF10-32), including 2 bushings and sealing rings, 0.5 mm bore, suitable for classical HPLC, 1 pc. A0148

Coupling to connect 2 capillaries with 1/16" OD (material: PEEK/PETP, thread: UNF10-32), including 2 bushings and sealing rings, 0.5 mm bore, suitable for classical HPLC, 5 pcs. A0149

Coupling to connect 2 capillaries with 1/16" OD (material: PEEK, thread: 10-32 UNF), including 2 one-piece PEEK fittings, 0.5 mm bore, suitable for classical HPLC, 1 pc. A0233

Coupling to connect 2 capillaries with 1/16" OD (material: PEEK, thread: 10-32 UNF), without fittings, 0.5 mm bore, suitable for classical HPLC, 1 pc. A0233-1

Coupling to connect 2 capillaries with 1/16" and 1/8" OD (material: PEEK, thread: 10-32 UNF, M8x1), including 2 one piece fittings (1x 1/16", 1x 1/8"), 1 mm bore, 1 pc. A1407

Coupling to connect 2 capillaries with 1/8" OD (material: PEEK, thread: M8x1), including 2 one piece fittings 1/8", 2 mm bore, suitable for preparative HPLC, 1 pc. A14071



A0117V1



A2512



A0845



A0480

Metal couplings

Coupling to connect 2 capillaries with 1/16" OD (material: titanium, thread: 10-32 UNF), including 2 bushings and ferrules, 0.5 mm bore, suitable for classical HPLC A0117V1

Coupling to connect 2 capillaries with 1/16" OD (material: stainless steel, thread: 10-32 UNF), including 2 bushings and ferrules, 0.5 mm bore, suitable for classical HPLC A0117

Coupling to connect 2 capillaries with 1/16" OD (material: stainless steel, thread: 10-32 UNF), including 2 bushings and ferrules, 0.5 mm bore, suitable for classical HPLC, 5 sets A0118

Coupling to connect 2 capillaries with 1/16" OD (material: stainless steel, thread: 10-32 UNF), including 2 bushings and ferrules, 0.5 mm bore, suitable for classical HPLC, 25 sets A0119

Coupling to connect 2 capillaries with 1/8" OD (material: stainless steel, thread: M8x1), including 2 bushings and ferrules, 2 mm bore, suitable for preparative HPLC A2512

Coupling to connect a capillary with 1/16" OD to a capillary with 1/8" OD (material: stainless steel, thread: M8x1, 10-32 UNF), 1 mm bore A2513

Coupling Dynaseal to connect a capillary with 1/16" OD to a capillary with 1/8" OD (material: stainless steel, thread: M8x1, 10-32 UNF), including Dynaseal bushings and ferrules (1x 1/16", 1x 1/8"), 1 mm bore A0485

Coupling Dynaseal to connect 2 capillaries with 1/8" OD (material: stainless steel, thread: M8x1), including 2 Dynaseal bushings and ferrules, 2 mm bore, suitable for preparative HPLC A0480



A58263



A58264



A58265



A58266

SST Swagelok® unions & reducing unions

| | |
|--|---------|
| Union to connect 2 capillaries with 1/4" OD, material: stainless steel, Swagelok® | A58263 |
| Reducer to connect a capillary with 3/8" OD to a capillary with 1/4" OD, material: stainless steel, Swagelok® | A58264 |
| Reducer to connect a capillary with 8 mm OD to a capillary with 1/4" OD, material: stainless steel, Swagelok® | A58265 |
| Reducer to connect a capillary with 1/8" OD to a capillary with 1/4" OD, material: stainless steel, Swagelok® | A58266 |
| Male connector to connect a 1/4" OD capillary to a 1/4" male NPT adapter (material: stainless steel) for A9861 | A58267 |
| Male connector to connect a 4 mm OD capillary to a 1/4" male NPT adapter (material: stainless steel) for A9861 | A58268 |
| Male connector to connect a 1/8" OD capillary to a 1/4" male NPT adapter (material: stainless steel) for A9861 | A58269 |
| Reducer to connect a capillary with 1/16" OD to a 1/8" OD pipe, material: stainless steel, Swagelok® | A58270 |
| Reducer to connect a capillary with 1/8" OD to a 1/4" pipe union, material: stainless steel, Swagelok® | A58271 |
| Reducer for 1/4" OD capillary to 1/8" OD pipe socket, material: stainless steel, Swagelok® | A582713 |
| Reducer to connect a 1/16" tube socket to 1/4" pipe union, material: stainless steel, Swagelok® | A58273 |
| Reducer to connect a capillary with 4 mm OD to a 1/8" pipe union, material: stainless steel, Swagelok® | A58282 |

Connectors



A2511



A0120



A58260



A58261

Metal T-connectors

| | |
|--|--------|
| T-connector to connect 3 capillaries with 1/8" OD (material: stainless steel, thread: M8x1), including 3 bushings and ferrules | A2511 |
| T-connector to connect 3 capillaries with 1/16" OD (material: stainless steel, thread: 10-32 UNF), including 3 bushings and ferrules | A0120 |
| T-connector to connect 3 capillaries with 1/8" OD (material: stainless steel, Swagelok®) | A58260 |
| T-connector to connect 3 capillaries with 1/4" OD (material: stainless steel, Swagelok®) | A58261 |
| T-connector to connect 3 capillaries with 1/4" OD (material: titanium, Swagelok®) | A58262 |



A150-1



A2511-1



A0150

Polymer T-connectors

T-connector to connect 3 capillaries with 1/16" OD (material: PETP/POM, thread: 10-32 UNF, coned), inclusive 3 bushings and sealing rings

A0150

T-connector to connect 3 capillaries with 1/16" OD (material: PEEK, thread: 10-32 UNF, coned), without bushings

A150-1

T-connector to connect 3 capillaries with 1/8" OD (material: PEEK, thread: M8x1, coned), including 2 one piece 1/8"-PEEK fittings

A2511-1



A0121



A1096

SST X-connectors

X-connector to connect 4 capillaries with 1/16" OD (material: stainless steel, thread: 10-32 UNF), including 4 bushings and ferrules

A0121

X-connector to connect 4 capillaries with 1/8" OD (material: stainless steel, thread: M8x1), including 4 bushings and ferrules

A1096

X-connector to connect 4 tubings with 1/4" OD (material: stainless steel, Swagelok®) for 1000 ml/ min Systems

A58272



A0151

Polymer X-connectors

X-connector to connect 4 capillaries with 1/16" OD (material: PEEK, thread: 10-32 UNF), including 4 one-piece fittings

A0151



A5800



A5805



A5804

Pressure release valves

| | |
|--|---------|
| Pressure Release Valve Kit for AZURA® pump P 2.1L and 80P (25 to 50 bar), 1/8", stainless steel, cross piece titanium | A5800 |
| Pressure Release Valve Kit for AZURA® pump P 2.1L and 80P (3,4 to 24 bar), 1/8", stainless steel, cross piece titanium | A5801 |
| Pressure Release Valve for AZURA® pump P 2.1L and 80P (without spring), 1/4", stainless steel | A5802 |
| Back-Pressure Regulator/pressure release valve for 1/8" and 1/16" OD tubing, 134 µl volume, PEEK, provides a constant back-pressure of 1.4 bar (20 psi), contains pressure release valve tee and fittings for 1/8" and 1/16" | A5804 |
| Back-Pressure Regulator/pressure release valve for 1/16" OD tubing, 134 µl volume, PEEK, provides a constant back-pressure of 0.3 bar (5 psi), contains Y assembly and fittings | A5804-1 |
| Back-Pressure Regulator/pressure release valve kit for 1/16" OD tubing, stainless steel, provides a constant back-pressure of 3 bar, contains pressure release valve tee and fittings for 1/16" | A5805 |
| Back-Pressure Regulator/pressure release valve for 1/16" OD tubing, stainless steel, provides a constant back-pressure of 52 bar, contains pressure release valve tee and fittings for 1/16" | A5805-1 |
| Spring for pressure release valve, 25 - 50 bar | M1070 |
| Spring for pressure release valve, 3.4 - 24 bar | M1080 |



A70087



A70088



A70084

Backpressure regulators

| | |
|--|--------|
| Backpressure Regulator for 1/16" OD tubing, 10-32 threads, PEEK, Range 1-20 bar (15-300 psi) | A70087 |
| Backpressure Regulator for 1/16" OD tubing, 10-32 threads, PEEK, Range 20-103 bar (300-1500 psi) | A70088 |
| Backpressure Regulator for 1/16" OD tubing, 10-32 threads, stainless steel, Range 90-300 bar (1300-4200 psi) | A70084 |
| Spare membranes for Backpressure Regulators A70084, A70087, A70088 | A70082 |

Capillaries and Start up kits



A0130



KNAUER Capillaries, straight




Capillaries 1/16", SST

| | |
|--|-------|
| Stainless steel, 1/16" OD, 0.1 mm ID, 300 cm length, 1 pcs. | A0130 |
| Stainless steel, 1/16" OD, 0.25 mm ID, 300 cm length, 1 pcs. | A0131 |
| Stainless steel, 1/16" OD, 0.5 mm ID, 300 cm length, 1 pcs. | A0132 |
| Stainless steel, 1/16" OD, 0.7 mm ID, 300 cm length, 1 pcs. | A0133 |
| Stainless steel, 1/16" OD, 1 mm ID, 300 cm length, 1 pcs. | A0134 |
| Stainless steel, 1/16" OD, 0.1 mm ID, 10 cm length, 10 pcs. | A0123 |
| Stainless steel, 1/16" OD, 0.1 mm ID, 20 cm length, 10 pcs. | A0124 |
| Stainless steel, 1/16" OD, 0.1 mm ID, 30 cm length, 10 pcs. | A0125 |
| Stainless steel, 1/16" OD, 0.25 mm ID, 10 cm length, 10 pcs. | A0126 |
| Stainless steel, 1/16" OD, 0.25 mm ID, 20 cm length, 10 pcs. | A0127 |
| Stainless steel, 1/16" OD, 0.25 mm ID, 30 cm length, 10 pcs. | A0128 |







Capillaries 1/16", Titanium

| | |
|---|-------|
| Titanium, 1/16" OD, 0.7 mm ID, 50 cm length, 1 pcs. | A0506 |
|---|-------|

Capillaries 1/4", SST

| | |
|--|----------|
|  Stainless steel, 1/4" OD, 4.6 mm ID, 100 cm length, straight, 1 pcs. | A01322-4 |
|  Stainless steel, 1/4" OD, 4.6 mm ID, 150 cm length, straight, 1 pcs. | A01322-5 |
|  Stainless steel, 1/4" OD, 4.6 mm ID, 200 cm length, straight, 1 pcs. | A01322-6 |

Capillaries 1/8", SST

| | |
|--|---------|
| Stainless steel, 1/8"OD, 1.6 mm ID, 150 cm length, oval bent, 1 pcs. | A0639 |
| Stainless steel, 1/8" OD, 2.2 mm ID, 150 cm length, oval bent, 1 pcs. | A0640 |
|  Stainless steel, 1/8" OD, 2.2 mm ID, 100 cm length, straight, 1 pcs. | A0640-4 |
|  Stainless steel, 1/8" OD, 2.2 mm ID, 150 cm length, straight, 1 pcs. | A0640-5 |
|  Stainless steel, 1/8" OD, 2.2 mm ID, 200 cm length, straight, 1 pcs. | A0640-6 |
|  Stainless steel, 1/8"OD, 1.6 mm ID, 100 cm length, straight, 1 pcs. | A0639-4 |
|  Stainless steel, 1/8"OD, 1.6 mm ID, 150 cm length, straight, 1 pcs. | A0639-5 |
|  Stainless steel, 1/8"OD, 1.6 mm ID, 200 cm length, straight, 1 pcs. | A0639-6 |



Analytical K-Connect start-up kits, SST, 1/32" capillaries with fitting sleeves for 1/16" connections

| | |
|---|-------|
| Set of precut capillaries 0.1 mm, fittings, connectors and adapters, red | AZF40 |
| Set of precut capillaries 0.18 mm ID, fittings, connectors and adapters, yellow | AZF50 |
| Set of precut capillaries 0.45 mm, fittings, connectors and adapters, black | AZF60 |
| Set of precut capillaries 0.1 and 0.18 mm ID, fittings, connectors and adapters, red/yellow | AZF80 |



Capillaries, AZURA® Analytical K-Connect, SST 1/32" with fitting sleeve for 1/16" connections

| | |
|--|-------|
| Stainless steel, 0.1 mm ID, 150 mm length, red | AZF41 |
| Stainless steel, 0.1 mm ID, 300 mm length, red | AZF42 |
| Stainless steel, 0.1 mm ID, 400 mm length, red | AZF43 |
| Stainless steel, 0.1 mm ID, 700 mm length, red | AZF44 |
| Stainless steel, 0.1 mm ID, 900 mm length, red | AZF45 |
| Stainless steel, 0.18 mm ID, 150 mm length, yellow | AZF51 |
| Stainless steel, 0.18 mm ID, 300 mm length, yellow | AZF52 |
| Stainless steel, 0.18 mm ID, 400 mm length, yellow | AZF53 |
| Stainless steel, 0.18 mm ID, 700 mm length, yellow | AZF54 |
| Stainless steel, 0.18 mm ID, 900 mm length, yellow | AZF55 |
| Stainless steel, 0.45 mm ID, 150 mm length, black | AZF61 |
| Stainless steel, 0.45 mm ID, 300 mm length, black | AZF62 |
| Stainless steel, 0.45 mm ID, 400 mm length, black | AZF63 |
| Stainless steel, 0.45 mm ID, 700 mm length, black | AZF64 |
| Stainless steel, 0.45 mm ID, 900 mm length, black | AZF65 |

AZURA® Capillary start-up kit, SST

| | |
|--|-------------------------|
| AZURA® start-up kit 1/16" stainless steel, capillary kit | A9849 |
| AZURA® start-up kit 1/16", stainless steel, semi-prep, capillary kit | A9849-1 |
| AZURA® start-up kit 1/8" stainless steel, capillary kit | A9850 |
| AZURA® start-up kit 1/16" stainless steel, 0.25 mm ID precut capillaries | AZF70 |
| AZURA® Accessory kit for ScaleUp system, 1/16" (0.25, 0.5, and 0.7 mm ID), stainless steel | A9850-1 |
| AZURA® Start-Up Kit 1/4" HPG, stainless steel, set of capillaries and fittings | A9850-2 |
| AZURA® Start-Up Kit 1/4" LPG, stainless steel, Set of capillaries and fittings | A9850-3 |



A50041



A9849-2



A70501

AZURA® Capillary start-up kits for special HPLC systems

| | |
|--|---------|
| AZURA® GPC Cleanup Start up kit, Tefzel-(ETFE) tubing, OD 1/16", ID 0.7 mm | A50041 |
| AZURA® Accessory kit for ScaleUp system, 1/16" (0.25, 0.5, and 0.7 mm ID), stainless steel | A9850-1 |
| AZURA® Capillary Start-up kit for educational system 1/16", stainless steel | A9849-2 |
| AZURA® Start-up Kit PEEK, for Analytical HPLC System, up to 5 ml/min or 300 bar | A70501 |

Tubing

Articles grouped under the expression "by the meter" can be shipped in the desired length, by simply ordering it multiple times. E.g. ordering 3.4 x A2528 will result in capillary with a length of minimum 3.4 meters.



Note: If you need tubings with an exact length please contact KNAUER directly.

Tubing start-up kits for FPLC

| | |
|---|---------|
| AZURA® FPLC Start-up kit, PEEK, 1/16" for 10 ml/min FPLC systems | A70500 |
| AZURA® FPLC Start-up kit, transparent FEP, 1/16" for FPLC systems up to 10 ml/min and 20 bar | A70500A |
| AZURA® FPLC Start-up kit, PEEK/Tefzel, 1/16" for 50 ml/min FPLC systems | A70600 |
| AZURA® FPLC Start-up kit, FEP/PEEK, 1/8" for 100 ml/min - 500 ml/min FPLC systems | A70300 |
| AZURA® FPLC Start-up kit, PEEK, 1/8" for 100 ml/min - 500 ml/min FPLC systems, up to 100 bar. | A70300A |
| AZURA® FPLC Start-up kit, PEEK/Tefzel, 1/16" for FPLC systems up to 100 ml/min | A70300B |
| AZURA® FPLC Start-up kit, 1/4" for 1000 ml/min FPLC systems | A70400 |
| AZURA® Start-Up Kit 1/4" HPG, PFA, set of capillaries and fittings | A9850-4 |



A2522

Tubing 1/16" OD, PEEK, by meter

| | |
|--|-------|
| 0.13 mm ID, variable length, max. pressure 420 bar, red striped | A2522 |
| 0.18 mm ID, variable length, max. pressure 400 bar, yellow striped | A2523 |
| 0.25 mm ID, variable length, max. pressure 385 bar, blue striped | A2524 |
| 0.50 mm ID, variable length, max. pressure 350 bar, orange striped | A2525 |
| 0.75 mm ID, variable length, max. pressure 240 bar, green striped | A2526 |
| 1.00 mm ID, variable length, max. pressure 165 bar, grey striped | A2527 |
| 1.40 mm ID, variable length, max. pressure 50 bar, black striped | A2528 |

Tubing 1/8" OD, PEEK, by meter

| | |
|---|-------|
| 0.75 mm ID, variable length, max. pressure 345 bar, natural | A2541 |
| 1.59 mm ID, variable length, max. pressure 220 bar, natural | A2540 |
| 2.00 mm ID, variable length, max. pressure 165 bar, natural | A2542 |



A0182-1

Tubing 1/16" OD, Tefzel™, by meter

| | |
|--|----------|
| 0.25 mm ID, variable length, max. pressure 185 bar | A0182-1 |
| 0.75 mm ID, variable length, max. pressure 115 bar | A0183-1 |
| 1.0 mm ID, variable length, max. pressure 85 bar | A04781-1 |



A0478-1

Tubing 1/8" OD, ETFE, by meter

| | |
|--|---------|
| 1.6 mm ID, variable length, max. pressure 70 bar | A0478-1 |
|--|---------|

Tubing, various OD, PTFE, by meter

| | |
|---|----------|
| 0.45 mm ID, variable length, max. pressure 150 bar, 1.6 mm (1/16") OD | A0152-1 |
| 0.9 mm ID, variable length, 1.6 mm OD | A04782-1 |
| 1.45 mm ID, variable length, max. pressure < 10 bar, 2 mm OD | A0153-1 |
| 1.5 mm ID, variable length, max. pressure 35 bar, 3.2 mm (1/8") OD | A0732-1 |
| 2 mm ID, variable length, 1/8" OD | A0873-1 |
| 3 mm ID, variable length, max. pressure 20 bar, 4 mm OD | A0154-1 |
| 7 mm ID, variable length, 9 mm OD | A1099-1 |
| 1.6 mm ID, variable length, 1/8" OD, black, antistatic | A3306 |
| 4.4 mm ID, variable length, 1/4" OD, black, antistatic | A3307 |

Tubing, various OD, PFA, by meter

| | |
|--|----------|
| PFA tubing, 1/4" OD, 4 mm ID, translucent, max. pressure 15.4 bar, variable length | A31891 |
| PFA tubing, 1/8" OD, 1.6 mm ID, translucent, variable length | A31892 |
| PFA tubing, 1/4" OD, 4.8 mm ID, translucent, max. pressure 15.4 bar, variable length | A31891-1 |
| PFA tubing, 1/8" OD, 2.4 mm ID, translucent, variable length | A31892-1 |

Tubing, various OD, FEP

| | |
|---|---------|
| 2.1 mm ID, 300 cm length, 1/8" OD (FEP tubing) | A9869 |
| 0.81 mm ID, 300 cm length, FEP tubing, 1/16" OD | A9869-1 |

Inline filters



A3381



A00161



B2



A00164-1

Inline filters, SST, for HPLC

| | |
|---|----------|
| Inline Filter (prep.) 5-10 µm, stainless steel, max. flow rate 1000 ml/min (for 1/8" tubing) | A3381 |
| Replacement frit for A3381 5-10 µm, stainless steel, max. flow rate 1000 ml/min | A3381-1 |
| Inline Filter, PEEK body, Stainless steel frit, 1/16", to protect your column, with 2 µm pore size, 3 pcs., easily connected directly to any column | A00161 |
| UHPLC/HPLC precolumn filter, universal, 0.5 µm titanium frit, set of 5, stainless steel body, up to 1034 bar | B2 |
| Inline Filter, stainless steel, frit 0.5 µm, 0.2 µl, for 1/16" capillaries, 0.25 mm bore, up to 1375 bar | A00164 |
| Frit 0.5 µm, 0.2 µl for Inline Filter, stainless steel with 0.25 mm bore up to 1375 bar, 5 pcs. | A00164-1 |



A3378



A3378-1



A00162

Inline filters, biocompatible, for FPLC

| | |
|--|---------|
| Inline Filter, PEEK/Titanium, 1/16", biocompatible, to protect your column, with 2 µm pore size titanium frit | A3378 |
| Inline Filter, PEEK/Titanium, 1/16", biocompatible, to protect your column, with 10 µm pore size titanium frit | A3379 |
| Replacement Frits 2 µm for Inline Filter, PEEK/Titanium, biocompatible | A3378-1 |
| Replacement Frits 10 µm for Inline Filter, PEEK/Titanium, biocompatible | A3379-1 |
| Inline Filter, PEEK body, Titanium frit, 1/16", to protect your column, with 0.5 µm pore size, 3 pcs., easily connected directly to any column | A00162 |
| Inline Filter, PEEK body, Titanium frit, 1/16", to protect your column, with 2 µm pore size, 3 pcs., easily connected directly to any column | A00163 |



A5811

Shut-off valves

| | |
|--|-------|
| Shut-off valve, PEEK, 1/16", including connectors (1/4-28 flat bottom) | A5811 |
| Shut-off valve, PEEK, 1/8", including connectors (1/4"-28 flat bottom) | A5812 |



A1980

Adapters

| | |
|--|-------|
| Luer Adapter to 10-32, ETFE, female Luer to male 10/32 threads for injection, simply screw the adapter in the port of your injection valve | A1980 |
|--|-------|

Safety-caps



A59257



A59257-1



A59259



A59258

Safety caps sets for AZURA analytical systems

| | |
|--|----------|
| for isocratic systems, incl. filters, bottles and fittings | A59257 |
| for LPG systems, incl. filters, bottles and fittings (4 pcs.) | A59257-1 |
| Eluent waste kit for all AZURA® Analytical systems, incl. filter, waste can and cap | A59258 |
| Safety Cap set for AZURA Preparative systems, for one eluent line, incl. filter, bottle and fittings | A59259 |
| Waste Cap set for AZURA Preparative systems, incl. filter, canister and fittings | A59259-1 |



A59260



A59261



A59231



A59234

Safety caps

| | |
|--|--------|
| Eluent Safety Cap, GL45 Thread, 2 ports, 1/4"-28 connection, including air valve and fittings | A59260 |
| Eluent Safety Cap Prep, GL45 Thread, 2 ports for 1/4" tubing, including air filter and fittings | A59261 |
| Eluent Safety Cap Prep, GL45 Thread, 2 ports for 1/4" tubing, including air filter and fittings | A59262 |
| Eluent Safety Cap Filter, spare part, 6 months usable | A59263 |
| VICI Cap, GL45 Thread, 3 ports, 1/4"-28 connection, including O-ring EPDM, nuts and ferrules | A59231 |
| VICI Safety Cap with stopcocks, GL45 Thread, 3 ports, 1/4"-28 connection, including O-ring EPDM, nuts and ferrules | A59234 |
| VICI Safety Cap with stopcocks, GL45 Thread, 4 ports, 1/4"-28 connection, including O-ring EPDM, nuts and ferrules | A59235 |
| VICI Waste Cap, GL45 Thread, 3 ports 1/4"-28 connection, 1 x 10M x 1 for barbed hose adapter, including O-ring EPDM, nuts and ferrules | A59236 |
| VICI Cap, GL45 Thread, 2 ports, 1/4"-28 connection, including O-ring EPDM, nuts and ferrules | A59232 |
| VICI Cap, GL45 Thread, 2 ports, 1/4"-28 connection, including O-ring EPDM, nuts and ferrules | A59230 |
| VICI Safety Cap with stopcocks, GL45 Thread, 2 ports, 1/4"-28 connection, including O-ring EPDM, nuts and ferrules | A59233 |



A59240



A59241



A59246



A59248

Safety caps accessories

| | |
|---|--------|
| VICI Safety Air Inlet Valve with 4 mm filter, fit any VICI cap or VICI safety cap | A59240 |
| VICI Safety Air Inlet Valve with 15 mm filter, fit any VICI cap or VICI safety cap | A59241 |
| VICI Safety Exhaust Filter filled with absorbent, fit any VICI cap or VICI safety cap | A59242 |
| VICI Safety Exhaust Filter with detector, filled with absorbent, fit any VICI cap or VICI safety cap | A59243 |
| O-ring FEP coated for sealing all VICI caps or VICI safety caps, improved chemical resistance | A59244 |
| VICI 1/4-28 flangeless nuts, PPS, for 1/16" tubing, for VICI caps, 10 pcs. | A59245 |
| VICI 1/4-28 flangeless nuts, PPS, for 1/8" tubing, for VICI caps, 10 pcs. | A59246 |
| VICI inverted ferrules, ETFE, for 1/16" tubing, suitable for A59245, for VICI caps, 10 pcs. | A59247 |
| VICI inverted ferrules, ETFE, for 1/8" tubing, suitable for A59246, for VICI caps, 10 pcs. | A59248 |
| VICI plugs, PEEK, 1/4"-28, 1 pcs., for closing unused ports for VICI caps | A59249 |
| VICI Barbed hose adapter for 1/8" tubing, for VICI caps | A59251 |
| Cellulose filter, 0.2 µm, 4 mm diameter for VICI Safety Air Inlet Valve, fit any VICI cap or VICI safety cap | A59252 |
| Cellulose filter, 0.2 µm, 15 mm diameter for VICI Safety Air Inlet Valve, fit any VICI cap or VICI safety cap | A59253 |
| VICI Barbed hose adapter for 8 mm ID tubing, for VICI caps | A59254 |
| for basic solutions in IC, fit any VICI cap or VICI safety cap | A59255 |
| AZURA® Tubing kit with cap and solvent filter (A3375, stainless steel, 10 µm), suitable for all analytical HPLC systems | A9650 |

Lab equipment

Tools



X0219



X0003



X0030

Wrenches & tightening tools

| | |
|--|--------|
| Torque wrench basic tool, 1-25 Nm, without plug-in head | X0219 |
| Open-jaw plug-in head for Torque wrench X0219, 1-17 mm (for 100 - 1000 ml pump head in-/outlet and LPGblock) | X0220 |
| Open-jaw plug-in head for Torque wrench X0219, 1-10 mm (for Smartline I pump heads) | X0221 |
| Open-jaw plug-in head for Torque wrench X0219, 1-13 mm (for 10 - 50 ml Smartline II/ AZURA® pump heads in-/outlet) | X0222 |
| Double open-end wrench, 1/4" and 5/16" | X0003 |
| Double open-end wrench, 8/10 mm | X0030 |
| Double open-end wrenches, 2 pc., 1/4" and 5/16" | A0138 |
| Tightening tools for PEEK fittings, blue, 1/16" fittings 1/4" hex head nut (10-32 threads) | A25030 |
| Tightening tools for PEEK fittings, green, 1/32" fittings 3/16 hex head nut (6-40 threads) | A25031 |



A0569



A0851



A0809



A9865

Capillary and tube cutter

| | |
|--|-------|
| Tube cutter, suitable for all tubes | A0569 |
| Capillary cutter for PEEK capillaries and tubings with OD up to 4 mm | A0851 |
| Metal capillary cutting pliers for 1/16" capillaries | A0809 |
| Metal capillary cutter for 1/8" capillaries | A9865 |



A0137



A9864



A9870

Capillary graters and benders

| | |
|--|-------|
| Capillary grater for degrating of 1/16" stainless steel capillaries, can also be used to remove column filters | A0137 |
| Capillary grater for degrating of 1/8" stainless steel capillaries | A9864 |
| Tube bender for 1/8" and 3/16" tubings with an bend radius of 90° | A9870 |



A1033

Tool kits for AZURA® systems

| | |
|---|---------|
| Tool Kit AZURA® for systems with PEEK or pre-cut capillary kits | A1033 |
| Tool Kit AZURA® for 1/16" systems (stainless steel) | A1033-1 |
| Tool Kit AZURA® for 1/8" systems (stainless steel) | A1033-2 |

Racks



A70010



A70011



Application example (devices not included in the scope of delivery)



A9860

LC racks - space saving solution for AZURA system setup

The Benchtop Racks area solution to install AZURA® L systems at space-limited sites, especially in cold rooms.

| | |
|---|--------|
| Benchtop rack: AZURA® S 300 x 160 x 210 mm (WHD), designed to place an AZURA® S device with a height of 129 mm beneath it | A70016 |
| Benchtop rack: AZURA® L low 480 x 190 x 420 mm (WHD), designed to place AZURA® S or low AZURA® L devices with a height of 150 mm beneath it | A70010 |
| Benchtop rack: AZURA® L high 480 x 430 x 420 mm (WHD), designed to place the Foxy fraction collector or AZURA® L devices beneath it | A70011 |
| Benchtop rack: special manufactured with customized dimensions | A70015 |
| Product Riser AZURA®: Set of 4 feet that lift the device to a height of 28 mm for easy handling of the waste tube of the drainage system | A9860 |


AZURA® Mounting brackets



A9853



A9854-2



A9854-3



A9853-5

| | |
|---|---------|
| Mounting bracket AZURA® L Bio for manual KNAUER injection valve, pH-flowcell and a prepacked column | A9854-1 |
| Mounting bracket AZURA® L for Vici valve drives | A9853-2 |
| Mounting bracket AZURA® S for manual KNAUER injection valve | A9854-2 |
| Mounting bracket AZURA® L for columns with 25 - 29 mm AD | A9853-3 |
| Mounting bracket AZURA® L for AZURA® Valve Unifier VU 4.1 (both-sided) or AZURA® UVD 2.1S and AZURA® CM 2.1S (left-sided on AZURA® L) | A9854-3 |
| Mounting bracket AZURA® L for KNAUER flow cells | A9853-5 |
| Mounting bracket AZURA® L for prep sample loop | A9853-6 |
| Mounting bracket AZURA® L for Hypershear mixing chambers | A9853-8 |
| Mounting bracket AZURA® L for KNAUER manual injection valves | A9853 |

Column holders



A9847



A1319



A70190



A0070A

LC column holder/multi column base

| | |
|---|--------|
| Column holder: Magnetic clip, for all KNAUER columns with 3, 4 and 4.6 mm ID, compatible with all AZURA® devices | A9847 |
| Glass column holder, stand, plate and 2 clamps, can hold one glass column in the dimensions of 10-40 mm ID | A1319 |
| Multi Column Base Bio 60 x 40 x 130cm (w x d x h) for up to 3 MPLC columns with conn. for cooling device | A70190 |
| Multi Column Base including bosshead and clamps, serves as a holder for up to 3 columns with inner diameter up to 50 mm, especially made for preparative column solutions | A0070A |



A4364



A4368



A2820



A2820A

Accessories for LC column holder

| | |
|---|---------|
| 3-finger clamps, long shaft, finger with silicone coating, clamp width 12-100 one piece | A4364 |
| 3-finger clamps, short shaft, finger with silicone coating, clamp width 12-100, one piece | A4364-1 |
| Clamp for Multi Column Base, short shaft, to fix an HPLC column or other accessories to the Multi Column Base, for up to 20 mm ID columns | A4368 |
| Clamp for Multi Column Base, long shaft, to fix an HPLC column or other accessories to the Multi Column Base, for up to 20 ID columns | A2820A |
| Bosshead squared for Multi Column Base, used in combination with clamps with a long shaft on the Multi Column Base | A2820A |

Installation accessories



A1071



A9862

| | |
|---|-------|
| HPLC Standard accessory kit | A1071 |
| Installation Box Kit, Box for small parts, KNAUER file folder and support sticker | A9862 |

Standards for Performance Verification (PV)



A PV procedure is recommended for testing newly installed AZURA® systems as well as for regularly monitoring the system performance.

This table gives an overview of the needed PV document, PV standard and separation column for a specific AZURA® system.

| Backpressure range | Type of detection | Flow cell path length [mm] | Injection: Sample loop volume [µl] | PV document | Article no. of PV standard | Article no. of HPLC column |
|----------------------------------|-------------------|----------------------------|--|--|----------------------------|----------------------------|
| UHPLC systems (max. 1000 bar) | UV, DAD | 10 | 1 - 20 | VPV-001: Analytical HPLC, UV detection | A01260-3 | 10BE181E2F |
| | UV, DAD | 10 | 21 - 100 | VPV-001: Analytical HPLC, UV detection | A01260-2 | 10BE181E2F |
| | UV, DAD | 50 | 1 - 20 | VPV-001: Analytical HPLC, UV detection | A01260-2 | 10BE181E2F |
| | UV, DAD | 50 | 21 - 100 | VPV-001: Analytical HPLC, UV detection | A01260-1 | 10BE181E2F |
| | FLD | all | 1 - 20 | VPV-004: Analytical HPLC, FL detection | A01262-2 | 10BE181E2F |
| | FLD | all | 21 - 100 | VPV-004: Analytical HPLC, FL detection | A01262-3 | 10BE181E2F |
| | RID | all | all | VPV-002: Analytical HPLC, RI detection | A01261-1 | 05WE184E2J |
| HPLC Plus systems (max. 862 bar) | UV, DAD | 10 | 1 - 20 | VPV-001: Analytical HPLC, UV detection | A01260-4 | 15WE181E2J |
| | UV, DAD | 10 | 21 - 100 | VPV-001: Analytical HPLC, UV detection | A01260-3 | 15WE181E2J |
| | UV, DAD | 50 | 1 - 20 | VPV-001: Analytical HPLC, UV detection | A01260-3 | 15WE181E2J |
| | UV, DAD | 50 | 21 - 100 | VPV-001: Analytical HPLC, UV detection | A01260-2 | 15WE181E2J |
| | FLD | all | 1 - 20 | VPV-004: Analytical HPLC, FL detection | A01262-1 | 15WE181E2J |
| | FLD | all | 21 - 100 | VPV-004: Analytical HPLC, FL detection | A01262-2 | 15WE181E2J |
| | RID | all | all | VPV-002: Analytical HPLC, RI detection | A01261-1 | 05WE184E2J |
| | ECD in PAD mode | n/a | 20, 100 | VPV-106: AZURA systems with ECD | A01132 | n/a |
| | ECD in DC mode | n/a | 20, 100 | VPV-106: AZURA systems with ECD | A01273-2 A01273-3 | n/a |
| | UV (normal phase) | 10, 50 | 10, 20, 100 | VPV-009: AZURA HPLC systems in normal phase mode | n/a | 15WE000E2J |
| RID (normal phase) | all | all | VPV-009: AZURA HPLC systems in normal phase mode | n/a | 15WE000E2J | |
| Preparative HPLC systems | UV, DAD | ≤ 2 | all | VPV-007: Preparative HPLC, UV detection | A01264-1 | 05JE181E2J |
| | UV, DAD | > 2 | all | VPV-007: Preparative HPLC, UV detection | A01264-2 | 05JE181E2J |
| | RID | all | all | VPV-008: Preparative HPLC, RI detection | A01265-1 | 05IE184E2J |
| FPLC systems | UV, DAD | all | all | VPV-003: AZURA FPLC systems | A01261-1 | 05WE184E2J |

Mobile Control (Chrom) for Windows 10

With the hand-held Mobile Control and Mobile Control Chrom software you have your AZURA® devices and systems at your fingertips. Remotely control and monitor your devices and enjoy the touchscreen-optimized user interface. Choose Mobile Control as an easy-to-use and cost-effective device control solution!

Mobile Control provides full access to AZURA® devices. Change device settings, set operating parameters, automate device control or check the system status... Mobile Control features all functionalities of a device display.

Do you want to acquire data without the overhead of an advanced chromatographic data system? **Mobile Control Chrom** features data acquisition from AZURA® detectors in addition to full device control.

Why to use Mobile Control (Chrom) software

Only pay for what you use: Mobile Control features basic functions to operate AZURA® devices and systems. Mobile Control can operate dedicated applications which do not require a highly developed and cost-intensive Chromatographic Data System (CDS).

Save space: Mobile Control runs on a tablet. Especially in labs with little space avoiding a desktop PC with keyboard and monitor can be a decisive factor. The touch-optimized user interface allows device control using just your fingers.

Save time: Mobile Control convinces due to an intuitive user interface and a clearly structured menu function. The training period is minimal in comparison to a complex CDS.

Free updates: With every release new features are available in Mobile Control. You can download the current version for free.

Free trial: To evaluate if Mobile Control holds up to your expectations, you can download the software and test the free trial option. Perfect for those who'd like to try before they buy.

Customized software design: Mobile Control is made by KNAUER and can be adapted to the requirements of our OEM partners



This software supports a wide range of instruments:

www.knauer.net/softwarecontrol



For PC hardware & periphery see p. 82

Specifications

Software

| | |
|------------------------------|---|
| Software name | Mobile Control - display software for AZURA® devices without data acquisition Mobile Control Chrom - display software for AZURA® devices with data acquisition |
| Operating system | Windows 10 |
| Software version | Mobile Control v5.5.36, Data Viewer v3.5.37 |
| Supported instruments | Consider release notes (downloads below) |
| Field of application | Display software, device control |

Expandability

| | |
|-------------------------------|-----------------------|
| Stand-alone | yes |
| Multi-user environment | yes |
| Report functions | yes |
| Special features | with tablet and mount |



Free demo version:

www.knauer.net/mobilecontrol

Ordering details:

Software

| | |
|--------|--|
| A9607 | Mobile Control without data acquisition including tablet and mount |
| A9608 | Mobile Control Chrom with data acquisition including tablet and mount |
| A9610 | Mobile Control without data acquisition |
| A9612 | Mobile Control Chrom with data acquisition |
| A9613 | Mobile Control Chrom with data acquisition and column test option |
| A96131 | Coming soon - Mobile Control Chrom with data acquisition and fraction collection option |
| A96132 | Coming soon - Mobile Control Chrom with data acquisition and fraction collection option including tablet and mount |
| A9614 | Upgrading Mobile Control to Mobile Control Chrom gaining data acquisition |

Accessories

| | |
|-----------|--|
| A96181 | USB-LAN ADAPTER Network adapter USB 2.0 ⇔ 10/100 Ethernet for tablets |
| A64809 | WiFi router, 8x LAN GBit RJ-45 ports, 1x WAN GBit RJ-45 port |
| A64809INT | WiFi router, 8x LAN GBit RJ-45 ports, 1x WAN GBit RJ-45 port, power plug UK, US or AUS |
| A64811 | Single device WLAN router for Mobile Control - 1x RJ45, 10/100 MBit; WLAN |
| A9617 | Mobile Control Mount - flexible tablet mount for tablets |

ClarityChrom®

KNAUER ClarityChrom® is a powerful, yet easy-to-use chromatography software (or chromatography data system, CDS) for instrument control, data acquisition and data processing. ClarityChrom is designed for smaller laboratories. It is an economical solution compared to other more complex chromatography software while still offering FDA 21 CFR Part 11 compliance.

ClarityChrom comes as a complete package with LC control and including autosampler control. It is scalable from 1 up to 4 systems; depending on the desired instruments. The built-in fractionation option as well as the optional extensions as SST for automated system tests, PDA for 3D (UV spectra) data handling, GPC analysis, MS and GC control cover a wide range of the requirements for a CDS on a modern lab. KNAUER additionally offers a more advanced fractionation with the KNAUER FRC control module.

ClarityChrom supports all KNAUER devices that can be controlled by software. Please refer to the instrument support list in the Support section of our website, the download link can be found below. Beside this, devices and systems from more than 45 manufacturers can be controlled. Additionally, data acquisition can also be performed with any detector providing a voltage output by simply connecting a KNAUER IFU 2.1 interface box or any other supported A/D converter.


The system suitability (SST) extension automates the calculation of system suitability parameters for system validation and calculates up to 12 parameters and compares the results with the limits the user has set.

The PDA extension allows to acquire and process 3D data from a photo diode array detector (KNAUER PDA detectors are fully supported). The PDA extension provides peak purity analysis and peak identification by spectral library search in self-made or commercial spectra libraries.

The SEC/GPC extension provides interactive and automated gel permeation chromatography analysis, including recalibration and GPC reporting, as well as simplifies the retrieval of GPC data. The GPC extension allows flow rate and multi-detector delay corrections and includes Narrow, Broad and Broad on Narrow calibrations.

ClarityChrom comes with some basic fractionation functionality. The KNAUER-exclusive KNAUER FRC control module for ClarityChrom adds more drivers of several fraction collectors and supports the peak recognition by level and/or slope as well as fractionation by time. Also more advanced functionality as solvent recycling, manual fractionation and rack view with detailed fraction information and chromatogram links are available. The functionality corresponds exactly to the KNAUER preparative functionality of discontinued ClarityChromPrep.

ClarityChrom offers all the necessary operations for an analytical lab. Moreover, the preparative version adds fractionation options to this feature list and allows more flexibility in the lab. ClarityChrom is the best solution for all laboratories searching for an up-to-date and robust software with support of devices from many manufacturers to be flexible in instrumentation but also meet the requirements for modern laboratories.



This software supports a wide range of instruments:
www.knauer.net/en/supp_cc



For PC hardware & periphery
see p. 82

Specifications

Software

| | |
|-----------------------|--|
| Software name | ClarityChrom |
| Extensions / Licenses | PDA / 3D UV, System suitability, Fraction collection, SEC/GPC, Mass spectrometry |
| System architecture | 32-bit CDS |
| Operating system | Windows 10, Windows 8.1, Windows 7, all 32- and 64-bit |

| | |
|--------------------------------------|---|
| Expandability | |
| Stand-alone | Workstation version, max. 4 systems controlled by one computer, max. 3 LC systems, max. 2 systems with PDA or 1 system with MS or special devices per computer |
| Client/server | No Client/Server functionality |
| Multi-user environment | Selectable system of user accounts with independently customizable behavior and appearance for individual users |
| Network environment | Easy offline data sharing (at the file level) among all stations in a local network |
| Fields of application | Analytical and preparative HPLC, GPC/SEC, GC, MS |
| Supported instruments | All KNAUER devices are supported, driver for devices from many other manufacturers are available |
| Instrument connection | Supports RS-232, Ethernet, PCI interface card, A/D-D/A interface |
| Recommended PC hardware | Pentium 2 GHz, 4 GB RAM, 80 GB free hard disk space, separate graphics card if one PC should control more than one system, USB for dongle, connectors as LAN, RS-232 etc. for device control |
| Graphics capabilities | Multiple chromatogram view and overlay, PDA view |
| Integration | 27 integration parameters (peak width, threshold, tangent slope ratio etc.) integration parameters programmable in time, automatic re-integration |
| Calculation types | with/without calibration (int./ext. standard method) |
| Security and GLP | Installation qualification test of the software; FDA 21 CFR Part 11 conformance, validation with virtual detector |
| Instrument control | method-based instrument control, Instrument status display and Direct-Control mode, |
| Calibration | 6 types of calibration curves, up to 20 levels, reference peaks, groups, unlimited number of standards (peaks), LOD, LOQ |
| Chromatogram operations | Overlay view, custom labels and settings, also applying mathematical operations to chromatograms |
| Automation | Sequences, automatic launch of selected commands or applications immediately following chromatogram acquisition - Post run, Batch |
| Presentation of results | Integrated customizable table of results, columns with userdefined calculation, summary table, and export in text or database format |
| Calculations | Custom: 12 predefined mathematical operators, 15 basic and 4 summary functions, special: Kovats indexes for GC, determination of noise/drift, performance calculations |
| Data import and export | ASCII, AIA, dBase |
| Additional options/extensions | |
| FRC option | separate license option; Control of fraction collectors and KNAUER valve drives as fraction collector, fractionation per time/level/slope, rack info with filling level and chromatogram link |
| PDA option | separate license option; 3D chromatogram, peak purity analysis, spectrum search in self-made or commercial spectra library |
| GPC/SEC option | separate license option; molecular weight determination in size exclusion chromatography with various calibration methods |
| System suitability test | separate license option; automates the calculation of system suitability parameters for system validation |
| Note | Autosampler control included |

Ordering details:

Software

| | |
|-------|--|
| A1670 | ClarityChrom® single instrument license for one time base |
| A1674 | ClarityChrom® offline license for data evaluation |
| A1671 | ClarityChrom® additional instrument license on additional time base |
| A1676 | ClarityChrom® option for PDA data processing |
| A1677 | ClarityChrom® system suitability option |
| A1678 | ClarityChrom® option for GPC data processing |
| A1679 | ClarityChrom® option for MS data processing |
| A1682 | ClarityChrom® KNAUER FRC control module for preparative HPLC |
| A1681 | Upgrade for one system from former version to ClarityChrom® |
| A1687 | Upgrade for former ClarityChrom® Prep to latest ClarityChrom® with KNAUER FRC control module |
| A1690 | 30-day trial version of ClarityChrom |
| A1675 | ClarityChrom® university package one offline license |

PurityChrom®



PurityChrom is a chromatography software especially designed for the area of preparative purifications and FPLC applications. PurityChrom provides a user-friendly and clearly structured interface. The **system visualization** offers a graphical representation of the purification system and allows easy handling even of complex flow processes. Furthermore, each device which is displayed in the fluidic scheme can be manually controlled, giving the opportunity to optimize, change and adapt your conditions during the run.

You have the option to create a **method** based on volume, column volume or time. There is also a possibility to pause or to change the method parameters during the run, which gives you complete control over your chromatography process. In PurityChrom you can define important functions in your method with **variables**. This allows you to write methods that can be adapted more flexibly to a specific sample or column, just before the run with only one click. In combination with the **sequence table** a quick and easy method scouting provides you with the best method for your purification problem in less time.

For **fractionation**, you can use a fractionation valve as well as a fraction collector.

Current guidelines and regulations like 21 CFR part 11 are entirely supported. Please check for more information about supported devices the Release Notes of the latest PurityChrom version. With an unlimited number of **free offline licenses**, you can write methods and evaluate runs on any computer of your choice, without blocking the system.

The **basic version** is limited to 3 data channels and the control of eight devices (excl. autosampler). The **upgrade version** (A2652) supports 8 data channels and an unlimited number of devices including an autosampler. The **3D option** (A2654) allows the support of an diode-array detector and the **MS option** (A2655) the usage of a mass spectrometer.



This software supports a wide range of instruments:

www.knauer.net/softwarecontrol



For PC hardware & periphery see p. 82

Specifications

| | |
|--------------------------------------|--|
| Software | |
| Software name | PurityChrom® 5 |
| Operating system | Windows XP, Windows Vista, Windows 7, Windows 8, Windows 10 (English or German only) |
| Expandability | |
| Stand-alone | License for controlling one system |
| Fields of application | FPLC and Prep LC |
| Instrument connection | supports RS-232, Ethernet, A/D-D/A interface |
| Recommended PC hardware | CPU/Memory: Pentium III or higher with at least 1 GHz; 32 or 64 bit; at least 512 MB RAM (Windows XP) and 2 GB (Windows Vista and higher); Graphics: Screen with minimal resolution 1024 x 768 Connectors and Slots: USB for license dongle; COM, USB or LAN according to connected instruments |
| Integration | Real-time analysis of peaks, automatic or manual integration and baseline correction |
| Security and GLP | FDA 21 CFR Part 11 conformance |
| Automation | via sequences and autosampler control files |
| Presentation of results | Individual report configuration as pdf or csv |
| Calculations | Column performance calculations according to DAB |
| Data import and export | Comma Separated Value, AIA/ANDI, ChromStar Slice |
| Special features | User administration |
| Additional options/extensions | |
| FRC option | Included |
| FRC features | Control of fraction collectors and KNAUER valve drives as fractionation valve, fractionation per Time/Level/Slope, rack info with filling level and chromatogram link |
| PDA option | Special licence option; no 3D presentation |
| Note | For autosampler control the upgrade licence is needed |

Ordering details:

Software

| | |
|-------|---|
| A2650 | Basic License for one system |
| A2652 | Extends the Basic License to an unlimited number of controllable devices and 8 data channels, adds autosampler and stacked injections support |
| A2654 | 3D option for a diode-array detector (DAD) |
| A2655 | Mass spectrometry (MS) option for supporting the mass spectrometer 4000 MiD® |
| A2656 | PurityChrom® Maintenance and Support including free updates and 5 hours Software support by KNAUER |

Purity Chrom® MCC / MCC PLUS

PurityChrom® MCC is a special version of our purification software PurityChrom® and is optimized to be used with continuous chromatography systems e.g. SMBC systems. PurityChrom® provides a very user-friendly and clearly structured interface. The system visualization offers a graphical representation and allows easy handling even of complex flow processes. Furthermore each device which is displayed in the fluidic scheme can be manually controlled, giving the opportunity to optimize, change and adapt your conditions also during the run. The new PurityChrom MCC SMB parameter wizard helps you to generate new SMB methods and optimize your parameters while the process is running. With the integrated starting point calculator, you can easily generate you SMB method with the adsorption isotherms of your substances. There is also a possibility to pause your method during a run. The hold function provides you with complete control over your chromatography process.

PurityChrom® MCC Plus is a special software extension enabling monitoring of up to 16 data channels and controlling of up to 8 independent pumps without gradient formation. Accordingly, the software can manage complex, preparative purification systems with an enhanced number of multiple devices. For example, in comparison to other PurityChrom® software packages, the combination of a multiwavelength detector and more than one single UV detector is enabled and up to 8 flowmeters can be controlled in one system.



This software supports a wide range of instruments:

www.knauer.net/softwarecontrol



For PC hardware & periphery see p. 82

Specifications

Software


| | |
|----------------------|--|
| Software name | PurityChrom® MCC / PurityChrom® MCC PLUS |
| Operating system | Windows 10 |
| Field of application | SMB, prep LC |

Expandability

| | |
|-------------------------|---|
| Stand-alone | License for controlling one system |
| Multi-user environment | Selectable system of user accounts with independently customizable behavior and appearance for individual users |
| Instrument connection | supports RS-232, Ethernet, A/D-D/A interface |
| Recommended PC hardware | CPU/Memory: Pentium III or higher with at least 1 GHz; 32 or 64 bit; at least 512 MB RAM (Windows XP) and 2 GB (Windows Vista and higher); Graphics: Screen with minimal resolution 1024 x 768 Connectors and Slots: USB for license dongle; COM, USB or LAN according to connected instruments |
| Security and GLP | FDA 21 CFR Part 11 conformance |
| Automation | via control files |
| Presentation of results | Individual report configuration as pdf or csv |
| Data import and export | Comma Separated Value, AIA/ANDI |
| Special features | controlling of up to 4 independent pumps without gradient formation |

Ordering details:

Device

| | |
|---|---|
| A2659 | PurityChrom® MCC: Software solution to control and monitor AZURA® multi column chromatography systems - SMB |
|  A2657 | PurityChrom® MCC PLUS: Software solution for complex preparative systems without gradient formation |

OpenLAB® CDS EZChrom Edition

OpenLAB CDS EZChrom Edition is the next generation of chromatography data systems and the successor of ChromGate CDS. OpenLAB CDS EZChrom Edition provides chromatography data acquisition, processing and control of GC and LC chromatographs and is used in chromatography operations ranging from single user/single instrument to multi-user/multi-instrument laboratories. It provides support of devices from KNAUER and many other manufacturers.

The basic workstation license can only be installed on one PC and allows for control and data acquisition from one system. The license includes System Suitability, Fraction Collector Control and one year Software Maintenance Agreement (SMA).

The system suitability option allows for test if the system is suitable for particular analysis by testing several parameters as resolution, peak asymmetry and theoretical plates.

The KNAUER fraction collector control option includes the drivers of several fraction collectors and supports fractionation by time, the peak recognition by level and/or slope, also with spectral confirmation. Collet Slices allows for setting a desired volume for each fraction, within the defined fraction vial volume. In the fraction collector configuration the delay volume and the fraction vial volume can be defined. This ensures that the target substance will be collected in the fraction vial and the fraction vial will not overflow. The pump flow rate, which is required for calculation of delay and fraction vial fill level, will be automatically read from the pump's method setup. If a chromatogram of your separation already exists, the required fractionation commands can be derived directly from the chromatogram with a double mouse-click. The rack view gives an overview of the already collected fraction, their volume and retention time. The manual fraction control and the option to use the KNAUER electric valves for fractionation gives you more flexibility. The combination of virtual detector and virtual fraction collector allows for optimizing the fractionation settings from an existing chromatogram of your separations without any physically existing device and, therefore, without the loss of solvent or target substance.

OpenLAB EZChrom Edition and EZChrom Elite are registered trademarks of Agilent Technologies, Inc.



This software supports a wide range of instruments:

www.knauer.net/softwarecontrol



For PC hardware & periphery see p. 82

Specifications

Software

| | |
|------------------------------|--|
| Software name | OpenLAB CDS EZChrom Edition |
| Extensions / Licenses | Fraction collection, System suitability, PDA / 3D UV |
| System architecture | 32-bit CDS |
| Operating system | Depends on CDS version. Latest version, supported by KNAUER drivers, is A.04.09. It runs on Windows 10 Prof./Enterprise, 64-bit and Windows 7 Prof., 32- and 64-bit. |

Additional options/extensions

| | |
|--------------------------------|---|
| FRC option | always included, for preparative HPLC, adds tools for detector controlled fraction collection, solvent and peak recycling, stacked injection, rack view with information about RT and volume |
| FRC features | fractionation can be controlled by time (volume), level, slope including AND/OR combination of these criteria, spectra comparison, local maximum and local minimum, slices, full manual control of fractionation during a run |
| PDA option | 3D chromatogram, peak purity analysis, spectrum search in self-made or commercial spectra library (must be converted in OpenLAB spectral library format) |
| GPC/SEC option | license is discontinued |
| System suitability test | license always included, automates the calculation of system suitability parameters for system validation |

Ordering details:

Software

| | |
|----------|---|
| A2600-1 | OpenLAB® CDS EZChrom Edition workstation for one system with SMA and 4x System Suitability |
| A2610-1 | OpenLAB® CDS EZChrom Edition 3D option for UV detectors MW-1, 2550 and 2600 |
| A2611-1 | OpenLAB® CDS EZChrom Edition 3D UV Option for DAD DAD6.1L, DAD2.1L, PDA-1, S2850 |
| A2618-01 | OpenLAB® CDS EZChrom Edition drivers for 80LT, 85LT, 90LT, 100LT and LC from Sedere |
| A2602-1 | OpenLAB® CDS EZChrom Edition Instrument Control License |
| A2614-1 | OpenLAB® CDS EZChrom Edition for distributed systems - please ask for desired configuration |

Chromeleon™ 7.2 Drivers

Thermo Scientific™ Dionex™ Chromeleon™ is one of the most wide-spread chromatography data systems. Its intuitive handling benefits laboratory workflow and the highly developed algorithms simplify data processing. It offers a broad range of third-party drivers and can be easily used with existing HPLC systems. KNAUER offers drivers for a lot of its devices.

Disclaimer: KNAUER Wissenschaftliche Geräte GmbH is solely responsible for development, testing and support of Thermo Scientific™ Dionex™ Chromeleon™ Chromatography Data System driver software for KNAUER instruments and therefore solely liable for damages associated with the use of this driver software.

Specifications

Computer requirements

| | |
|---------------------------------|--|
| Operating system | Windows 10 Enterprise or Professional Edition. Windows 8.1 Professional, 64-bit; Windows 7 SP1 Professional, Enterprise, 64-bit, (32-bit version is not recommended); Windows Vista SP2 Business, Ultimate, 32-bit (Vista is not recommended) |
| CPU (recommended) | 3 GHz Intel Core i7 or better |
| Memory RAM (recommended) | 8 GB |
| Free Hard Disk Space | 120 GB available, for system with PDA detectors |
| Optical Drive | DVD |
| Display (recommended) | 1280 x 1024, 32-bit color |
| USB Ports | 1 port for USB license key |
| Ethernet Port | 1 port for router (for system connection) |

Ordering details:

Drivers

| | |
|---------|---|
| A1783-2 | Thermo Scientific™ Dionex™ Chromeleon™ "7.2" Driver CD (AZURA® only) |
| A1783-3 | Shimadzu LC Driver for Chromeleon 7.2 Shimadzu CBM-20A required. Only for CBM-20A available. |
| A1783-4 | Sedex Driver for Chromeleon 7.2; For Sedex 85LT / 90LT; Instrument Controller Class 3 necessary |
| A1783-5 | Sedex Driver for Chromeleon 7.2; For Sedex FP / LC / 100LT; Instrument Controller Class 3 necessary |
| A1783-6 | Teledyne Foxy R2 Driver for Chromeleon™ 7.2 |

Enterprise

| | |
|---------|--|
| A1791-1 | Thermo Scientific™ Dionex™ Chromeleon™ 7.2 Enterprise. Complete Software Package incl. Secure Client Lizenz: Data Client, Instrument Operation, Report Designer Pro, Compliance Tools |
| A1792-1 | Thermo Scientific™ Dionex™ Chromeleon™ 7.2 Enterprise license - enables full control of one 3rd Party LC Instrument, includes Instrument Controller, Spectral, Fraction Collection and one Class 3 license |

Workgroup

| | |
|---------|--|
| A1780-2 | Thermo Scientific™ Dionex™ Chromeleon™ 7.2 Bundle Workstation. Complete Software Package incl. License Dongle |
| A1782-2 | Thermo Scientific™ Dionex™ Chromeleon™ 7.2 Instrument Controller Option - Instrument Class 3. Max 2 LC per workstation |
| A1783-8 | Thermo Scientific Dionex Chromeleon™ 7.2 Instrument Controller Option- Fraction collection License |
| A1784-2 | Thermo Scientific™ Dionex™ Chromeleon™ 7.2 Instrument Controller Option - PDA License |
| A1787-2 | Thermo Scientific™ Dionex™ Chromeleon™ 7.2 Instrument Controller Option - MS License |

PC Hardware & periphery

Desktop PCs

| | |
|---|--------|
| Desktop PC (SFF) for OpenLAB® and Chromeleon™ with 24" monitor, English edition Windows 10 Prof. 64-bit English, Intel® Core™ i7, 8 GB RAM, 256 GB SSD, two network cards | A13121 |
| Desktop PC (SFF) for OpenLAB® and Chromeleon™ with 24" monitor, German edition Windows 10 Prof. 64-bit German, Intel® Core™ i7, 8 GB RAM, 256 GB SSD, two network cards | A13111 |
| Desktop PC (SFF) for PurityChrom® and ClarityChrom® with 24" monitor, English edition Windows 10 Prof. 64-bit English, Intel® Core™ i5, 8 GB RAM, 256 GB SSD, two network cards | A13120 |
| Desktop PC (SFF) for PurityChrom® and ClarityChrom® with 24" monitor, German edition Windows 10 Prof. 64-bit German, Intel® Core™ i5, 8 GB RAM, 256 GB SSD, two network cards | A13110 |
| Laptop for OpenLAB®, Windows 10 Prof. German, Intel® Core™ i5, 8 GB RAM, 500 GB HDD, German edition | A13113 |
| Laptop for PurityChrom® and ClarityChrom®, Windows 10, min. Intel® Core™ i3, 8 GB RAM, 256 GB SSD, German edition | A13112 |
| Microsoft Surface Pro for PurityChrom®, Windows 10, Core m3, 4GB RAM, 128 GB SSD, German edition with keyboard, Surface Pen and docking station | A13114 |



Configuration on request

| | |
|---|--------|
| Tell us your requirements and we will figure out the matching CDS workstation. We offer the complete CDS installation and promise you a smooth operation. | A13130 |
|---|--------|



A64809

Network devices

| | |
|---|-----------|
| WiFi router, 8x LAN GBit RJ-45 ports, 1x WAN GBit RJ-45 port | A64809 |
| WiFi router, 8x LAN GBit RJ-45 ports, 1x WAN GBit RJ-45 port, power plug UK, US or AUS | A64809INT |
| 8-port LAN GBit Switch NetGear GS108GE, 8x RJ-45, GBit, Auto MDI-X | A3119 |
| 8-port LAN GBit Switch NetGear GS108GE, 8x RJ-45, GBit, Auto MDI-X, power plug UK, US or AUS | A3119INT |
| 5-port LAN Fast Ethernet Switch NetGear FS105, 5x RJ-45 10/100 MBit | A3126 |
| 5-port Fast Ethernet switch NetGear FS105 5x RJ-45, 100 MBit, power plug UK, US or AUS | A3126INT |
| 16-port LAN GBit Ethernet Switch NetGear GS316 16x RJ-45, 10/100 MBit, Auto MDI-X, power plug UK, US or AUS | A3129 |



AZB00XA

IT accessories

| | |
|--|----------|
| VSCOM USB 4 COM 4 x RS-232 DE9 on USB | A3114 |
| AZURA® Interface Box IFU 2.1 LAN, A/D converter, 4 channels | AZB00XA |
| NEW Ethernet Eventbox for 12 digital inputs and outputs each; only supported under PurityChrom® & PurityChrom MCC Plus® | AZB01 |
| NEW Input cable for Ethernet Eventbox (5 m, M3 plug, open ends with wire end ferrules) | AZB01-01 |
| NEW Output cable for Ethernet Eventbox (3 m, hollow plug, open ends with wire end ferrules) | AZB01-02 |
| NEW Cable for connection of an air sensor to an Ethernet Eventbox (2 m, 2-pole and 3-pole plug) | AZB01-03 |
| RS-232 f/f cable 9-pol nullmodem | A0895 |
| RS-232 m/f cable 9-pol | A0884 |
| APC Smart UPS 1500 VA, uninterruptible power supply for up to 8 devices | A3121 |



A3121

Power cables

| | |
|--|---------|
| Power cable for Europe, 2 m, with rubber connector type C13, 230 V | M1642 |
| Power cable for Switzerland, 2 m, with rubber connector type C13, 230 V | M1597 |
| Power cable for UK, 2.5 m, with rubber connector type C13, 230 V | M1278 |
| Power cable for USA, 2 m, with rubber connector type C13, 115 V | M1651 |
| Power cable, 1.5 m, with rubber connector for UPS APC Smart connector | M2561 |
| Power Supply: Distribution Box 24 V for 6 devices like air sensor, external pressure sensor, IFU 2.1 LAN | AZS80SA |
| NEW EU power cable with 4 cold-device plugs and cover caps | A12345 |



AZS80SA

KNAUER Services

Application Services

With profound application knowledge of analytical and preparative HPLC and FPLC, our team is at your service around the world.

Our experts are pleased to receive your inquiries and requests and will offer attractive customized solutions.



HPLC method development

Qualify, quantify or purify

Do you plan to separate substances by HPLC in order to qualify, quantify or even purify without spending too much time in developing a suitable method? We offer an application and method development service and support you to select a suitable system for your lab.

According to your specifications we prepare an efficient HPLC or FPLC method including advice for an appropriate sample preparation.

HPLC method transfer & optimization

For optimized quality and speed

Do you intend to perform your analyses faster, more efficient and cost effective? We are happy to support you with our profound expertise and experience in liquid chromatography. The team assists in transferring LC applications and methods.

1. Method transfer

We investigate the transfer of your method to one of our HPLC systems. Especially complex separations can cause trouble when transferring them to a different system.

We ensure continuous and consistent quality after the transfer.

2. Method optimization

Using ultra-pure solvents in HPLC can increase the expenses of an analysis substantially. A shift from classic HPLC columns to smaller inner diameters and smaller particle size could cut costs enormously since considerably less solvent is required. We optimize and transfer your LC analyses in order to obtain identical, or even better and faster results, reduce eluent consumption and operating costs.

Rent-an-expert

Get professional assistance in your lab

Some of the numerous chromatographic challenges are better solved in your own lab with your own HPLC equipment. Just order a specialist for your assistance in your lab.

In order to develop the best procedure for your HPLC/FPLC or even purification challenge, we will together compile a concept with you in advance.

In addition, we offer a range of column screening services. If you are unsure which column you need, simply ask our sales department for a quote. We offer analytical column screening support for the most common phases (C18, C18A and C18H) and specialty phases (Eurocat screening or chiral screening).

Chiral Column Screening Services

Chiral column screening and/or method development and optimization

As most chiral separations are not predictable, KNAUER offers a screening service to find the best suiting Eurospher II Chiral column for your chiral separation task.


- Column screening with all available Eurospher II Chiral columns
- Optional: Method optimization
- Results will be reported completely and send as a report
- Method parameters and column specifications will be handed over directly

Eurokat column screening for the analysis of carbohydrates

Not sure which column separates your saccharides best? We offer a screening service for Eurokat columns that are recommended for the separation of sugars and all types of carbohydrates.

- Column screening with all available Eurokat columns
- Optional: Method optimization
- Results will be reported completely and send as a report
- Method parameters and column specifications will be handed over directly



 **Note:** Details and requirements must be discussed previously with KNAUER's application specialists.

Find more information about chiral columns at: www.knauer.net/chiralcscreening

KNAUER Academy

KNAUER has been successfully leading courses for many years for its customers, dealers and sales staff. Our main goal is to familiarize every participant with the latest chromatographic technologies in small groups with practical examples. We offer HPLC courses for beginners and advanced users. In individual courses, participants can receive specialized knowledge, e.g. in UHPLC, FPLC or preparative HPLC. Take part in one of the regularly offered courses or book an individual training on special topics.

Workshops at KNAUER in Berlin or on site (see dates online: www.knauer.net/academy)

Ordering information online or upon request: Tel. +49 30 8097270, E-Mail academy@knauer.net

Update due to Covid-19

We currently only offer online trainings for our most popular courses HPLC basic training and FPLC basic training. Our regulary on-site trainings are on hold due to the Corona pandemic. Please contact us if you are interested in a individual training for your team - we will try to make it possible.

| HPLC Workshops | Description | Article number |
|---------------------------------|---|----------------|
| HPLC Basics (1 day) | Practical work in small groups on compact HPLC systems from installation to system performance verification. | AL0520 |
| HPLC Troubleshooting (1 day) | Participants gain theoretical as well as practical knowledge in troubleshooting detectors, pumps, autosamplers and columns. | AL0510 |
| HPLC Method Development (1 day) | Learn HPLC method development from the beginning and become an HPLC method development pro with our training. | AL0511 |

| Online HPLC Training | Description | Article number |
|--|--|----------------|
| Online HPLC basic training (1/2 day - German language) | Online workshop on HPLC systems from installation to System Performance Verification. Learn HPLC from the beginning and become a HPLC pro with our training. | AL0520-1 |

| FPLC Workshop | Description | Article number |
|---------------------------------------|--|----------------|
| FPLC Basics & Troubleshooting (1 day) | Learn protein purification from the beginning and become a FPLC pro with our training. | AL0580 |

| Online FPLC Training | Description | Article number |
|--|--|----------------|
| Online FPLC training (1/2 day - German language) | This course teaches the basics of FPLC in theory and practice and thus contributes to a better understanding of chromatographic separation techniques. The participant will learn to purify proteins and to develop methods. | AL0580-1 |

| Prep Workshop | Description | Article number |
|---|--|----------------|
| Preparative HPLC Basics & Troubleshooting (1 day) | Learn preparative LC from the beginning and become a preparative LC pro with our training. | AL0570 |

| SMB Workshops | Description | Article number |
|---------------------------------|---|----------------|
| SMB Basics (1 day) | The participants will gain basic knowledge on SMB chromatography. | AL0590 |
| SMB Method Development (2 days) | In this course, SMB method development is explained and practiced using a sample application. | AL0501 |

| Software Workshops | Description | Article number |
|---------------------------------|---|----------------|
| ClarityChrom® Software (2 days) | Deepen your knowledge and improve your analyses through better software skills. The courses provide step by step explanation of the software and all components including advanced functions. | AL0551 |
| OpenLAB® EE Software (2 days) | | AL0550 |
| PurityChrom® Software (1 day) | | AL0552 |

| Maintenance Trainings | Description | Article number |
|---|---|----------------|
| Maintaining KNAUER instruments (individual duration, at KNAUER) | Learn how to perform preventative maintenance on KNAUER equipment by yourself (at KNAUER in Berlin, on site or online). | WM0001 |
| Maintaining KNAUER instruments (at customer) | | WM0002 |

| Individual Workshops (individual duration, at KNAUER in Berlin or on site) | | |
|--|--|--|
| Academy Individual Workshop (individual duration) | Do you wish a different course date or a customized workshop for your department? We will gladly help to achieve your individual goals (at KNAUER in Berlin or on site). | |

Research

Scientific research generates new results and knowledge for industry and society. Currently, KNAUER is involved in different research projects. Obviously, we mostly focus on activities where we can efficiently contribute with our expertise in HPLC technology.

With our research commitments, we intend to generate new knowledge in the field of chromatography as well as learn even more about our own products.

Are you looking for a competent partner in scientific research projects? Do not hesitate to contact us: academy@knauer.net

Compliance

Qualification



Note: Standard procedure for IQ and OQ can be handled differently in individual cases for devices.

Installation Qualification (IQ)

The customer may request the IQ, which is free of charge. In case of a request, the Technical Support of KNAUER or from a provider authorized by KNAUER performs this functionality test during the installation.

The IQ is a standardized document including:

- confirmation of flawless condition at delivery
- check if the delivery is complete
- certification on the functionality of the device

Operation Qualification (OQ)

The Operation Qualification includes an extensive functionality test according to KNAUER standard OQ documents. The Operation Qualification is a standardized document. It is not part of the delivery, please contact the Technical Support in case of request.

The OQ includes the following:

- definition of customer requirements and acceptance terms
- documentation on device specifications
- device functionality check at installation site

Test intervals: To make sure that the device operates within the specified range, the device should be tested regularly. The test intervals depend on the use of the device.

Execution: The test can be carried out either by the Technical Support of KNAUER or from a provider authorized by KNAUER (for a fee).



| Instrument | IQ Document |
|-----------------|--------------------------------|
| all instruments | VIQ-Installation-Qualification |

| Instrument / Software | OQ Doc. |
|---|---------------|
| AZURA® Assistant ASM 2.1L, ASM 2.2L | VOQ-ASM |
| AZURA® AS 6.1L, AS 3950, PLATINblue AS-1 | VOQ-AS |
| AZURA® CM 2.1S | VOQCM21SA |
| AZURA® CT 2.1 Column Thermostat | VOQCT21 |
| AZURA® DAD 6.1L, DAD 2.1L, MWD 2.1L | VOQ-DAD |
| AZURA® RID 2.1L, Smartline S2300 | VOQ-RID-2.1L |
| AZURA® UVD 2.1S, UVD 2.1L | VOQ-Detectors |
| Flow cells | VOQ-Flowcells |
| Fraction collectors | VOQ-FRC |
| Osmometer K-7400 | VOQ-K7400 |
| Osmometer K-7400S | VOQ-K7400S |
| Pumps AZURA®, Smartline, BlueShadow, Platinblue | VOQ-Pumps |
| PurityChrom® | VOQ-PUC |
| RF20A/RF20Axs | VOQ-RF20 |
| System OQ for analytical systems | VOQ-Sys-01 |
| Valves | VOQ-Valves |
| Impigement Jets Mixing Skid | VOQ-IJM |

Performance Verification (PV)

Definition: The document Performance Verification (PV) is part of the quality management system of KNAUER. The Performance Verification includes a qualification test of an AZURA® LC system and must be purchased from the manufacturer. The PV is a standardized KNAUER document and includes:

- Documentation on device specifications
- All necessary method parameters to perform the PV


Goals: The system runs reliably within the documented specifications and the PV is a summary of the results with comments and evaluations.

Target group: The test can be carried out either by the Technical Support of KNAUER, from a provider authorized by KNAUER or by the customer.

| System | Document |
|---|---------------------------------------|
| AZURA® analytical systems with UV detector used in reversed phase mode | VPV-001-AZURA-UV |
| AZURA® analytical systems with RI detector used in reversed phase mode | VPV-002-AZURA-RID |
| AZURA® FPLC systems | VPV-003-AZURA-FPLC |
| AZURA® analytical systems with FLD detector used in reversed phase mode | VPV-004-AZURA-FLD |
| AZURA® SMB Lab and Pilot systems | VPV-005-AZURA-SMB |
| AZURA® preparative systems with UV detector used in reversed phase mode | VPV-007-AZURA-Prep |
| AZURA® preparative systems with RI detector used in reversed phase mode | VPV-008-AZURA-Prep-RID |
| AZURA® systems with UV or RI detector used in normal phase mode | VPV-009-AZURA-HPLC-RI-UV-normal-phase |
| AZURA® systems with ECD detector and flow cell with GC or Au working electrode. | VPV-106-ECD |

Material certification

Upon request customized material certification for all wetted parts with varying degrees of complexity from manufacturer statement (only material) to full documentation (e.g. material certification 3.1, FDA compliance statements).

 **Note:** Retrospective material certification is not possible.

FAT / SAT

The factory acceptance test (FAT) refers to the functional test that is performed upon completion of the manufacturing process to prove the equipment has the same specification and functionality that indicated in the datasheet, specification and purchase order. We are experienced in establishing such test procedures with you before your equipment is shipped.

The acceptance of the equipment at your site (site acceptance test, SAT) is also possible: A technician comes to you and ensures that everything works to your utmost satisfaction. In addition, we can integrate the equipment into the existing production environment, if necessary.



Capillary labeling

Complex HPLC systems with a myriad of valves and variable flow paths can be somewhat confusing. We offer professional capillary labeling upon request, to aid end-users in everyday use.

Support

We are committed to provide the best quality support with experienced staff and technical expertise. All standard user instructions, helpful video tutorials, and a structured section of frequently asked questions is freely accessible on our web page www.knauer.net.

If you need further support, our friendly Support team is happy to help you via e-mail, phone or Team Viewer. They will work with you personally until all issues are resolved.

Contact

Do you have questions about the installation or the operation of your instrument or software?

Support in Germany

(Austria & Switzerland on case-to-case basis):

Phone: +49 30 809727-111 (workdays 9-17h CET)

Email: support@knauer.net

International Support:

Contact your local KNAUER partner for support:

www.knauer.net/en/Support/Distributors-worldwide

Worldwide Technical Services

Our highest goal is to keep your laboratory work as effective and productive as possible. Therefore, we not only pay attention to the highest quality in the development and production of our components and instruments, but also stand by your side after the purchase. With our wide range of services, we are ready to meet any demands to your full satisfaction.

KNAUER offers worldwide quality service of all products, purchased from KNAUER or our authorized partners. All KNAUER Service technicians have completed a specialized service training in the KNAUER headquarter in Berlin, Germany. They are ready to help on site ensuring efficient operation and minimized downtime.

Installation & Instruction

Our experienced KNAUER Service technicians can ensure the proper set-up of your instruments. Get in contact whether you want to use a single device, install a complete system or update your chromatography data system.

KNAUER installations always include introduction in proper handling of the devices as well as tips for self-maintenance and imparting of necessary software knowledge.

On request you may add an IQ, OQ, PV or PQ for compliance (see page 86).

Maintenance

Preventive maintenance has proven to be very successful in ensuring the highest availability of HPLC equipment. Unforeseeable failures of individual system components are thus almost impossible, production processes and laboratory capacities can be planned safely.

We offer maintenance services customized to your needs. You may either ship your instruments to the nearest KNAUER Service facility or contact your local dealer for on-site service of an authorized KNAUER Service technician.

Repair

KNAUER still repairs and maintains the following product lines: the current AZURA®, the former Smartline and PLATINblue devices and - to our best abilities - the Wellchrom equipment which was introduced in the 90s.

If you discover any malfunction of your device, don't worry, we will repair it for you! Please contact your local dealer for shipment matters or ask for an on-site visit of our skilled KNAUER service technicians.

Development Services

Software development

How does your software limit you?

Many devices rely on some kind of software to run and interact with you, either internal software (firmware) or drivers and application software on your PC.

Development of firmware for HPLC devices like

- UHPLC and HPLC pumps
- UV, PDA, RI, detectors
- Autosamplers
- Valves
- Column ovens
- Fraction collectors

Development of device drivers for

- OpenLAB® CDS
- Chromeleon™
- HyStar
- ClarityChrom® (Clarity based)



KNAUER software support for firmware, drivers and software solutions

To provide the most useful tools for your daily work, our team of software engineers combines its expertise in developing firmware, instrument control drivers, as well as application software. KNAUER also has a long experience in customizing instrument operation and in developing drivers for various OEM customers.

Let us know about your software challenges - we will program a solution!

Hardware development

KNAUER has a long experience in customizing scientific equipment according to your needs. With on-site hardware designers, mechanical production and assembly, we can provide tailor made products under certain conditions. Contact us for more information.

Storage of instruments and systems

At times equipment must be removed from your laboratory or you are forced to order equipment before your laboratory is up and running. We can offer storage facilities where your equipment can be stored for future use, giving you peace of mind knowing that you are protecting your investment.

Configuration of your PC

We strongly recommend ordering a KNAUER computer with your HPLC system. However, we understand that sometimes certain constraints do not allow this. We offer a PC configuration service of your PC, in order to assure a safe and reliable installation.

 **Note:** We cannot guarantee installation on a non-KNAUER PC.

Power cable overview

Allocation of power plug types to devices

Every device is supplied with a power plug of the AZURA® series (cold-device plug) in the suitable country-specific version (see Table 2).

Exception of allocation (Table 1)

| Device | Power plug type |
|--|--|
| <ul style="list-style-type: none"> BlueShadow Pump 40P BlueShadow Detector 40D/50D Smartline Degasser (article no. A5328) Osmometer | Smartline series (see Table 2) |
| <ul style="list-style-type: none"> Router Switch | Power plug is supplied. For outside Europe, a suitable adapter is supplied (see Table 2). |
| Degasser (article no. AZE03, AZE03-1, AZE02-1) | Power plug is supplied for US, UK, Europe, Australia. |
| <ul style="list-style-type: none"> Pressure Control (article no. AZG10) Pressure Sensor (article no. AZG10-1) Airsensor (article no. A70092, A70093, A70082) Interface Box (article no. AZB00XA) | Power distributor (article no. AZS80SA) and accessories kit with 1x power plug (article no. F1518) is needed. The distributor can provide power for up to 6 devices. Only one power distributor per system is required. <ul style="list-style-type: none"> Power plug for China: Article no. M3027D Power plug for Australia: Article no. M3027C |
| <ul style="list-style-type: none"> Gas Sensor (article no. A70111) Leak Sensor (article no. A70112) | Power distributor (article no. A70110) , includes 1x power plug |
| Tablet for Mobile Control | The Tablet includes an european power cable. Order M1279/M1277 in addition for an US/UK power cable. |

Overview of country-specific power plugs, routers and switches

If no suitable adapter is available for a specific country, contact the responsible distributor:
www.knauer.net/en/Support/Distributors-worldwide

Overview (Table 2)

| Power plugs/ routers/ switches | Article no. USA | Article no. UK | Article no. CH | Article no. Europe | Article no. Argentina |
|---|----------------------------------|----------------------------------|----------------|--------------------|-----------------------|
| Power plug AZURA® series (cold-device plug) | M1651 | M1278 | M1597 | M1642 | M3233 |
| Power plug Smartline series | M1279 | M1277 | M1479-1 | M1479 | - |
| Router (power plug incl.): Router WLAN, 8x LAN | A64809INT Adapter: M0447V2 | A64809INT Adapter: M0447V1 | - | A64809 | - |
| Switch (power plug incl.): Switch 8x LAN | A3119INT Adapter: M0447V2 | A3119INT Adapter: M0447V1 | - | A3119 | - |
| Switch (power plug incl.): Switch 5x LAN | A3126INT Adapter: M0447V2 | A3126INT Adapter: M0447V2 | - | A3126 | - |



Note: For connecting multiple devices, we provide a special power plug for up to 4 AZURA® devices (Europe), Article no. A12345.

Allocation interfaces to devices

Currently, PCs from KNAUER have no serial interface (RS-232). Thus, to operate the following devices, you must install a serial interface (USB-4COM, Article no. A3114):

- Sedex 85 LT
- Osmometer (only with software)
- Shimadzu RF-20A/Axs, ordered from KNAUER, comes with an RS-232 adapter card PCI-e x1 for desktop computers since this device does not work reliably with USB adapters.



Note: If the tablet for Mobile Control should be connected via LAN and not WLAN, the USB-to-LAN Adapter (article no. A96181) is required.

You find the driver on the KNAUER website: www.knauer.net/en/usb-lan-adapter

Detail overview of devices by power plug type

| AZURA® series (cold-device plug) | Smartline series |
|--|--|
| All devices of AZURA® series | Analytical Pumps 40P (Article no. APC30xx) |
| PCs and monitors | UV Detector 40D/50D |
| Preparative pumps BlueShadow 80P (Article no. APD20xx) | Smartline Degasser (Article no. A5328) |
| Micro devices | Osmometers |
| - BlueShadow Pump 10P/20P | |
| - BlueShadow Detector 10D | |
| - Degasser 20DG (Article no. AZE02) | |
| Detectors | |
| - RF20A (Article no. A59200) | |
| - RF20AXS, CBM-20A (Article no. A59201) | |
| - Gabi Star | |
| - HERM flumo, HERM LB500 | |
| - Sedex85LT, Sedex90LT, Sedex100LT, Sedex LC (Article no. A0754-x) | |
| - CHIRALYSER-MP | |
| Autosamplers | |
| AZURA® Column Thermostat CT 2.1 (Article no. A05852) | |
| Fraction collectors | |
| - Foxy® R1/R2 (Article no. A59100/A59102/A591021) | |
| - LABOCOL Vario-4000 (Article no. A591022/ A591024) | |
| External pressure sensor (Article no. AZG10-2) | |



Document no. V1662, version 2.3. Last accessed on 2020/02/17.

This document is subject to technical changes. For an up-to-date version, please visit:


www.knauer.net/en/cableoverview

System configurator

HPLC/UHPLC by KNAUER

MAKE YOUR PRESELECTION

- UHPLC**
 (SST, max. 1000 bar)
- HPLC**
 (SST, max. 862 bar)
- Bio-Inert**
 (metal-free, max. 400 bar)



ELUENT SELECTION & DELIVERY

- 5 ml/min binary gradient pump P 6.1L (UHPLC)
- 5 ml/min quaternary gradient pump P 6.1L (UHPLC)
- 10 ml/min binary gradient pump P 6.1L
- 10 ml/min quaternary pump P 6.1L
- x solvent selection valve (6 further inlets)

SAMPLE INJECTION

- Manual injection valve
- Autosampler AS 6.1L
- Autosampler AS 6.1L cool/heat

COLUMN SELECTION & THERMOSTAT

- 2 column selection
- 8 column selection
- Column thermostat
- Column kit HPLC
- Column kit UHPLC
- Eluent pre-heating cartridge 0.1 mm ID UHPLC
- Eluent pre-heating cartridge 0.18 mm ID HPLC

DETECTION

- UV/VIS single wavelength
- UV/VIS multiple wavelength
- Conductivity
- Refractive index
- Light Scattering
- A/D-converter (integration of further detectors)
- DAD 2.1L
- DAD 6.1L
- Fluorescence Detector RF-20 A
- Fluorescence Detector RF-20 Axs
- ECD 2.1

ACCESSORIES

- 0.1 mm tubing
- 0.18 mm tubing
- PEEK tubing
- x Back pressure regulator
- Workstation (Windows)

FLOW CELLS FOR UV-DETECTOR

- 10 mm/10 µl Pressure proof
- 10 mm/2 µl LightGuide®
- 50 mm/6 µl LightGuide®
- 3 mm/2 µl (up to 100 ml/min) Pressure proof

SOFTWARE

- ClarityChrom®
- OpenLAB®
- Chromeleon™
- Mobile Control

COMMON APPLICATIONS


- Reversed phase
- Normal phase
- other...
- System Qualification

System configurator

Preparative HPLC by KNAUER

MAKE YOUR PRESELECTION

 SST

 Titanium


BUFFER SELECTION & DELIVERY

- 10 ml/min binary gradient pump P 6.1L
- 10 ml/min quaternary pump P 6.1L
- 50 ml/min binary gradient pump P 6.1L
- x 100 ml/min pump P 2.1L
- x 250 ml/min pump P 2.1L
- x 500 ml/min pump P 2.1L
- x 1000 ml/min pump P 2.1L
- Ternary gradient module for pump P 2.1L
- Binary gradient module for pump P 2.1L
- x solvent selection valve (6 further inlets)

SAMPLE INJECTION

- Injection valve
- Sample pump module
- Sample selection valve: x inlets
- Autosampler AS 6.1L
- Autosampler AS 6.1L cool/heat

COLUMN SELECTION & THERMOSTAT

- Column selection (two columns or one bypass)

DETECTION

- UV/VIS single wavelength
- UV/VIS multiwave length
- DAD 2.1L
- Fluorescence Detector RF-20 A
- Conductivity
- pH
- Refractive index
- Light Scattering
- A/D-converter (integration of further detectors)

FRACTION COLLECTION

- Fractionation valve
- Foxy fraction collector with fixed rack types
- Labocol fraction collector with individual rack types
- Rack for fraction collector
- Flow splitter

ACCESSORIES

- | | | | | |
|-----------------------------|-----------------------------|--|--|--|
| x Airsensor main pump | x Airsensor feed pump | <input type="checkbox"/> Pressure control (2 pressure sensors) | x Back pressure regulator | <input type="checkbox"/> AZURA Organizer |
| x Tubing 1/16" | x Tubing 1/8" | x Tubing 1/4" | <input type="checkbox"/> Workstation (Windows) | |

SOFTWARE

- | | | |
|--|---|---------------------------------------|
| <input type="checkbox"/> ClarityChrom® | <input type="checkbox"/> OpenLAB® | <input type="checkbox"/> PurityChrom® |
| <input type="checkbox"/> Chromeleon™ | <input type="checkbox"/> Mobile Control | |

COMMON APPLICATIONS


- | | |
|---|---|
| <input type="checkbox"/> Reversed phase | <input type="checkbox"/> Normal phase |
| <input type="checkbox"/> other... | <input type="checkbox"/> System Qualification |

System configurator

Bio purification by KNAUER

METHOD

- SEC**
Size Exclusion Chromatography
- AC**
Affinity Chromatography
- IEX**
Ion-Exchange Chromatography
- HIC**
Hydrophobic Interaction Chromatography



BUFFER SELECTION & DELIVERY

- 10 ml/min binary gradient pump P 6.1L
- 10 ml/min quaternary pump P 6.1L
- 50 ml/min binary gradient pump P 6.1L
- x 100 ml/min pump P 2.1L
- x 250 ml/min pump P 2.1L
- x 500 ml/min pump P 2.1L
- x 1000 ml/min pump P 2.1L
- Ternary gradient module for pump P 2.1L
- Binary gradient module for pump P 2.1L
- x Buffer selection valve (8 further inlets)

SAMPLE INJECTION

- Multi-Injection valve
- x Injection valve
- Sample pump module
- Sample selection valve: x inlets
- Biocompatible Autosampler AS 6.1L

COLUMN SELECTION

- Column selection valve up to 50 ml/min (5 columns, one bypass, reverse flow)
- Column selection (two columns or one bypass)
- Column selection high flow (5 columns, one bypass)
- Column selection high flow (7 columns, one bypass, reverse flow)

DETECTION

- UV/VIS single wavelength
- UV/VIS multiwave-length
- Conductivity
- pH
- Fluorescence
- Refractive index
- Light Scattering
- Analog integration of further detectors

FRACTION COLLECTION

- Outlet valve
- Foxy fraction collector with fixed rack types
- Labocol fraction collector with individual rack type
- Rack for fraction collector

COLUMNS & MEDIA

- | | | |
|--|---|--|
| <input type="checkbox"/> SEC : Desalting ml | <input type="checkbox"/> AC : Protein A ml | <input type="checkbox"/> IEX : DEAE - Weak anion exchange ml |
| <input type="checkbox"/> SEC : SEC 75 ml | <input type="checkbox"/> AC : Protein G ml | <input type="checkbox"/> IEX : CM - Weak cation exchange ml |
| <input type="checkbox"/> SEC : SEC 200 ml | <input type="checkbox"/> AC : Ni-NTA ml | <input type="checkbox"/> IEX : Q - Strong anion exchange ml |
| | <input type="checkbox"/> AC : Glutathione ml | <input type="checkbox"/> IEX : SP - Strong cation exchange ml |

ACCESSORIES

- | | | | | |
|------------------------------|------------------------------|--|--|--|
| x Air sensor main pump | x Air sensor feed pump | <input type="checkbox"/> Pressure control (2 pressure sensors) | x Back pressure regulator | <input type="checkbox"/> AZURA Organizer |
| x Tubing 1/16" | x Tubing 1/8" | x Tubing 1/4" | <input type="checkbox"/> Workstation (Windows) | |

KNAUER BlueShadow Pumps and Detectors



Versatile stand alone instruments for your lab and production systems

KNAUER BlueShadow pumps & detectors are the ideal choices for upgrading your existing LC, reaction system or process instruments.

BlueShadow Pumps 40P and 80P



BlueShadow Pump 40P



BlueShadow Pump 80P

Pumps from the BlueShadow line can be integrated into every existing LC system, but they can also be used for high-pressure dosing applications. KNAUER dosing pumps are highly accurate two-piston pumps for applications in the chemical and pharmaceutical industries as well as in research and method development.

They pump and dose aqueous and organic liquids, aggressive media or liquid gases. The metering pumps impress with their high chemical resistance, excellent flow rate precision and low pulsation of the pumped medium in a wide range of applications.

BlueShadow Detectors 40D and 50D



BlueShadow Detector 40D



BlueShadow Detector 50D

Detectors from the BlueShadow line are spectrophotometers that can be used for LC applications, reaction monitoring, and other applications. They are offering excellent technical specifications in a highly flexible and compact design.

The flow cells are easily accessible, can be changed quickly and cover flow rates from 10 $\mu\text{l}/\text{min}$ up to 10 l/min. With the unique fiber optics design of the BlueShadow 40D, the flow cell can also be separated from the detector and directly placed in the stream of the product flow.

For more information on the BlueShadow devices please refer to the KNAUER website or our High-Pressure Dosing Pumps and Accessories selection guide: www.knauer.net/dosingpump-psg.

KNAUER GMP Services NEW

KNAUER Services for Good Manufacturing Practice for biopharmaceutical industry

KNAUER provides equipment for downstream processing in the biopharmaceutical industry such as skids for the formulation of lipid nanoparticles, or chromatographic systems for mRNA purification or continuous chromatography. KNAUER provides a wide range of services to support our customers and to ensure that GMP requirements are met.

KNAUER's GMP services are based on our hardware- and software-solutions; encompassing product safety, quality control and the training of personnel. Risk management, in relation to GMP, is covered by the user:



Further information:
www.knauer.net/gmp

Product safety:

Documentation on the compliance of materials used for wetted parts is an important requirement for product safety. In the bio-pharmaceutical industry, potentially harmful substances must be avoided in liquids for clinical, cosmetic or food applications. Therefore, any materials of the liquid flow path that come into contact with the final product must meet certain criteria. According to our end user's requirements KNAUER can provide compliance with the order (EN 10204-2.1), certificates of compliance on the materials used for wetted parts, and further documentation from the supplier such as 2.1 certificates.

Overview of KNAUER options:

Certificates are available for KNAUER products and selected third party products. Contact sales@knauer.net.

| Type of certificate/statement | Unit of quantity | Article number |
|--|---|----------------|
| Declaration of Compliance with order (EN 10204-2.1) | for 1 order | A0000TDCOO |
| Declaration of Compliance (EN 10204-2.1) with Certificate and compliance of material of wetted parts | for 1 article with less than 5 components | A0000COMS |
| Declaration of Compliance (EN 10204-2.1) with Certificate and compliance of material of wetted parts | for 1 article with 5 or more components | A0000COM |
| Declaration of Compliance (EN 10204-2.1) with Certificate and compliance of material of wetted parts for one pump | for 1 article | A0000COMP |
| Declaration of Compliance (EN 10204-2.1) with Certificate and compliance of material of wetted parts for one pump head | for 1 article | A0000COMPK |
| Declaration of Compliance (EN 10204-2.1) with Certificate and compliance of material of wetted parts for one valve | for 1 article | A0000COMV |
| Declaration of Compliance (EN 10204-2.1) with Certificate and compliance of material of wetted parts for one detector | for 1 article | A0000COMD |
| TSE/BSE Statement: Customized order- and article related | for 1 article | A0000TDTSE |
| Stepfile per device without functional groups | for 1 article | A0000IDSTE |
| Documentation on compliance of material of wetted parts: Compliance with the order (EN 10204-2.1); Certificate of compliance on material of wetted parts; Documentation on references (supplier information of material) and wetted parts (certificates such as 2.1) | customized | A0000TD |

Conversion tables

Dimensions

| mm | inches | inches | mm |
|-------|--------|--------|------|
| 0.10 | .004" | 1/32" | 0.8 |
| 0.12 | .005" | 1/16" | 1.6 |
| 0.15 | .006" | 1/8" | 3.2 |
| 0.25 | .010" | 3/16" | 4.8 |
| 0.40 | .016" | 1/4" | 6.4 |
| 0.50 | .020" | 3/8" | 9.5 |
| 0.75 | .030" | 1/2" | 12.7 |
| 1.00 | .040" | 1" | 25.4 |
| 1.50 | .060" | | |
| 2.00 | .080" | | |
| 4.60 | .180" | | |
| 6.00 | .236" | | |
| 6.40 | .253" | | |
| 7.00 | .276" | | |
| 10.00 | .400" | | |

Tubing volume/length

| Tubing ID | µl/cm | µl/in |
|-----------|-------|-------|
| .004" | 0.08 | 0.21 |
| .005" | 0.13 | 0.32 |
| .010" | 0.51 | 1.29 |
| .015" | 1.14 | 2.90 |
| .020" | 2.03 | 5.15 |
| .025" | 3.17 | 8.04 |
| .030" | 4.56 | 11.58 |
| .040" | 8.11 | 20.59 |
| .060" | 18.24 | 46.33 |
| .070" | 24.83 | 63.06 |
| .085" | 36.61 | 92.99 |

Pressure

| MPa | bar | psi |
|-----|-------|--------|
| 5 | 50 | 725 |
| 10 | 100 | 1 450 |
| 20 | 200 | 2 901 |
| 30 | 300 | 4 351 |
| 40 | 400 | 5 802 |
| 50 | 500 | 7 252 |
| 60 | 600 | 8 702 |
| 70 | 700 | 10 153 |
| 80 | 800 | 11 603 |
| 90 | 900 | 13 054 |
| 100 | 1 000 | 14 504 |
| 110 | 1 100 | 15 954 |
| 120 | 1 200 | 17 405 |
| 130 | 1 300 | 18 855 |
| 140 | 1 400 | 20 306 |
| 150 | 1 500 | 21 756 |
| 160 | 1 600 | 23 206 |
| 170 | 1 700 | 24 657 |
| 180 | 1 800 | 26 107 |
| 190 | 1 900 | 27 558 |
| 200 | 2 000 | 29 008 |

Temperature

| °C | °F | °C | °F | °C | °F |
|-----|-----|-----|-----|-----|-----|
| -40 | -40 | 65 | 149 | 170 | 338 |
| -35 | -31 | 70 | 158 | 175 | 347 |
| -30 | -22 | 75 | 167 | 180 | 356 |
| -25 | -13 | 80 | 176 | 185 | 365 |
| -20 | -4 | 85 | 185 | 190 | 374 |
| -15 | 5 | 90 | 194 | 195 | 383 |
| -10 | 14 | 95 | 203 | 200 | 392 |
| -5 | 23 | 100 | 212 | 205 | 401 |
| 0 | 32 | 105 | 221 | 210 | 410 |
| 5 | 41 | 110 | 230 | 215 | 419 |
| 10 | 50 | 115 | 239 | 220 | 428 |
| 15 | 59 | 120 | 248 | 225 | 437 |
| 20 | 68 | 125 | 257 | 230 | 446 |
| 25 | 77 | 130 | 266 | 235 | 455 |
| 30 | 86 | 135 | 275 | 240 | 464 |
| 35 | 95 | 140 | 284 | 245 | 473 |
| 40 | 104 | 145 | 293 | 250 | 482 |
| 45 | 113 | 150 | 302 | 255 | 491 |
| 50 | 122 | 155 | 311 | 260 | 500 |
| 55 | 131 | 160 | 320 | 265 | 509 |
| 60 | 140 | 165 | 329 | 270 | 518 |

Terms & Conditions

1. Definition of terms

The following terms and conditions apply to every order received by KNAUER and every delivery of goods. This holds as well in case of contradictory buying conditions of the purchaser. Exceptions are only valid when confirmed by KNAUER in writing. Purchase orders are only binding if confirmed by KNAUER in writing.

2. Payment

Deliveries are due and payable, net, within 30 days of invoice date or in advance. Deductions are not allowed. Foreign deliveries must be paid by irrevocable letter of credit or in advance. All bank and transfer fees must be paid by the customer. The consequences arising out of delay are due to statutory provisions. Payments are due irrespective of an eventual notice of defect, except such defects are evidently justified.

3. Delivery

Delivery dates are not binding unless expressly stated in the contract as binding dates. Delay in delivery requires a written reminder and an adequate additional grace period set by the customer. KNAUER is only liable for claims for damages under the requirements of no. 6.

4. Claims

Condition for any warranty claim is the immediate inspection of the goods upon delivery, and complaint towards and damage assessment together with the carrier, and an immediate written complaint to KNAUER. The complaint must be made within five workdays in case of visible defects or losses.

5. Risk liability

Delivery is made at the customer's own risk. As soon as the goods leave KNAUER's plant the risk of accidental loss, destruction or deterioration passes to the customer.

6. Warranty and damages

6.1. Warranty claims

The warranty begins with receipt of the goods. If commissioning has been ordered, after commissioning. In the case of delayed commissioning, the warranty begins at the latest four weeks after receipt of the goods unless the supplier is responsible for delayed commissioning.

The warranty for osmometers and liquid chromatography instruments is limited to two years, excluding glass breakage, damages due to stoppage and consumable materials such as membranes, light bulbs, columns, bushings, gaskets and valves. KNAUER's liability shall be restricted to the replacement of defective material or repair only. Transportation costs are borne by the customer. In case of failure of replacement or repair the customer may demand a reduction in price or cancellation of the contract with respect to the defective material. The customer has to inspect the goods delivered immediately and shall immediately give written notification of any defects to KNAUER, in case of non-obvious defects within 10 working days after delivery at the very latest.

6.2. Claims for damages

The liability of KNAUER shall be restricted to intentional acts and acts of gross negligence and compensation shall only be due for direct, foreseeable damages. Liability for breach of a material, essential duty of the contract, liability because of personal injury, liability according to the stipulations of the German Law on Product Liability and liability for the lack of the condition of the contract goods guaranteed by KNAUER remain unaffected.

7. Third party rights on industrial or other intellectual property

KNAUER shall not be liable for the infringement of third party rights founded on industrial or other intellectual property caused by the use of the delivered goods. The customer is fully responsible for the products manufactured with the goods. In particular KNAUER is not obliged to indemnify and hold harmless the customer from all claims raised by third parties based on the infringement of their industrial or intellectual property rights by the use of the goods.

8. Property rights

The ownership of the goods shall remain with KNAUER until payment in full for all our claims resulting from our business relation is received. In case of improper treatment of the goods or in case of default KNAUER may demand the return of the delivered goods. This demand entails resignation of the contract only if KNAUER declares it explicitly.

Resellers are allowed to sell the goods to third parties in due course of the business. The customer herewith assigns his resale claims against third parties to KNAUER.

9. Export

Instruments and products delivered by KNAUER may not be exported to a country other than of the customer's headquarters without KNAUER's prior written permission.

10. Place of settlement and court of jurisdiction

The place of performance is Berlin. Proper venue for all claims is the competent local court at KNAUER's principal place of business - Berlin. KNAUER reserves the right to sue the customer at his principal place of business.

This agreement shall be governed by the laws of the Federal Republic of Germany excluding the UN-Convention on the International Sale of Goods (CISG).

KNAUER Wissenschaftliche Geräte GmbH
Hegauer Weg 38
14163 Berlin, Germany

These terms and conditions apply since June 1, 2016


KNAUER Brochures



AZURA® Analytical HPLC/UHPLC
(Document no. V7852US)



AZURA® Bio purification
(Document no. V7855US)




AZURA® Preparative HPLC
(Document no. V7820US)



AZURA® SMB systems
(Document no. V7741US)



Freezing point osmometry
(Document no. V7716US)




KNAUER Dosing Pump Selection Guide
(Document no. V7866US)



KNAUER LNP Flyer
(Document no. V7720US)




KNAUER OEM Brochure
(Document no. V7712US)



KNAUER Column Selection Guide
(Document no. V7803US)



All KNAUER brochures:
www.knauer.net/brochures



KNAUER

Wissenschaftliche Geräte GmbH
 Hegauer Weg 38
 14163 Berlin, Germany

+49 30 809727-0
 +49 30 8015010 (Fax)

info@knauer.net
 www.knauer.net

CEO

Alexandra Knauer, CEO
 Carsten Losch, CEO

Commercial register

Berlin-Charlottenburg
 Register No.: 93 HRB 15674
 VAT-ID-No.: DE136737469
 EORI Number DE 2620448
 DUNS Nr. 31-790-0785